

# Indoor Seating Comparison Chart



SPECIFICATIONS	JENX CORNER SEAT			DISCOVERY LUNAR	SPEX DISCOVERY HI-LO BASE	COBRA BASE		NEAT ZEAT	JENX ATOM	JENX JUNIOR+ SEAT			JENX MULTISEAT		GRAVITY CHAIR			
	Size 1	Size 2	Size 3			Size 1	Size 2			Size 1	Size 2	Size 3	Size 4	Size 5	Sml	Med	Lrg	XLrg
DESCRIPTION	Popular as the child's first seating system at floor level.			Kindy/Nursery Base for the Discovery Seat	High & Low level seating for the Discovery Seat	High & Low level seating for the Bingo or Spex Wonderseat		Modular indoor seating for activity/feeding	A modular seating system for users 6 months - 5 years	Comfortable, contoured and stable modular seating system, with outstanding levels of adjustment and a wide variety of support options.			Larger version of the Jenx Junior+ Seat, designed to suit older children and adolescents.		Supportive alternative seating and foam positioning system to augment postural control for a wide range of individuals. Also used as a 'hugging' chair for individuals who benefit from this type of sensory input.			
APPROX AGE OF USER*	6 mths-2 yrs	1-6 yrs	4-10 yrs	1-5 yrs	1-5 yrs	1-5 yrs	3-15 yrs	6 mths - Adolescent	6 mths - 5 yrs	18 mths - 5 yrs	3-8 yrs	5-12 yrs	8-14 yrs	10-18+ yrs	Refer to Measurements			
SEAT DEPTH (cm)	35	37	46	14-30	14-30	16-34	22-45	Kindy: 9-31 / Std: 18-36	9-34	21-30.5	24-32	27.5-39	28.5-45	35.5-52	26	32	40	45
SEAT WIDTH (cm)	33	39	49	13-27	13-27	18-30	23-35	40	18-27	15-37.5			20-44		27 front 20 cm back	30.5 cm front 23 cm back	37 cm front 27 cm back	40 cm front 30 cm back
SHOULDER HEIGHT (cm)	31-42	44-62	50-68	24-39	24-39	28-46	34-63	N/A	27-38	31-41	33-43	39-52	46.5-59	52.5-66	60cm total back	69cm total back	74cm total back	80cm total back
FOOTPLATE HEIGHT RANGE (cm)	N/A			14-28	14-28	16-37	16-40	8-56	16.5-32.5 (Fixed Angle Footplate) 11.5-36 (Flip Up Footplate)	19.5-27	25-36.5	25-40.5	25-47		24 cm to floor	26 cm to floor	28 cm to floor	30 cm to floor
MAX LOAD (kg)	25 kg	35 kg	55 kg	25	25	50		75	25	50			100		50	70	100	130
APPROX PRODUCT WEIGHT (kg)	Approx. 16			Base only 2.5	Manual: Base only 15.7 kg Electric: Base only 17.3 kg	11.3	11.5	Variable	13	~20kg (variable)			Variable		4	5.2	7.4	9.6
TILT IN SPACE	✗			40°	40°	45°		✗	✓	-10° - +30°			-10° - +30°		Use wedge underneath			
RECLINE	✗			90°-160°	90°-160°	90°-160°		✗	-15° - +15°	-15° - +15°			Adjustable Angle		✗			
GROWTH ADJUST SEAT	✓			✓	✓	✓		✓	✓	✓			✓		✗			
SEAT HEIGHT (cm)	N/A			22-32	Manual: 38.5 - 74.8 Electric: 38.5 - 63.5	27-61		Variable	28.5-75	29-76	33-76	36.5-76	33-57		N/A			
HI-LOW ACTIVATION	N/A			Tool Adjustment	Motorised or Manual Gas Assist	Hand Operated Gas Assist		N/A	Hydraulic Foot Pump	Hydraulic Foot Pump			Hydraulic Foot Pump		N/A			
TRAY	✓			✓	✓	✓		✓	✓	✓			✓		✓			
POSTURAL SUPPORT LEVEL	Mild/Moderate			Mild/Moderate/Complex	Mild/Moderate/Complex	Moderate/Complex		Mild/Moderate	Mild/Moderate/Complex	Mild/Moderate/Complex			Mild/Moderate/Complex		Mild/Moderate/Complex			

\* Indicator only, for accurate product selection client measurements are the best guide.

Medifab is as diligent as possible in compiling this information. Medifab does not guarantee the correctness and completeness of the information provided. Equally, Medifab does not guarantee that this information is up to date. For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Medifab. Product specifications are subject to change without notice.