Houdini 13 Harness

- Can be fitted on wheelchairs and some easy chairs as well as using in vehicles
- All buckles are **behind the seat** so the user cannot release them.



Model	Code	Age Group	Buckle
13J	4604-1302-017	3 to 8 years	Anti-Escape
13H	4604-1303-017	7 yrs to Adult	Anti-Escape

Houdini 22 Harness

- Can be used on **single** or bench seats
- Fit to the user before entering vehicle
- Good option for **taxi's** and also suitable for hip spicas.





TAXI

Model	Code	Age Group	Buckle
22C	4604-2211-017	3 to 8 years	Anti-Escape
22D	4604-2212-017	7 yrs to Adult	Anti-Escape

Houdini 27 Harness

- Centre-release buckle to provide postural support, anti-escape to prevent removal
- Attachment straps fasten horizontally around seat
- Suitable for single backrest seats and air travel*





Model	Centre-Release Buckle	Anti-Escape	Age Group
27A	4604-2702-017	4604-2702-026	3 to 8 yrs
27B	4604-2703-017	4604-2703-026	7 yrs to Adult

MODEL 13 BENEFITS

Suitable for use with people who are confused or have challenging behaviour, as all the buckles are behind the seat back.

In vehicles, this harness is suitable for use on a **single seat only**, and requires a gap between the vehicle seat and backrest for installation.

Not a safety harness: must be worn in conjunction with vehicle safety belt.

Upgrade all 3 buckles (neck strap and 2 on waist band) Pen/key release 4604-1305-031

MODEL 22 BENEFITS

Consists of a removable torso section which is attached to the individual prior to entering the vehicle.

Once the individual is seated, two buckles are engaged to provide safe travel.

Not a safety harness: must be worn in

conjunction with vehicle safety belt.

Additional seat attachment strap for use in another vehicle 4604-2222-017

Upgrade all 4 buckles (crotch strap, upper neck strap and 2 on waist band): Pen/key release 4604-2215-031 Magnetic key release 4604-2216-029

MODEL 27 BENEFITS

Designed to offer postural support and a certain degree of restraint as it is a full 5-point harness. Identical to Houdini 28; however straps fasten horizontally (28's fasten vertically).

Only suitable for single-backed seats.

Not a safety harness: must be worn in conjunction with vehicle safety belt.

Upgrade buckle on waist belt: Pen/key release 4604-2805-031 Magnetic key release 4604-2806-029

* Signed approval must be obtained from the airline prior to travel.



Houdini 28 Harness

- Medifab's **most popular** model of Houdini Harness
- Attachment straps fasten vertically around seat
- Centre-release buckle to provide postural support, anti-escape to prevent removal.





Model	Code	Age Group	Buckle
28	4604-2803-017	Junior to Adult	Centre-Release
28SB	4604-2803-026	Junior to Adult	Anti-Escape

Houdini 31 Harness

- Most secure harness, with twin anti-escape buckles
- Can be used on deep 'bucket' single or bench seat.
- Can be fitted to some regular booster seats - if straps can be attached to backrest





Model	Code	Age Group	Buckle
31A	4604-3101-017	3 to 8 yrs	Anti-Escape
31B	4604-3102-017	7 yrs to Adult	Anti-Escape

BUCKLE OPTIONS



Anti-Escape

Successfully prevents those with inquisitive fingers, challenging behaviour or learning difficulties from releasing the buckle. Operates in reverse to normal, as it is disengaged by raising a spring-loaded lever and will not press release.



Pen/Key Release

Requires a tool such as a pointed pen or key to prod the buckle release mechanism to open. Medifab supply a release tool with each of these buckles



Magnet Release

Magnetically released buckles require a particular magnet tool to disengage the buckle cap. Most advanced secure buckle option.

MODEL 28 BENEFITS

Designed to offer postural support and a certain degree of restraint as it is a full 5-point harness. This harness is the same as the Model 27, however the seat attachment straps go vertically around the vehicle seat.

This is suitable for **single or bench seats**. For seats that 'narrow' towards the top of the backrest, the Model 27 is more suitable.

Not a safety harness: must be worn in conjunction with vehicle safety belt.

Upgrade buckle on waist belt: Pen/key release 4604-2805-031 Magnetic key release 4604-2806-029

MODEL 31 BENEFITS

This harness offers more support due to a front PVC pad that is attached around the individual's torso with two straps fitted with anti-escape buckles.

Not a safety harness: must be worn in conjunction with vehicle safety belt.

Upgrade both buckles on waist band: Pen/key release 4604-3105-031 Magnetic key release 4604-3106-029

CROSS STRAP

Fit to a harness to hold shoulder straps together & prevent them slipping off.



PVC with side-release buckle 4604-3112-019 **PVC with anti-escape** 4604-3004-026

PVC with pen release 4604-3004-025



Cross Strap Padded Sleeve (side release buckle) Sml 4604-3117-017 Lge 4604-3118-017



Designs and styles may vary.