


# Jenx Atom

Modern, child-centred styling, combined with a comprehensive range of support options make Atom the ideal seating system to develop correct seating posture and function.

A modular seating system for users 6 months - 5 years



Different **base options** with a variety of **fixed or adjustable heights** allow the flexibility for each child to enjoy activities from floor level to breakfast bar height

**Tool free adjustments**, cleverly designed wing knobs offer a comfortable grip to enable effective clamping with minimal effort

**Unrivalled seat depth growth** to accommodate early intervention requirements

**Open sides** ensures adjustments are easily made with the user remaining seated

**Memory foam ischial well** designed to nestle the pelvis. Provides comfort, pressure distribution and stability to promote correct pelvic alignment

### Y-Base

- Versatile and flexible, with lockable swivel castors
- High-low range of 25.8-75 cm seat height from floor
- Tilt from 30° recline - 10° prone
- Base footprint of 65 x 78 cm

3615-8740-000




## FUNCTIONAL BENEFITS



### Head Positioning

Through our innovative Multigrip Head Supports, we offer a number of solutions that adapt to changing head positioning needs.



### Trunk Positioning

All trunk positioning accessories are optional and are intended to offer just the right level of support.



### Leg Positioning

With a number of different control options windsweeping, abduction and adduction can be easily influenced or controlled using our standard leg positioning accessories.



### Foot Positioning

The child's feet can be precisely positioned using a number of standard accessories.

Code	Approx Age Range	Max User Weight	Tilt	Back Pad Height	Chest Width	Hip Width	Seat Depth	Lower Leg Length	Foot Support Angle
3431-0001-022	6mths - 5yrs	25 kg	15° - 15°	27-38cm	20-28cm	18-27cm	9-34cm	16.5-32.5cm	5° - 35°

# Accessories

## Wooden Base

High 3431-8720-000  
Low 3431-8710-000



## Multigrip Head Support

1 Link 3431-4601-022  
2 Link 3431-4602-022  
3 Link 3431-4603-022



## Shoulder Protractors

3431-2750-022



## Lumbar Pad

3431-2501-000



## Back Pad Infill

3431-2370-000



## Thoracic Support Pads

High-Rise 3430-3151-000  
High-Rise Flexible 3431-3251-022



## Arm Supports

3431-5501-000  
Arm Support Covers 3431-5510-022



## Waistcoat Anterior Chest Support

3615-6701-000



## Nursery Tray

3431-7101-022  
Infill 3431-7182-000  
Half Padding 3431-7442-000  
Full Padding 3431-7452-000



## Elbow Blocks

3730-7350-000



## Activity Frame

3615-7475-000  
Grab Rail 3615-7200-000



## 3 Point Lap Strap

3431-6251-000  
Spex 4 Point Lap Strap 1406-6420-027



## Abduction Block

0 3615-1650-000  
1 3615-1651-000



## Long Leg Sitting Board

3431-0281-022



## Leg Control Pads

3431-1751-000



## Fixed Angle Foot Support

3431-5302-000  
Flip Up Angle Foot Support 3431-5301-000



## Foot Support Pad

3431-5701-000



## Sandals

5356-6081-000



## Ankle Snugs

5356-6131-000



## Sandal Raising Blocks

5356-6091-000

