

# Jenx Therapy Bench

Ideal platform for a variety of different activities such as sitting and balance therapy, dynamic sitting balance and sitting posture improvement

Therapy aid designed to promote upright sitting and balance



Easily **height adjustable** using tightening bolt.


Easy to clean, antimicrobial, wipe-down vinyl cushions

**Interchangeable cushions** include the standard 1" flat cushion and the half round cushion

The simple design means that the Therapy Bench can be **easily moved** with just 1 hand

## ACCESSORIES


**Half Round Cushion**  
Size 1 7510-0501-054  
Size 2 7510-0502-054




**Pelvic Support**  
7510-1900-054



**Abduction Block**  
7510-1750-054



**Knee Blocks**  
Size 1 7510-1601-054  
Size 2 7510-1602-054



Size	Code	Max Load Weight	Bench Height from Floor	Bench Top Length	Bench Top Width	Flat Cushion Width	Flat Cushion Length	Flat Cushion Depth
Size 1 - Full Length	7510-0010-054	85kg	165-250mm	980mm	200mm	250mm	1000mm	25mm
Size 1 - Half Length	7510-0015-054	85kg	165-250mm	480mm	200mm	250mm	500mm	25mm
Size 2 - Full Length	7510-0020-054	85kg	235-395mm	980mm	250mm	300mm	1000mm	25mm
Size 2 - Half Length	7510-0025-054	85kg	235-395mm	480mm	250mm	300mm	500mm	25mm