

# BINGO Mini



## EARLY INTERVENTION ADAPTIVE STROLLER

Distributed by  
**medifab**

AUSTRALIAN  
SCRIPTING AND ORDERING GUIDE

# Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

# Product Details

## Bingo Mini

<b>Order Code, Size 0</b>	2880-0008-000
<b>Conversion Kit from Mini to Size 1 (stroller and Cobra)</b> Converts BINGO Mini frames to accept the BINGO Evo size 1 seat unit, extending the growth of the solution (size 1 seat unit must be purchased separately).	2880-8080-086

Seat Unit		Mobility Base with swivelling wheels	
Seat width	13-25cm	Push handle height	75-110cm
Seat depth	12-27cm	Overall width	63cm
Lower leg length	16-36cm	Wheel diameter front (swivel)	19cm
Back support height	37-51cm	Wheel diameter rear	29cm
Height of slots for shoulder straps	25.5-39.5cm	Tilt in space (forward and rear facing)	45°
Seat to floor height (dependant on tilt position of seat)	50cm	Folded size mobility base (incl. rear wheels)	80 x 63 x 38cm
Backrest recline	80-180°	Folded size mobility base (excl. rear wheels)	80 x 63 x 31cm
Weight of seat unit	7.8kg	Weight mobility base (swivelling front wheels)	9.8kg
Load capacity of seat unit	25kg	Load capacity mobility base	60kg
Weight of seat unit and mobility base combined	17.6kg		

\*Including upholstery and standard accessories




### Stroller package includes:

- Matte black foldable stroller frame with black wheels
- Removable seat unit
- Tie-downs for transport
- Removable hip supports (width & depth adjustable)
- Fabric bag (max. 3kg)
- Adjustable suspension set
- Receiver for grab rail or therapy tray
- Elevating foot support hanger (angle adjustable 85-180°)
- **Note: Upholstery not included in package.**



Items outlined in blue are required as standard package

## Upholstery Colour Options

Berry Soft		Turquoise Soft		Black	
Seat Upholstery (required)					
2880-8800-347		2880-8800-349		2880-8800-346	

Berry Soft		Turquoise Soft		Black	
Hip Support Upholstery (required)					
2880-1074-347		2880-1074-349		2880-1074-346	

## Upholstery Upgrades

### Upholstery Pressure Relief Foam Upgrades

Seat  
2880-0391-000

Back  
2880-2392-000

## Head Supports

### Head Support Pads

2880-4470-346



# Head Supports (continued)

Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
<b>Merino/Spacer Fabric</b> 2010-4110-099		Small (17cm Internal width)		Small (17cm Internal width)	
		<b>Merino/Spacer Fabric</b> 2010-4210-099		<b>Merino/Spacer Fabric</b> 2010-4310-099	
		<b>Wipe Down Cover</b> 2010-4210-209		<b>Wipe Down Cover</b> 2010-4310-209	
<b>Wipe Down Cover</b> 2010-4110-209		Large (21cm Internal width)		Large (21cm Internal width)	
		<b>Merino/Spacer Fabric</b> 2010-4211-099		<b>Merino/Spacer Fabric</b> 2010-4311-099	
		<b>Wipe Down Cover</b> 2010-4211-209		<b>Wipe Down Cover</b> 2010-4311-209	

# Lateral Supports

Standard Lateral Supports		Spex Narrow Straight Lateral Supports	
<b>Size 0</b> 80mm (H) x 100mm (D) 2880-3009-346		<b>Size 0</b> 85mm (H) x 115mm (D) 2880-3007-100	
Spex Broad Lateral Supports		Spex Teardrop Lateral Supports	
<b>Size 1</b> 105mm (H) x 130mm (D) 2880-3007-310		<b>Size 0</b> 150mm (H) x 105mm (D) 2880-3007-500	

# Sacral Support

Sacral foam insert including Spacer	
On velcro attachment. Also includes foam spacer which adheres behind Sacral insert where required. 10cm (H) x 18cm (W) 2880-2510-000	

# Seat Options

Thigh Wedge		Trochanter Wedges	
Velcro attached under upholstery. 22.5cm (W) x 10.5cm (D). Height: 2cm tapering down to 1cm 2006-0505-000		Option to contour standard flat cushion. Velcro attached under upholstery. 2006-0660-000	

# Abductor Pommel

## Swing down abductor pommel

### Size 1

For increased abduction of the legs.  
9cm (D) x 6.5cm (W) x 9cm (H)  
2880-1651-346



# Hip Belts & Harnesses

## Spex 5-Point Harness

Must be worn when stroller is used for vehicle transportation.  
2006-6551-000



## Wonderseat Harness

Wonderseat Padded Centrepoint Shoulder Harness.

Size 0, Pad height 22cm  
1409-5660-017



## Wonderseat Hip Belt


Wonderseat 4-Point Padded Hip Belt.  
Pad width 24cm.  
1406-5420-027



# Foot Positioning

## Footplate Cover

Neoprene fabric protects feet from the hard footplate surface.  
2880-5305-000



## Footplate Spacer

Spacer block to reduce lower leg length by 4cm.  
2880-5321-000



## Ankle Footfasts

Strap inner circumference 15-23cm.  
2880-6101-000



# Tray & Grab Rail Options

## Therapy Tray

Cannot be used at the same time as Grab Rail.  
2880-7110-000



## Therapy Tray Infill

For smaller trunk cut out.  
2880-7180-000




## Grab Rail


Cannot be used at the same time as Therapy Tray.  
2880-7250-346






# Critical Care Options


IV Pole	Oxygen Bottle Holder
<p>With mounting bracket for attachment to mobility base. Fits all sizes. 2880-8275-000</p>	<p>With mounting bracket for attachment to mobility base. 120mm diameter for "C Type" bottle. Fits all sizes. 2880-8253-000</p>
<p><b>Extra Bracket</b> For QR mounting of IV pole to second base. (Cobra or other indoor base). 2880-8275-050</p>	

Chassis Mount Tray	
<p><b>Size 1</b> 2880-8058-000</p>	


## Canopy

Includes mini mounts and rain cover	
<p>2880-7501-346</p>	

## Hand Strap

Wrist safety tether strap	
<p>2005-6930-000</p>	

## Swivel Lock

To hold castors straight when brisk walking	
<p>2880-8440-000</p>	

# Cobra Indoor Base

## Cobra Hi-Low Indoor Frame with 75mm Twin Castors

The seat unit easily transfers out of the stroller base and into the Cobra indoor base for feeding and activity times. Lower down the floor for playing and activities, and raise up to a table for meal time! Tilt settings can still be adjusted in the same way as the stroller base.

**Size 0**  
2880-8720-000



# Separate Seat Unit

## Use with Cobra Base

Separate seat unit for use with the Cobra Base if a combined stroller and seat unit is not purchased. This is also a great option if two seat units are preferred to avoid swapping the seat unit between the stroller and Cobra bases. Must be purchased with upholstery.

**One Size**  
2880-0010-000



# Twin Bingo Stroller

Connection mount that allows two strollers to be connected one in front of the other for an additional sibling. The strollers are able to be easily detached and used separately where required, providing flexibility when travelling with or without siblings.

This is not available for use with a size 3 Bingo Evo. This can be used to connect either: bases of the same size (and corresponding sized Bingo Spex strollers), or the Bingo Mini and Evo Size 1. It is not suitable for use between any other size combinations.

Please fill out an additional script form to order the 2nd stroller unit with a unique configuration and accessories.

Please note that this can only be fitted at the time of ordering both stroller units (Medifab cannot retrofit this attachment).

2880-8630-000



## Attendant Brake

Recommended for purchase with Twin Stroller configuration due to the increased combined weight and inertia when both strollers are connected.

2880-8452-000

# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: BINGO Mini\_20240321\_AU