

Carrot XL

CARROT 3

Carrot XL caters for safety and the importance of keeping the user in a supported posture for travel. This larger sized car seat is suitable for older children and adult users who will not fit into the Carrot 3000.

Larger version of the popular Carrot 3000 Car Seat



Extendable shoulder protector wings hold and protect the body from side impacts and centrifugal forces

Ergonomical large head support

Special large anterior chest support design with easy length adjustment to support the shoulders of older users

Correct pelvis alignment with comfort padded hip belt

Low profile wide seat base

UPGRADE KIT

Upgrade an existing Carrot 3000 to XL.

Carrot XL Upgrade kit can be purchased to **upgrade an existing Carrot 3000 Car Seat** to an XL version, instead of purchasing a whole new car seat when the child grows out of the existing seat.



Kit includes:

- XL Low profile wide Seat Base
- Longer metal rods for base
- H-Anterior Chest Support
- Hip belt

4116-0030-259

Please note that the XL Upgrade Kit must be fitted by Medifab Ltd or an authorised Medifab distributor.

Package includes

- Low contour seat base
- Back support
- 50mm height extension
- Deep headrest
- Width adjustable shoulder protector wings
- 4 point H anterior chest support
- Package available with **Standard Buckle** (4116-0020-259), **Pen-Key Release Buckle** (4116-0020-025), or **Anti-Escape Buckle** (4116-0020-026)



FUNCTIONAL BENEFITS

- Equipped as standard with the low profile wide seat base, large size headrest, shoulder protector wings, and 50mm back support extension.
- Extendable shoulder protector wings hold and protect the body from side impacts and centrifugal forces.
- Tether options include additional security tether for the top of back support or ISOFIX latch connectors. These may be required by your local authorities.

If user is too small for the Carrot XL Car Seat, the Carrot 3000 Car Seat is available for smaller users.

The Carrot XL conforms to regulation (EU) 2017/745 on medical devices.

Code	Age Range	User Weight Range	User Height Range	Seat Depth	Seat Width	Shoulder Height	Seat Weight
Standard: 4116-0020-259	10 to Adult	36-75kg	140cm - 175cm	38 or 42cm	48cm	49.5 - 62cm	10.5kg
Pen-Key Release: 4116-0020-025							
Anti-Escape: 4116-0020-026							

Accessories

Tether Strap
4116-6953-017



ISOFIX Lower Straps
4116-6910-000



Pen-Key Release Harness Belt
4116-6125-256



Tilt Wedge
4116-0291-259



100mm Back Support Height Extension Infill Pad
4116-2482-259



Forearm Support Tray
4116-5544-259



Lateral Trunk Supports
4116-3250-090




Air Mesh Cushion Pad
4116-8880-099



Hipbelt (Seat Base)
Medium 4116-6532-027
Large 4116-6533-027



Crotch Belt
(requires customisation of seat base)
4116-6320-000



Airplane Fixing Strap
4116-6960-017

