

# Stroller Comparison Chart



SPECIFICATIONS	BINGO MINI	BINGO EVO			BINGO SPEX		PIXI		DURO		KOOL STRIDE	MULTIROLLER
		Size 1	Size 2	Size 3	Size 1	Size 2	Size 1	Size 2	Size 1	Size 2		
DESCRIPTION	A smart stroller base with ultimate early-intervention seating.	A versatile and modular stroller system for users with either basic or critical care needs.			A smart-looking stroller with superior postural support that enables parents to take their child outside and explore.		Light and nimble stroller with comfortable rigid seating.		Configurable stroller for everyday challenges, with a great range of support and positioning accessories.		Go off road and beyond!	Incredibly light, incredibly tough!
APPROX AGE OF USER*	4mths-4yrs	1-5yrs	3-12yrs	8-15yrs	2-5yrs	3-7yrs	3-7yrs	6-14yrs	3-12yrs	8-15yrs	8-21yrs	8-18
MAX LOAD (kg)	25	50			50		35	50	75		70	75
SEAT WIDTH (cm)	13-25	18-30	23-35	23-35	17-28	22-33	28	33	23-34	28-39	43	30-38
SEAT DEPTH (cm)	12-27	16-34	22-40	30-45	18-30	23-38	25 or 30	35 or 40	26-40	31-45	30	33
SHOULDER HEIGHT (cm)	25.5-39.5	28-46	34-55	42-63	36-44	40-48	54 or 60	64 or 70	32-40	42-62	60	29, 38, 48
FOOT SUPPORT HEIGHT RANGE (cm)	16-34	19-40	19-40	19-40	25-40	25-45	13-35	13-39	20-44	20-44	45	20 - 30
FRONT/REAR FACING	Both	Both			Both		Front		Front		Front	Front
INDOOR FRAME OPTIONS	✓	✓			✓		✗		✗		✗	✗
APPROX PRODUCT WEIGHT (kg)	17.5	15	16.1	18	29	30	12.5	15	19.2	21	15	9
TILT IN SPACE	45°	45°			45°		✗		Static 16° or 20°		✗	✗
RECLINE	80°-180° in 10° steps	90°-160°		90°-135°	90°-160°		90°-125°		78-120°		✗	✗
FOOT SUPPORT ELEVATION	✓	✓			✓		✗		✓		✗	✗
GROWTH ADJUST SEAT	✓	✓			✓		5cm growth bracket, seat depth and back height		✓		✗	✓
SWIVEL FRONT WHEELS	✓	✓			✓		✓		✓		✗	✓
ACTIVITY TRAY	✓	✓			✓		✓		✓		✗	✗
SUN CANOPY	✓	✓			✓		✓		✓		✓	✓
RAIN SHIELD	✓	✓			✓		✓		✓		✓	✗
BASKET STORAGE	✓	✓			✓		✓		✓		✓	✗
POSTURAL SUPPORT LEVEL	Moderate/Complex	Moderate/Complex			Complex		Moderate		Moderate		Mild	Mild
TYRE SOLID / PNEUMATIC	Solid	Solid			Solid		Combination		Solid		Pneumatic	Pneumatic
TRANSPORT APPROVED	✓	✓			✓		✓		✓		✗	✗

\*Indicator only, for accurate product selection client measurements are the best guide.