

Jenx Atom



MODULAR SEATING SYSTEM

Distributed by
medifab


SCRIPTING AND ORDERING GUIDE

October 2023


Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Measurements	User Measurements
A Hip width	
B Seat depth	
C Seat to top of shoulder	
D Seat to top of head	
E Chest width	
F Seat to axilla	
G Leg length	





Product Specifications	
Average age	6 Months-5 Years
Weight capacity	25kg
Seat width	18cm-27cm
Seat depth	9cm-34cm
Chest support width	20cm-28cm
Seat to top of backpad height	27cm-38cm
Foot support angle	-20° to 15°
Foot support height	11cm-36cm
Back support angle	-75° to 105°
Tilt range	30° posterior, 10° anterior



Seat Unit/Colour Choice

Base Options


Red	Pink	Y-Base
<p>Includes:</p> <ul style="list-style-type: none"> Padded seat with memory foam Shaped back pad Adjustable flip away hip support pads Push handle Upper back pad 	<p>Includes:</p> <ul style="list-style-type: none"> Padded seat with memory foam Shaped back pad Adjustable flip away hip support pads Push handle Upper back pad Arm supports Tray table support Shoulder protractors High rise flexi thoracic pads Thoracic pad covers 2 Link multigrip head support Flip up angle adjustable footplate 	<p>Black Y-base with lockable swivel castors Floor to seat height: 30-75cm Tilt range: 30° posterior 10° anterior</p> 
3431-0001-022	3431-0010-028	3615-8740-000
		

Head Supports

1 Link Multigrip	2 Link Multigrip		
<p>Single link wings for some degree of lateral support. 33cm W x 5cm D x 7.5cm H 5356-4601-000</p>	<p>2 Links in each wing for extra support and positioning. 47.5cm W x 5cm D x 12cm H 5356-4602-000</p>		
<p>Wipeable Cover (Optional) 5356-4671-000</p>	<p>Wipeable Cover (Optional) 5356-4672-000</p>		

Back Pad Options





Lumbar Pad		Back Pad Infill	
<p>Provides additional lumbar support. 3431-2501-000</p>		<p>A 50mm back pad to reduce the seat depth for smaller children. 3431-2370-000</p>	

Washable Cover	
<p>For Seat and Back Support. Comes in Red. 3431-2470-022</p>	

Shoulder & Trunk Supports

Shoulder Protractors	High-Rise Short Thoracic Pad	High-Rise Flexible Thoracic Pads
<p>Flip away when not required. Independently adjustable for width and angle. 13cm High x 15.5cm Deep 3431-2750-022</p>	<p>Short, slim trunk support pads with single-handed flip away action. No strap. 3430-3151-000</p>	<p>Slimline and super flexible trunk support pads with straps and single-handed flip away action. 8cm High x 18cm Deep 3431-3251-022</p>
		
	<p>Cover (optional) 3615-3171-022</p>	<p>Cover (required) 3431-3271-022</p>

Harnesses & Hip Belt

Waistcoat Harness	Wonderseat Harness	3 Point Lap Strap	4-Point Hip Belt
Provides additional trunk support to children with low tone. 3615-6701-000	Size 0, Pad height 22cm Padded Centrepoint Shoulder Harness. 1409-5660-017	Dual lapstrap with additional safety of a groin harness. 3431-6251-000	Wonderseat Padded Hip Belt, Pad width 24cm. 1406-5420-027
			

Leg Supports

Leg Positioning Pads	Long Leg Sitting Board	Abduction Block	Abduction Block
Ideal for controlling windsweeping or creating an abducted seating position. Required fixed angle foot support 3431-1751-000	(Red) For very young children or those needing hamstring stretch. Used in place of foot support. 3431-0281-022	Size 0 Length: 9cm Height above Seat: 8.5cm Overall height: 14.5cm Width Front: 9cm Width Rear: 7cm 3615-1650-000	Size 1 Length: 18cm Height above Seat: 7.5cm Overall height: 15.5cm Width Front: 10cm Width Rear: 8.5cm 3615-1651-000
			



Foot Supports

Flip-Up Angle Foot Support	Foot Support Pad
Effortlessly flips away to aid in transfers. 11.5cm to 36cm of Leg-Length Adjustment 3431-5301-000	Attaches to foot support for softer surface. Can only be used on the flip-up foot support. 3cm Thick 3431-5701-000
	



Foot Supports (cont.)

Sandals		Sandal Raising Blocks		Ankle Snugs	
<p>Ideal for positioning feet where a firmer hold on the foot is required. 5356-6081-000</p>		<p>Easily accommodate leg length discrepancies by giving sandals extra height where needed. 5356-6091-000</p>		<p>Features unique magnetic buckle fastening mechanism. Sandals required. 5356-6131-000</p>	
8cm W x 17cm L		4 x 1.5cm Blocks		Inner Circumference 12cm - 22.5cm	



Arm Supports

Arm Supports		Covers for Arm Supports	
<p>Required for using a Nursery Tray. 3431-5501-022</p>		<p>(Red) Provide soft top for arm supports when tray is not used. 3431-5510-022</p>	

Tray & Accessories

Nursery Tray		Nursery Tray Infill	
<p>A useful surface for work, feeding or play. Requires arm supports. 3431-7101-022</p>		<p>Designed to work with the Nursery Tray to reduce the cut-out for smaller children. 3431-7182-000</p>	
58cm W x 48cm D Internal Cutout - 31.5cm W x 17.5cm D			

Nursery Tray Padding		Elbow Blocks	
<p>Half Padding 3431-7442-000</p>		<p>Comes as a pair. Mounts to Nursery Tray corners for elbow support. 9.5cm High x 10cm Wide 3430-7350-000</p>	
<p>Full Padding 3431-7452-000</p>			

Activity Frame		Grab Bar	
<p>Great for encouraging head control and hand/eye co-ordination. Requires tray. 54cm W x 16cm D x 35cm H 3615-7475-000</p>		<p>Nursery Tray required. Bar dimensions: 54cm (W) x 9.5cm (D) x 7cm (H) 3615-7200-000</p>	

medifab

Distributed by Medifab Ltd

22 Detroit Drive
Rolleston 7675
New Zealand

0800 543 343

26 Pardoe Street
Devonport TAS 7310
Australia

1300 543 343

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Jenx Atom_20231024