

Jenx Multistander



**PRONE OR SUPINE
STANDING AID**


Distributed by
medifab

**AUSTRALIAN
SCRIPTING AND ORDERING GUIDE**

December 2023

Client Details

Client/User Name:	
Date:	
NDIS/Funding:	
User Measurements	
A Hip width	
B Seat depth	
C Seat to top of shoulder	
D Seat to top of head	
E Chest width	
F Seat to axilla	
G Leg length	



Medifab Notes

This image shows a single sheet of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Product Details

Product Specifications	Size 1	Size 1 - With Abduction	Size 2
Average age	9 Months – 6 Years	9 Months – 6 Years	3 – 13 Years
Weight capacity (kg)	35	35	60
Hip width (cm)	14 – 27.5	14 – 27.5	20 – 32
Hip height (cm)	25 – 44	25 – 44	37.5 – 60
Chest width (cm)	14 – 27.5	14 – 27.5	20 – 32
Chest height (cm)	49 – 85	50 – 89	75 – 120
Height range (cm)	62 – 116	60 – 116	110 – 155
Tilt range (°) (90 = upright)	10 – 90	10 – 90	0 – 85 (Powered) / 10 – 85 (Manual)
Frame width (cm)	61.5	61.5	67.5
Frame length (cm)	76	76	107
Abduction range (°)	-	30 (Either Leg)	-
Height of surface when horizontal (cm)	56	56	77.5

*Max ranges are achieved with additional accessories

Jenx Multistander Package Includes:

- Multistander manual or powered frame
- Base with 4 lockable swivel castors
- Knee block mounting brackets
- Angle adjustable footplate

Key

Blue	Suitable for Prone and Supine
Green	Suitable for Supine only
Yellow	Suitable for Prone only




Items outlined in blue are recommended as part of standard package


Power & Size Options

Manual		Battery Powered	
Size 1 5230-0001-000		Size 2 5230-0008-000	
Size 2 5230-0002-000			


Colour Options

Blue		Lime	
-066		-054	

Abduction Accessory

Only Available on Size 1	
Supports 0-30 degrees in each leg in 5 degree increments. Fitted prior to shipping. 5230-5000-000	



Support Boards

Size 1	Size 2	
PU chest & hip support boards. 450mm Wide 5230-2201-XXX	PU chest & hip support boards. 520mm Wide 5230-2202-XXX	

Head Supports




Moulded Head Support		Oval Head Support	
Moulded polyurethane headrest for children who require minimal to moderate head support. 280mm wide x 210mm high 5230-4101-XXX		A simple support for children who require minimal head positioning. 190mm wide x 100mm high 5356-4611-XXX	

Key Blue Prone and Supine Green Supine Only Yellow Prone Only

Multigrip - 1 Link	Multigrip - 2 Link	Multigrip - 3 Link	
Single link wings for some degree of lateral support. 33cm W x 5cm D x 7.5cm H 5356-4601-000	2 Links in each wing for extra support and positioning. 47.5cm W x 5cm D x 12cm H 5356-4602-000 	3 adjustable links in each wing for ultimate lateral headrest control and comfort. 62cm W x 5cm D x 12cm H 5356-4603-000	
Wipeable Cover (Optional) 5356-4671-000	Wipeable Cover (Optional) 5356-4672-000	Wipeable Cover (Optional) 5356-4673-000	

Pelvic Supports

Support Pad & Straps Small	Support Pad & Straps Medium	
A 100mm B 90mm C 90mm D 120mm E 25mm 5356-3221-XXX	A 100mm B 130mm C 90mm D 140mm E 25mm 5356-3222-XXX	


Flexi Thoracics (Pair) Size 1	Flexi Thoracics (Pair) Size 1 Long	Flexi Thoracics (Pair) Size 2	 
A 180mm B 80mm 5356-3333-000	A 215mm B 80mm 5356-3334-000 	A 220mm B 110mm 5356-3335-000	
Covers Required 5356-3363-XXX	Covers Required 5356-3364-XXX	Covers Required 5356-3365-XXX	

Hi-Rise Thoracic Laterals

Adjustable for height, width and angle.

A 90mm **B** 90mm
C 200mm **D** 120mm
E 25mm


5230-3251-XXX





Trunk Supports

Support Pad & Straps Small	Support Pad & Straps Medium	
A 100mm B 90mm C 90mm D 120mm E 25mm 5356-3221-XXX	A 100mm B 130mm C 90mm D 140mm E 25mm 5356-3222-XXX	

Trunk Supports (cont.)

Flexi Thoracics (Pair) Size 1	Flexi Thoracics (Pair) Size 1 Long	Flexi Thoracics (Pair) Size 2	
A 180mm B 80mm 5356-3333-000	A 215mm B 80mm 5356-3334-000	A 220mm B 110mm 5356-3335-000	
Covers Required 5356-3363-XXX	Covers Required 5356-3364-XXX	Covers Required 5356-3365-XXX	

Hi-Rise Thoracic Laterals	Multigrip Body Support		
Adjustable for height, width and angle. A 90mm B 90mm C 200mm D 120mm E 25mm 5230-3251-XXX		For Size 2 only. Provides extra support on shoulders or lumbar area. 420mm Wide 5230-2721-XXX	

Sternal Prompt
For Prone 5230-2910-000
Cover required 5230-8801-XXX



Waistcoat Harness

Size 2	Size 3	Size 4	Mounting Kit
Provides additional trunk and shoulder stability			Mounting bracket for Harness (Required for Size 1 Multistander with Abduction Only). 5230-part-056
A 5cm B 32cm C 14cm D 15.5cm 5230-6712-000	A 6.5cm B 37cm C 17.5cm D 18cm 5230-6713-000	A 8.5cm B 39.5cm C 20.5cm D 18.5cm 5230-6714-000	



Knee Supports

Size 1	Size 2	
Easy to attach and adjust for height, angle & depth		
A 100mm / B 45mm / C 75mm 5356-1751-XXX	A 120mm / B 60mm / C 95mm 5356-1752-XXX	


Tray Options

Elbow Blocks	
<p>Suits all sizes. 10cm W x 14cm H x 8.5cm D 5356-7351-000</p>	

Adjustable Tray - Size 1	Adjustable Tray - Size 2	
<p>A 650mm B 450mm C 370mm D (With Insert) 330mm / (Without Insert) 290mm 5230-7111-000</p>	<p>A 700mm B 530mm C 460mm D (With Insert) 410mm / (Without Insert) 360mm 5230-7112-000</p>	

Foot Positioning

Sandals			
<p>Size 1 170mm long x 80mm wide 5356-6081-000</p>	<p>Size 2 220mm long x 95mm wide 5356-6082-000</p>	<p>Size 3 275mm long x 105mm wide 5356-6083-000</p>	

Sandal Raiser Blocks			
<p>Size 1 Kit of 4x15mm blocks 5356-6091-000</p>	<p>Size 2 Kit of 4x15mm blocks 5356-6092-000</p>	<p>Size 3 Kit of 4x15mm blocks 5356-6093-000</p>	

Ankle Snugs		
<p>Size 1 120mm - 225mm Circumference 120mm High 5356-6131-000</p>	<p>Size 2 165mm - 295mm Circumference 120mm High 5356-6132-000</p>	



Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Multistander_20231208_AU