

Jenx Supine



SUPINE STANDING SYSTEM

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

August 2024

Client Details

Client/User Name:

Date:

DOB:

User Measurements

A: Hip Width

B: Axilla Height


C: Chest Width

D: Overall User Height

E: Heel to Knee

F: Hip Height

User Weight



Product Specifications	Size 1 - Manual Tilt	Size 2 - Powered Tilt	Size 3 - Powered Tilt
Code	5350-0001-XXX	5350-0002-XXX	5353-0003-XXX
Average Age	9 Months - 4 Years	3 - 11 Years	10 - 18 Years
Weight Capacity	25 kg	50 kg	100 kg
Hip Width	13cm - 31cm	13cm - 31cm	22cm - 44cm
Chest Width	13cm - 31cm	13cm - 31cm	22cm - 44cm
Height	50cm - 106cm	90cm - 146cm	107cm - 190cm
Tilt Range	10° to 90°	10° to 90°	5° to 90°
Frame Width	57.5cm	65cm	75cm
Frame Length	66cm	99.5cm	107cm
Height of surface when horizontal	54cm	73cm	83cm

*The lateral postural supports can be used at any available slot on the supine stander.

Key

BLUE Suitable for Size 1

GREEN Suitable for Size 2

YELLOW Suitable for Size 3


Jenx Supine Package Includes:

- Frame with 4 lockable swivel castors (manual for size 1, powered for size 2 and 3)
- PU boards
- Height-adjustable foot plate (one-piece for size 1 and 2, split for size 3)







* Sections of board can be removed and added back on as user grows

Colour Options



Blue	Lime
	
XXXX-XXXX-066	XXXX-XXXX-054

Head Supports

Moulded Head Support	Head Pad Cushion with Side Supports
<p>Moulded polyurethane headrest for children who require minimal to moderate head support. 28cm wide x 21cm high 5353-4103-XXX</p> 	<p>Use to provide simple head support. 5350-4650-XXX</p>  <p>*Depending on size of stander</p>

Multigrip - 1 Link	Multigrip - 2 Link	Multigrip - 3 Link	
5350-4601-000	<p>For Size 1 & 2 5350-4602-000</p> <p>For Size 3 5353-4602-000</p>	5350-4603-000	  <p>A: 13cm / 5.1" B - 1 Link: 8cm / 3.1" C - 1 Link: 8cm / 3.1" D - 1 Link: 8cm / 3.1"</p>
Wipeable Cover (Optional) 5356-4671-000	Wipeable Cover (Optional) 5356-4672-000	Wipeable Cover (Optional) 5356-4673-000	

Trunk/Hip Supports


Support Pad & With Derotational Straps				
Small	Medium	Large	X-Large	
<p>A 10cm B 9cm C 9cm D 12cm E 2.5cm</p> <p>5356-3221-XXX</p>	<p>A 10cm B 13cm C 9cm D 14cm E 2.5cm</p> <p>5356-3222-XXX</p>	<p>A 20cm B 16cm C 9cm D 17cm E 2.5cm</p> <p>5353-3223-XXX</p>	<p>A 20cm B 26cm C 9cm D 27cm E 4cm</p> <p>5353-3224-XXX</p>	 

Key BLUE Suitable for Size 1 GREEN Suitable for Size 2 YELLOW Suitable for Size 3

Trunk/Hip Supports (cont.)

Flexible Support Pads


Size 1	Size 1 - Long	Size 2	Size 3
A 18cm B 8cm	A 21.5cm B 8cm	A 22cm B 11cm	A 32cm B 14cm
5356-3333-000	5356-3334-000	5356-3335-000	5356-3336-000
Cover (required) 5356-3363-XXX	Cover (required) 5356-3364-XXX	Cover (required) 5356-3365-XXX	Cover (required) 5356-3366-XXX



Waistcoat Harness

Provides additional trunk and shoulder stability

Size 2	Size 3	Size 4	Mounting bracket for Harness (Required)
A 5cm B 32cm C 14cm D 15.5cm	A 6.5cm B 37cm C 17.5cm D 18cm	A 8.5cm B 39.5cm C 20.5cm D 18.5cm	5350-part-060
5230-6712-000	5230-6713-000	5230-6714-000	



Abduction Block

Used to keep the body in midline position

Arched Foam Pads (1 pad per set - for use with abduction block)


Qty

12.5cm H x 5cm W x 26cm D 5353-1603-066		13.5cm W x 13.5 cm T x 5.5cm D Provides comfort and ability to increase abduction block width. Supplied individually. 5353-1775-000	
--	--	---	--

Knee Supports

Standard Knee Support with Derotational Straps

Small Blocks with Strap	Medium Blocks with Strap	XL Blocks with Size 1 Strap	XL Blocks with Size 2 Strap	XL Blocks with Size 3 Strap
A 10cm B 9cm C 9cm D 12cm E 2.5cm Strap Length 40cm Strap Width 12.5cm	A 10cm B 13cm C 9cm D 14cm E 2.5cm Strap Length 50cm Strap Width 12.5cm	A 20cm B 26cm C 9cm D 27cm E 4cm Strap Length 40cm Strap Width 18cm	A 20cm B 26cm C 9cm D 27cm E 4cm Strap Length 60cm Strap Width 20cm	A 20cm B 26cm C 9cm D 27cm E 4cm Strap Length 80cm Strap Width 22cm
5350-1761-XXX	5350-1762-XXX	5353-1765-XXX	5353-1766-XXX	5353-1767-XXX



Key BLUE Suitable for Size 1 GREEN Suitable for Size 2 YELLOW Suitable for Size 3

Knee Supports (cont.)

Flexion Control Knee Support with Derotational Straps

Includes backing plate and foam arched pads to build up support behind knee for children with knee flexion contractures. Use with PU spacer blocks below.

Medium	XL with Size 1 Strap	XL with Size 2 Strap	XL with Size 3 Strap	
A 10cm B 13cm C 9cm D 14cm E 2.5cm Strap Length 50cm Strap Width 12.5cm	A 20cm B 26cm C 9cm D 27cm E 4cm Strap Length 40cm Strap Width 18cm	A 20cm B 26cm C 9cm D 27cm E 4cm Strap Length 60cm Strap Width 20cm	A 20cm B 26cm C 9cm D 27cm E 4cm Strap Length 80cm Strap Width 22cm	
5350-1772-XXX	5353-1775-XXX	5353-1776-XXX	5353-1777-XXX	

PU Spacer Block (1 block per set)

Qty

Used in conjunction with Flexion Control Knee Supports or Abduction Block (size 3 only) to add additional spacing.


2.5cm Thick

5353-2730-000




Foot Positioning

Sandals

Size 1	Size 2	Size 3	Size 4	
8cm (W) x 17cm (L)	9.5cm (W) x 22cm (L)	10.5cm (W) x 27.5cm(L)	11.5cm (W) x 32cm(L)	
5356-6081-000	5356-6082-000	5356-6083-000	5356-6084-000	

Sandal Raising Blocks

To accommodate leg length discrepancies or for that extra bit of height! Please note, size 3 Supine Stander has individual foot-plates that are able to be off-set in height.

Size 1	Size 2	Size 3	Size 4	
4 x 1.5cm	4 x 1.5cm	4 x 1.5cm	4 x 1.5cm	
5356-6091-000	5356-6092-000	5356-6093-000	5356-6094-000	

Key BLUE Suitable for Size 1 GREEN Suitable for Size 2 YELLOW Suitable for Size 3

Foot Positioning (cont.)

Ankle Snugs for Sandals

Size 1

To wear with
Size 1 Sandals

12 - 22.5cm Internal
Circumference
12cm High

5356-6131-000

Size 2

To wear with Size
2 & 3 Sandals

16.5 - 29.5cm Internal
Circumference
12cm High

5356-6132-000




*Please note there are no ankle snugs for size 4 sandals

Tray Options

Adjustable Tray

Size 2 includes tray infill to cover the bowl if required


Size 2 with bowl



A 70cm / 27.6"
B 54cm / 21.3"
C 25.5cm / 10"
D 36cm / 14.2"
E 50cm / 19.7"

5350-7111-000

Size 3





A 73cm / 28.7"
B 63cm / 24.8"
C N/A
D 40cm / 15.7"
E 44cm / 17.3"

5353-7113-000

Tray Extension Wings

Increase support for the elbows and prevents arms from slipping off the tray

5353-7170-000




Elbow Blocks

Used to encourage arms into a
midline position

10cm W x 8.5cm D x 14cm H

5356-7351-000



Key BLUE Suitable for Size 1 GREEN Suitable for Size 2 YELLOW Suitable for Size 3

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Jenx Supine_20240821