

Jenx Therapy Bench

Order in only.
Not stocked in
Australia.



**THERAPY
BENCH**

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

October 2023

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Size Options

Measurements	Size 1 Full Length	Size 1 Half Length	Size 2 Full Length	Size 2 Half Length
Code	7510-0010-054	7510-0015-054	7510-0020-054	7510-0025-054
Max load weight	85kg		85kg	
Bench height from floor (without cushion)	165-250mm		235-395mm	
Bench top length	980mm	480mm	980mm	480mm
Bench top width	200mm		250mm	
Flat cushion width	250mm		300mm	
Flat cushion length	1000mm	500mm	1000mm	500mm
Flat cushion depth	25mm		25mm	


Package includes:

- Height Adjustable Bench
- 1" Comfortable Wipe Down Flat Cushion



Accessories

Half Round Cushion (half length only; order 2 for full length bench)

Size 1	Size 2	
500mm (L) x 250mm (W) x 125mm (H)	500mm (L) x 300mm (W) x 150mm (H)	
7510-0501-054	7510-0502-054	

Pelvic Support	Abduction Block
<p>Lateral size 2 support pads with flip away action. Cannot be used with half round cushion.</p> 	<p>Available in one size and used to provide simple abduction to ensure the right position for the child.</p> 
7510-1900-054	7510-1750-054

Accessories (cont.)

Knee Blocks - Size 1	Knee Blocks - Size 2	
Knee block with fully adjustable knee cups which can be used to control legs and create abduction		
7510-1601-054	7510-1602-054	

Additional Client Notes

Notes

medifab

Distributed by Medifab Ltd

22 Detroit Drive
Rolleston 7675
New Zealand

0800 543 343

solutions@medifab.com

medifab.com



DISTRIBUTED BY: