

BINGO Evo Stroller

Larger supportive seating system with many features, including a wide range of accessories and options to meet the needs of different children.

Modular seating & positioning stroller



Multi-directional suspension **protects the child** by absorbing ground bumps, and is effective off road or on rough terrain

Easily adjust to rear facing configuration for improved close interaction with caregiver or monitoring of safety

Stepless seat depth adjustment

Ergonomic, height-adjustable push handle

Sun canopy option available to allow parents and caregivers to still get out and about, knowing the child is protected from the sun

Adjustable seat unit
Seat depth and width, back support height and lower leg length

Note: Cover shown is an optional extra.

TUBIE FUN COVERS

Personalise and protect your stroller with waterproof and washable covers from Tubie Fun!



ONE SEAT. TWO BASE OPTIONS

Indoor Cobra base option available. Height adjustable while the child is seated.



TILT

The special construction of the tilt mechanism allows a 45° tilt of the seat unit without compromising the overall stability.



COMPACT FOLDED SIZE



Easier to store in a car boot! The innovative folding mechanism can be released easily on the push bar.

BACK SUPPORT RECLINE

The operating elements of the tilt & back support recline are positioned centrally and allow for easy and secure handling.



BINGO Switch:

At the heart of this innovative new design lies an interchangeable front end, a unique feature which eliminates the need to purchase multiple bases. Use the standard wheels for everyday journeys or quickly switch in the third wheel for off-road expeditions, delivering a safe, secure experience across all terrain.

Size 1: 2884-0001-500 **Size 2:** 2884-0002-500

Swivelling Front Wheel Attachment:

Size 1: 2884-9571-000 **Size 2:** 2884-9572-000



Stroller Packages Include:

Removable seat, matte black mobility base with swivel front black mag wheels and solid tyres, tie-downs, hip guides, basket, adjustable suspension and 3 sided foot support.



COLOUR OPTIONS:



Turquoise Soft
(-349)



Berry Soft
(-347)



Black
(-346)

Size	Code	Seat Load Capacity	Overall Length	Overall Width	Seat Depth	Seat Width	Back Support Recline	Shoulder Height	Foot Support Height	Overall Weight	Indoor Base Seat Height	Indoor Base Footprint
1	2881-0008-000	50kg	92cm	63cm	16-34cm	18-30cm	90-160°	28-46cm	16-37cm	15kg	27-67.5*	61cm (W) x 71cm (L)
2	2882-0008-000			68cm	22-40cm	23-35cm	90-160°	34-55cm	16-40cm	16.1kg		
3	2883-0008-000			68cm	27-45cm	23-35cm	90-135°	42-63cm	16-40cm	18kg		

* Depending on size of castors

BINGO Evo Accessories

Occipital Head Support
Merino/Spacer Fabric
2010-4110-099
Wipe Down
2010-4110-209



Std Lateral Head Support,
Merino/Spacer Fabric
Small: 2010-4210-099
Large: 2010-4211-099



Standard Lateral Head Support, Wipe Down
Small: 2010-4210-209
Large: 2010-4211-209



Ext Lateral Head Support,
Merino/Spacer Fabric
Small: 2010-4310-099
Large: 2010-4311-099



Extended Lateral Head Support, Wipe Down
Small: 2010-4310-209
Large: 2010-4311-209



Std Lateral Trunk Support
1: 2880-3010-346
2: 2880-3020-346
3: 2880-3030-346



Spex Lateral Trunk Supports
Narrow Straight: 2880-3008-100
Broad: 2880-3008-310
Teardrop: 2880-3008-500



Additional Padding
Sacral Support: 2020-2510-060
Thigh Wedge: 2006-0505-000



Medial Knee Support
1: 2880-1651-346
2: 2880-1652-346
3: 2880-1653-346



Anterior Chest Support
1: 2880-6751-346
2: 2880-6752-346
3: 2880-6753-346



Wonderseat Padded Centrepoint Harness
Size 0: 1409-5660-017
Size 1: 1409-6661-017



5-Point Anterior Chest Support
Bingo: 2880-6511-346
Spex: 2006-6551-000



Hip Belt
Spex 2-Point: 1405-6220-027
Wonderseat 4-Point: 1406-5420-027



Groin Strap
1: 2880-6361-346
2-3: 2880-6362-346



Foot Support Spacer 4cm thick
1: 2880-5321-000
2: 2880-5322-000



Ankle Foot Fast
1: 2880-6101-000
2: 2880-6102-000



Therapy Tray
1: 2880-7111-000
2-3: 2880-7112-000
Infill: 2880-7180-000



Grab Rail
1: 2880-1651-346
2-3: 2880-1652-346




IV Pole:
2880-8275-000
Extra Bracket
2880-8275-050



Oxygen Bottle Holder:
2880-8253-000



Chassis Mount Tray
1: 2880-8058-000
2: 2880-8059-000



Separate Seat Unit
1: 2881-0011-000
2: 2882-0012-000
3: 2883-0013-000



Sun Canopy & Rain Shield:
2880-7500-346



Cobra Indoor Base
1: 2880-8711-000
2: 2880-8721-000



Head Support Pads
1 2880-4471-346
2 2880-4472-346



Swivel Lock:
2880-8440-000



Cover Set - Llamas
Size 1: 2880-9561-121
Size 2: 2880-9562-121



5 Point Harness Cover:
2880-6511-121

Cover Set - Planets
Size 1: 2880-9561-036
Size 2: 2880-9562-036



5 Point Harness Cover:
2880-6511-036