

# Multiroller

Made from quality German-manufactured parts, the Multiroller is a lightweight, high-work-capacity stroller ideal for users with low-level postural support needs who require mobilisation assistance over long distances.


**Incredibly light, incredibly tough!**

**Build to Last!**  
Despite its light weight the Multiroller is ultra-durable, with a high work capacity and the ability to hold up to 75kg – perfect for growing users and older occupants.



**Lightweight and Compact**  
At just 9kg, this sleek and stylish stroller is particularly light, with an easy-fold function making it simple to transport.

**Easy Access**  
The high seat-to-floor ratio of the Multiroller makes transfers a breeze.



**Interchangeable Front End**  
Stylish, stable and easy to maneuver, the Multiroller includes an innovative feature which allows the front-end to be changed out for an alternative wheel set – simply swap in the optional jogger wheel set and you've instantly got an off-road solution ready for any terrain!





The high seat-to-floor ratio makes transfers much easier than conventional recreational strollers, while the sling upholstery is both breathable and comfortable. A canopy and drum brakes come as standard.

ACCESSORIES

---

**Running Wheel Add On**  
2610-8416-000



**5 Point Harness**  
2610-6500-000



---

**Seat Depth Extension**  
(Increases seat depth to 38cm)  
2610-0281-000



**Infill for Smaller Users**  
2610-0393-000



---

**2 Point Hip Belt – Centre Release – Small**  
1405-6211-028



**2 Point Hip Belt – Pen Key Release – Small**  
1405-6211-025



---

**Anti-Tip**  
2610-8420-000



Standard Package

- Stroller frame with easy to fold function
- Upholstered seat
- Height adjustable flip up footplate
- Interchangeable front end stroller frame wheelset
- Canopy
- Drum Brakes



**It's Got the Look!**  
Everything from the sleek black frame to the stylish sling upholstery has been designed to ensure the Multiroller looks as good as it performs.

Code	Overall Width	Seat Depth	Seat Width	Folded size*	Push Handle Height	Shoulder Strap Height	Lower Leg Length	Weight Limit	Approx Stroller Weight
2610-1000-000	68cm	33cm	30 - 38cm	34cm (H) x 56cm (W) x 84cm (L)	86 - 120cm	29cm, 38cm, 48cm	20 - 30cm	75kg	9kg

\*with rear wheels removed