

# Simple Stuff Works



## SLEEPING & POSITIONING SYSTEM

Distributed by  
**medifab**

---

SCRIPTING AND ORDERING GUIDE

# Client details

Client/User name:			
DOB:		Date:	
Therapist:			

# Notes

**Further Information**

---



---



---



---



---



---



---

CHOOSE A **PACKAGE** TO USE ON **FOAM MATTRESS**


# Packages

Infant Package	Child Package	Adult Package
<p><b>For infants from 0 to 5 years of age</b></p> <ul style="list-style-type: none"> <li>• 2x cot Bed Mesh</li> <li>• 2x cot Toppers</li> <li>• 2x cot cotton sheets</li> <li>• 2x small Lateral Supports</li> <li>• 2x small Wedges</li> <li>• Small Horseshoe pillow</li> <li>• Simple Stabiliser Size 1</li> <li>• Travel bag</li> </ul> <p>7150-7001-C00</p>	<p><b>For children from 5 to 12 years of age</b></p> <ul style="list-style-type: none"> <li>• 2x single Bed Mesh</li> <li>• 2x single Toppers</li> <li>• 2x single cotton sheets</li> <li>• 2x medium Lateral Supports</li> <li>• 2x small Wedges</li> <li>• 2x small Lateral Supports</li> <li>• Simple Stabiliser Size 2</li> <li>• Small Terry Towelling Pillow</li> <li>• Travel bag</li> </ul> <p>7150-7002-S00</p>	<p><b>For adults from 12+ years of age</b></p> <ul style="list-style-type: none"> <li>• 2x single Bed Mesh</li> <li>• 2x single Toppers</li> <li>• 2x single cotton sheets</li> <li>• 2x medium Lateral Supports</li> <li>• 2x medium Wedges</li> <li>• 2x large Lateral Supports</li> <li>• Large Terry Towelling Pillow</li> <li>• Travel bag</li> </ul> <p>7150-7003-S00</p>




OR REFER TO **PAGE 3** TO **CONFIGURE** YOUR OWN

# Sheet Sets

Cot Bed	Single Bed	Single Bed
<p>Set Includes:</p> <ul style="list-style-type: none"> <li>Stabilising mesh set (includes 2 mesh sheets; 1 bottom, 1 top)</li> <li>2 x cot temperate toppers</li> <li>2 x cot fitted sheets</li> </ul> <p>7150-7021-C00</p>	<p>Set Includes:</p> <ul style="list-style-type: none"> <li>Stabilising mesh set (includes 2 mesh sheets; 1 bottom, 1 top)</li> <li>2 x single temperate toppers</li> <li>2 x single fitted sheets</li> </ul> <p>7150-7022-S00</p>	

# Sheets


Stabilising Mesh (Set of 2)		
<p><b>Cot Bed</b> Base layer: 1500 mm (L) x 700 mm (W) Top layer: 1800 mm (L) x 1300 mm (W) 7150-7011-C00</p>	<p><b>Single Bed</b> Base layer: 1800 mm (L) x 1000 mm (W) Top layer: 1800 mm (L) x 1500 mm (W) 7150-7012-S00</p>	

Temperate Topper (Set of 2)		
<p><b>Cot Bed</b> 1000 mm (L) x 590 mm (W) 7150-7051-C00</p>	<p><b>Single Bed</b> 1000 mm (L) x 1000 mm (W) 7150-7052-S00</p>	

Fitted Sheets			
<p><b>Cot Bed</b> 1400 mm (L) x 700 mm (W) 7150-7121-C00</p>	<p><b>Single Bed/King Single</b> 1900 mm (L) x 900 mm (W) 7150-7122-S00</p>	<p><b>Double/Queen Bed</b> 2030 mm (L) x 1530 mm (W) 7150-7123-S00</p>	

Flat Sheet
<p><b>Single</b> 1620 mm (L) x 2600 mm (W) 7150-7126-S00</p>

# Wedges


Wedges (pair)		
With padded covers.		
<p><b>Small</b> A: 220 mm B: 110 mm C: 110 mm 7150-7351-000</p>	<p><b>Medium</b> A: 360 mm B: 80 mm C: 90 mm 7150-7352-000</p>	

# Lateral Supports

## Padded Lateral Supports




Also available in an option for air mattresses, designed to grip directly to the mattress surface.

Small (Foam)	Medium (Foam)	Large (Foam)
7150-7201-000	7150-7202-000	7150-7203-000
Small (Air)	Medium (Air)	Large (Air)
7150-7211-000	7150-7212-000	7150-7213-000
A: 180 mm B: 120 mm C: 180 mm	A: 220 mm B: 140 mm C: 220 mm	A: 280 mm B: 220 mm C: 280 mm



# Simple Stabiliser

## Supports hips and thighs in neutral supine alignment

Size 1	Size 2	Size 3	Size 4
H: 110mm D: 110mm W: 420mm	H: 190mm D: 350mm W: 460mm	H: 250mm D: 400mm W: 650mm	H: 280mm D: 400mm W: 650mm
To fit approx leg circumference of 250mm. Suggested age range 2 - 4yrs.	To fit approx leg circumference of 350mm. Suggested age range 2 - 10yrs.	To fit approx leg circumference of 550mm. Suggested age range 5 - 13yrs.	To fit approx leg circumference of 600mm. Suggested age range 13 - adult.
<b>With breathable inner &amp; terry cotton outer</b> 7150-7341-000	<b>With breathable inner &amp; terry cotton outer</b> 7150-7342-000	<b>With breathable inner &amp; terry cotton outer</b> 7150-7343-000	<b>With breathable inner &amp; terry cotton outer</b> 7150-7344-000
<b>With wipeable incontinence inner &amp; terry cotton outer</b> 7150-7345-000	<b>With wipeable incontinence inner &amp; terry cotton outer</b> 7150-7346-000	<b>With wipeable incontinence inner &amp; terry cotton outer</b> 7150-7347-000	<b>With wipeable incontinence inner &amp; terry cotton outer</b> 7150-7348-000
			
<b>Spare terry cotton cover</b> 7150-7221-000	<b>Spare terry cotton cover</b> 7150-7222-000	<b>Spare terry cotton cover</b> 7150-7223-000	<b>Spare terry cotton cover</b> 7150-7224-000
<b>Spare wipeable inner cover</b> 7150-7225-000	<b>Spare wipeable inner cover</b> 7150-7226-000	<b>Spare wipeable inner cover</b> 7150-7227-000	<b>Spare wipeable inner cover</b> 7150-7228-000

# Side Lying Leg Support

Assists with pelvic and thigh positioning in side-lying

Size 1	Size 2	Size 3
<b>A:</b> 270mm <b>B:</b> 380mm <b>C:</b> 280mm <b>D:</b> 80mm <b>E:</b> 130mm <b>F:</b> 80mm	<b>A:</b> 400mm <b>B:</b> 560mm <b>C:</b> 340mm <b>D:</b> 100mm <b>E:</b> 200mm <b>F:</b> 90mm	<b>A:</b> 400mm <b>B:</b> 610mm <b>C:</b> 400mm <b>D:</b> 110mm <b>E:</b> 150mm <b>F:</b> 70mm
With breathable inner & terry cotton outer (standard) 7150-7971-000	With breathable inner & terry cotton outer (standard) 7150-7972-000	With breathable inner & terry cotton outer (standard) 7150-7973-000
With wipeable incontinence inner & terry cotton outer 7150-7981-000	With wipeable incontinence inner & terry cotton outer 7150-7982-000	With wipeable incontinence inner & terry cotton outer 7150-7983-000
Spare terry cotton cover 7150-7991-000	Spare terry cotton cover 7150-7992-000	Spare terry cotton cover 7150-7993-000

# Foot Support


Provides positioning & skin protection for the feet (pair)

<b>Size 1</b> <b>A:</b> 210mm (max) <b>B:</b> 200mm 7150-5351-000	<b>Size 2</b> <b>A:</b> 250mm (max) <b>B:</b> 220mm 7150-5352-000	<b>Size 3</b> <b>A:</b> 300mm (max) <b>B:</b> 300mm 7150-5353-000	
--	--	--	--

# Pillows


## Roll Pillow

<b>Small</b> A: 150mm B: 400mm	<b>Large</b> A: 180mm B: 500mm
<b>With breathable inner &amp; terry cotton outer</b> 7150-7851-000	<b>With breathable inner &amp; terry cotton outer</b> 7150-7852-000
<b>With wipeable incontinence inner &amp; terry cotton outer</b> 7150-7855-000	<b>With wipeable incontinence inner &amp; terry cotton outer</b> 7150-7856-000
<b>Spare terry cotton cover</b> 7150-7857-000	<b>Spare terry cotton cover</b> 7150-7858-000




## Neck Support Pillow

(Suitable for air & foam mattresses). 300mm x 350mm x 100mm 7150-4050-000
<b>Spare terry cotton cover</b> 7150-4055-000



## Horseshoe Pillow

<b>Small</b> A: 460mm B: 190mm	<b>Large</b> A: 500mm B: 200mm
<b>With breathable inner &amp; terry cotton outer</b> 7150-7711-000	<b>With breathable inner &amp; terry cotton outer</b> 7150-7712-000
<b>With wipeable incontinence inner &amp; terry cotton outer</b> 7150-7721-000	<b>With wipeable incontinence inner &amp; terry cotton outer</b> 7150-7722-000
<b>Spare terry cotton cover</b> 7150-7731-000	<b>Spare terry cotton cover</b> 7150-7732-000



## Terry Towling Pillow

<b>Small</b> L: 480mm W: 350mm 7150-7501-000	<b>Large</b> L: 680mm W: 450mm 7150-7502-000
<b>Spare terry cotton cover</b> 7150-7511-000	<b>Spare terry cotton cover</b> 7150-7512-000



# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: SimpleStuffWorks\_20241216