

medifab

Shaping Better Lives



Juni Chair
Indoor Seating

Confidence,
engagement,
learning



Whether in the classroom, nursery or at home, better posture while seated creates a more positive experience for your child. It can help develop fine motor skills, increase confidence and the ability to engage and learn.

The Juni chair comes in four colours and the contemporary design helps to promote social inclusion.

Juni's seat and backrest provide unique contouring to promote better posture, give higher levels of comfort and help your child feel more secure and stable in the chair.


It is our vibrant chair for children who require minimal levels of postural support when seated.

We understand that your child's needs change during each stage of growth and development. That's why the seat's height and depth, as well as armrest height, can easily be adjusted.



Base Package includes seat unit and hip belt.

Features



Colour Options

Juni comes in four colours (Red, Blue, Lime and Black) and the contemporary design helps to promote social inclusion.

Tailor Juni to your child's needs

With exceptional adjustment of Juni, along with a great range of optional accessories, your child's Juni can be tailored to them, supporting them to achieve their goals and more!





Contoured seat base and back support as standard

Gentle contouring provides support and gives stability for sitting, promoting engagement and focus on the task at hand.



Low seat to floor height

Juni's low seat to floor height allows children to be sat with their feet on the floor, increasing support and allowing for more concentration on their play and learning.



High weight rating

Juni has a max user weight of 30-60kg. Its sturdy frame helps to weight and prevent chair movements.

Accessories



Side Cushions

Reduce seat width using side supports. The thickness of the cushion equates to the amount reduced.

For size 1 and 2 Juni:

25mm 3209-0001-025

50mm 3209-0001-050

75mm 3209-0011-050

For size 3 and 4 Juni:

25mm 3209-0003-025

50mm 3209-0003-050

75mm 3209-0013-050



Standard Skis

Helps to stabilise and prevent chair movements.

Size 1 3209-5021-000

Size 2 3209-5022-000

Size 3 3209-5023-000

Size 4 3209-5024-000



Mobile Skis

Makes the Juni chair easier to move, ideal for pushing under tables!

Size 1 3209-5051-000

Size 2 3209-5052-000

Size 3 3209-5053-000

Size 4 3209-5054-000



Sliding Footboard

Add footboard to mobile skis to support users legs

Size 1 3209-5071-000

Size 2 3209-5072-000

Size 3 3209-5073-000

Size 4 3209-5074-000



Adjustable and Foldable Table

Adjustable height work surface, can be easily folded for storage, a great companion to the Juni chair.


Size 1 3209-7191-000

Table height 31-45cm

Size 2 3209-7192-000

Table height 40-60cm

Specifications



Code	Size 1	Size 2	Size 3	Size 4
Black	3209-1001-099	3209-1002-099	3209-1003-099	3209-1004-099
Blue	3209-1001-060	3209-1002-060	3209-1003-060	3209-1004-060
Red	3209-1001-022	3209-1002-022	3209-1003-022	3209-1004-022
Lime	3209-1001-054	3209-1002-054	3209-1003-054	3209-1004-054
Seat Height (cm)	16-22	20-28	24-34	30-40
Seat Width (cm)	24	28	32.5	36.5
Seat Depth (cm)	18-24	22-30	26-36	32-42
Back Support Height (cm)	26	29	32.5	40
Load capacity of seat unit (kg)	30	40	50	60
Product weight (kg)	4	6	8	10

medifab

Distributed by Medifab Ltd

22 Detroit Drive
Rolleston 7675
New Zealand

26 Pardoe Street
Devonport TAS 7310
Australia

0800 543 343

1300 543 343

solutions@medifab.com

medifab.com

DISTRIBUTED BY:

Flyer: Juni Chair_20230914