

Zoomi+



INDOOR SEATING

SCRIPTING AND ORDERING GUIDE

Distributed by
medifab

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Order Code	Size 1
Black	3212-1001-099
Blue	3212-1001-060
Red	3212-1001-022
Lime	3212-1001-054

Seat Unit Measurements	
Seat height (cm)	81
Seat width (cm)	43
Seat depth (cm)	15-26 (increased adjustment possible with accessories)
Foot support height (cm)	26-60
Back support height (cm)	27
Load capacity of seat unit (kg)	85
Product weight (kg)	16

1 Ranges given are recommendations only - further adjustment range is possible.
 2 Foot board can be removed, leaving option of foot flat to floor.

Base package includes:

- Seat Unit
- Contoured Padded Seat
- High Back Support
- Arm Supports
- 5 Point Harness
- Foot Support
- Quick Release Tray
- Mobile Ski's



Accessories

Transfer Step

Add on accessory to facilitate a standing transfer.

Black
3212-5901-099

Blue
3212-5901-060

Red
3212-5901-022

Lime
3212-5901-054



Extended Slotted Footboard

Black
3212-5021-099

Blue
3212-5021-060

Red
3212-5021-022

Lime
3212-5021-054



Extended Contoured Padded Seat

Upgraded seat depth for bigger users, addition 10cm seat depth.

Black
3212-8991-099

Blue
3212-8991-060

Red
3212-8991-022

Lime
3212-8991-054



Wonderseat Padded Centrepoint Shoulder Harness

Size 0
Pad height 22cm
1409-5660-017

Size 1
Pad height 28cm
1409-5661-017



Spex 2-Point Hip Belt

X-small
2 point slim line, dual pull
with neoprene sleeve padding.
1405-6220-027



medifab

Distributed by Medifab Ltd

22 Detroit Drive
Rolleston 7675
New Zealand

0800 543 343

solutions@medifab.com

medifab.com



DISTRIBUTED BY: