

# Carrot 3000




## CAR SEAT

### AUSTRALIAN SCRIPTING AND ORDERING GUIDE

Distributed by  
**medifab**

# Client Details



Client/User Name:	
Date:	
NDIS/Funding:	
<b>User Measurements</b>	
A Hip width	
B Seat depth	
C Seat to top of shoulder	
D Seat to top of head	
E Chest width	
F Seat to axilla	
G Leg length	







# Seat Specifications

Age Range	User Weight Range	Seat Depth	Shoulder Height	Seat Weight	Testing Data
3-15 years	13.6-49kg	25-40 cm	32-47 cm 45-60 cm (when shoulder protector wings are added)	8.7kg	Tested to ECE R44/04 

# Basic Packages

Includes:	Small	Large
<ul style="list-style-type: none"> <li>Carrot seat unit</li> <li>5-point harness</li> <li>Headrest wings</li> <li>Tee-pad</li> <li>Free recline backrest</li> </ul>	<p>4116-0010-259</p> 	<p>Larger chest pad and deep headrest wings 4116-0000-259</p> 

# Behavioural Packages

Includes:	Small	Large - larger tee pad and deep headrest wings
<ul style="list-style-type: none"> <li>Carrot seat unit</li> <li>5-point harness</li> <li>Headrest wings</li> <li>Tee-pad</li> <li>Free recline backrest</li> </ul>	<p>Small package with anti-escape buckle. 4116-0010-026</p> 	<p>Large package with anti-escape buckle. 4116-0000-026</p> 
	<p>Small package with pen-key buckle. 4116-0010-025</p> 	<p>Large package with pen-key buckle. 4116-0000-025</p> 

### Swivel Base Attachment

#### ISOfix essential

(Standard adjustment of mandatory footprop is 34cm-40cm, if this is too short please contact Medifab Customer Service for Quote)

4116-0180-104



### Swivel Base Protection Plate

Foam plate to protect car seat cover, mounts under swivel base  
4116-0190-000




## Seat Depth Growth

For growth, please select one of the following 3 rodded options, and associated extensions below, if required:  
(Rodded components cannot be used in combination with each other).

### Option 1

#### Rodded 50mm Seat Depth Extension

4116-0281-259



#### 50mm Growth Pad

4116-0285-259



### Option 2

#### Rodded Waterfall Seat Depth Extension

4116-0283-259




#### 50mm Growth Pad\*

4116-0285-259



#### 100mm Growth Pad\*

4116-0289-259



### Option 3

#### Rodded Foot Support

##### Small

Lower leg (D): 25-30 cm

4116-5241-250

##### Medium

Lower leg (D): 30-40 cm

4116-5242-250

##### Large

Lower leg (D): 40-50 cm

4116-5243-250




#### 50mm Growth Pad\*

4116-0285-259



#### 100mm Growth Pad\*

4116-0289-259



#### Adaptor

For use of foot support with swivel base  
4116-5120-000



\*50mm and 100mm growth pads cannot be used in combination with each other.

# Back Support Height Growth






50mm Height Extension	100mm Height Extension	Shoulder Protector Wings
<p>50mm back support height extension infill pad. 4116-2481-259</p>	<p>50mm back support height extension infill pad. 4116-2482-259</p>	<p>(For children over 45cm shoulder height). 4116-2753-259</p>

\*50mm and 100mm infill pads cannot be used in combination with each other

# Installation Options




Tether Strap	ISOfix lower straps	Airplane Fixing Strap
<p>4116-6953-017</p>	<p>No installation option removes the need for the vehicle safety belt to be used, this must be used in every situation. 4116-6910-000</p>	<p>For securing to aircraft seat. 4116-6960-017</p>

# Cover Options

Wipe-down		
<p>Wipe-down seat base. 4116-0476-209</p>		<p>Tee-pad, Standard. 4116-6553-209</p>
<p>Wipe-down backrest. (Main unit only). 4116-2476-209</p>		<p>Tee-pad, Large. 4116-6558-209</p>
		
Removable Incontinence Cover		Air Mesh Cushion Pad
<p>Attaches using Velcro. Unable to be attached to wipe down covers. 4116-0475-209</p>		<p>Attaches using Velcro. Unable to be attached to wipe down covers. 4116-8880-099</p>
		

# Accessories (continued)



Tray / Forearm Support		Tilt Wedge		Upgraded Lateral Supports	
------------------------	--	------------	--	---------------------------	--

4116-5534-259		Provides 20° tilt. 4116-0291-259		Firm padding for added trunk control. 4116-3250-090	
---------------	---	-------------------------------------	--	--	---

Seat Belt Holder		Shoulder Strap Retainer	
------------------	--	-------------------------	--

Designed to keep seat belt away from the neck. 4116-6790-017		<b>Small</b> 4116-6721-000	
		<b>Medium</b> 4116-6722-000	
		<b>Large</b> 4116-6723-000	

Medial Thigh Support		Hinge Locking Belt	
----------------------	--	--------------------	--

Detachable & Washable. Attaches to chest support and crotch strap. 4116-1510-000		For easy carrying. 4116-8170-000	
---	--	-------------------------------------	---

Hip Belt		Accessory Cross Strap	
----------	--	-----------------------	--

<b>Small</b> 4116-6531-027		<b>Anti-Escape</b> 4604-3004-026	
<b>Medium</b> 4116-6532-027		<b>Pen Release</b> 4604-3004-025	
<b>Large</b> 4116-6533-027		<b>Side Release</b> 4604-3112-019	



# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: Carrot\_20231215\_AU