

# Duro Stroller

Whether navigating city footpaths or the outdoors terrain, Duro is a stable stroller solution with a great range of support and positioning accessories that will enable you to embrace the outdoors.

## Configurable rehab stroller for everyday challenges



Extremely agile with **great manoeuvring** due to compact mobility base

Adjustable back angle

2 options for seat angle

Adjustable knee angle

Generous **growth capacity** provides long-term use for growing users

Individually customisable to suit requirements, with proven accessories from the **BINGO Evo Stroller**

Very stable stroller option



Super compact folding makes the Duro easy to transport and stow in the boot of your car!



Reclining back support means you can relax back in the Duro when being out and about gets tiring!

## Colour options: *(Upholstery must be selected in addition to Stroller)*

**Turquoise**  
Size 1:  
2780-8801-349  
Size 2:  
2780-8802-349




**Black**  
Size 1:  
2780-8801-346  
Size 2:  
2780-8802-346


## ACCESSORIES

Standard package includes: hip supports, swivel front wheels, flip-up foot support, soft accessory bag, tie-down kit and retainer for accessories.

**Head Support Pads**  
1: 2880-4471-346  
2: 2880-4472-346



**Therapy Tray**  
1: 2880-7111-000  
2-3: 2880-7112-000  
Infill: 2880-7180-000



**Occipital Head Support Spacer:** 2010-4110-099  
**Wipe Down:** 2010-4110-209



**Std Lateral Head Support Spacer:** 2010-4210-099  
**Wipe Down:** 2010-4210-209



**Std Lateral Trunk Support**  
1: 2880-3010-346  
2: 2880-3020-346  
3: 2880-3030-346




**Spex Broad Lateral Supports:** 2880-3008-310  
**Spex Teardrop Lateral Supports:** 2880-3008-500



**Spex Centrepoint Harness**  
XS: 1409-6660-017  
S: 3780-4301-009



**5-Point Harness:** 2880-6511-346



**Canopy**  
Size 1: 2780-7501-346  
Size 2: 2780-7502-346



Size	Code	Weight Limit	Seat Width	Seat Depth	Seat Angle	Back Support Height	Back Support Angle	Lower Leg Length	Wheel Diameter	Folded Size (L x W x H)
1	2780-0001-052	75 kg	23 - 34 cm	26 - 40 cm	16° or 20°	45 - 63 cm	78 - 120°	20 - 44 cm	20 cm front, 30 cm rear	84x63x32 cm
2	2780-0002-052	75 kg	28 - 39 cm	35 - 45 cm	16° or 20°	62 - 83 cm	78 - 120°	20 - 42 cm	20 cm front, 30 cm rear	84x68x38 cm