


# SPEX CONSTRUCTA SHAPE CUSHION

## SPEX CONSTRUCTA SHAPE CUSHION PARTS (EXPLODED)

Use this checklist to identify the parts you have or the parts which need replacing.



### COVER

- A** OUTER COVER
- B** INNER COVER
- C** OVERLAY WITH COVER

### CONTOURS

- D** ISCHIAL INSERT
- E** LATERAL PELVIC CONTOURS
- F** LATERAL THIGH CONTOURS
- G** MEDIAL CONTOUR

### BASE

- H** POSITIONING BASE
- POSITIONING BLOCKS
- I** TROCHANTER BLOCKS
- J** 5CM THIGH BLOCK
- K** THIGH BLOCK

### ACCESSORIES


- G** STRATEGIC POSITIONING KIT

Note: Pictures are for reference and are not shown to scale. Refer to the Spex Spare Parts Script Form for more information.

# SPEX CONSTRUCTA FLEX CUSHION

## SPEX CONSTRUCTA FLEX CUSHION PARTS (EXPLODED)

Use this checklist to identify the parts you have or the parts which need replacing.



### COVER



- A** OUTER COVER
- B** INNER COVER
- C** OVERLAY WITH COVER

### CONTOURS

- D** ISCHIAL INSERT
- E** LATERAL PELVIC CONTOURS
- F** LATERAL THIGH CONTOURS
- G** MEDIAL CONTOURS

### BASE

- H** FLEX LEG WITH POSITIONING BLOCKS
- POSITIONING BLOCKS
- I** THIGH BLOCK
- J** 5CM THIGH BLOCK
- K** LEG WEDGE
- L** TROCHANTER BLOCKS
- M** REAR UNDERPANEL PIVOT
- N** PELVIC ADJUSTMENT KIT

Note: Pictures are for reference and are not shown to scale.  
Refer to the Spex Spare Parts Script Form for more information.