

FLUX Walker



VERSATILE POSTERIOR WALKING AID

AUSTRALIAN
SCRIPTING AND ORDERING GUIDE

Distributed by
medifab

Client Details

Client/User Name		Notes
Date		
NDIS/Funding		

Measurements

Product specifications	Size 1	Size 2	Size 3
Standard grab rail grip height	48-61cm	57-81cm	69-87cm
Standard grab rail grip width	41cm	46cm	51cm
Universal grip grab rail grip height	43-68cm	51-88cm	62-94cm
Universal grip grab rail grip width	29-41cm	34-46cm	39-51cm
Total frame width	64cm	69cm	74cm
Total frame length	65cm	80cm	88cm
Total load capacity	60kg	75kg	75kg
Front / rear wheel diameter	15/20cm	15/20cm	15/20cm

Package & Size Options

Package contents	Size 1	Size 2	Size 3
	5641-0001-000	5642-0001-000	5643-0001-000
<ul style="list-style-type: none"> • Standard Handlebar • Front Swivel Wheels • Polyurethan Tyres • Bumper Wheels • Friction Brakes • Anti-Reverse Mechanism 		 (Shown with universal handlebar)	 (Shown with universal handlebar)
	Luminous Red	Sky Blue	Gentian Blue

Accessories

Anti-tippers - without castors

5640-8421-000



Anti-tippers with castors

5640-8425-000



Fold Up Seat

For size 1-2 walker

5640-0201-000

For size 3 walker

5640-0204-000



Hand Disk

Shown with Universal Handlebar upgrade

5640-5550-000



Swivel Lock

To hold castors straight when brisk walking

2880-8440-000



Driving Brake

Includes Wheel Lock

5640-8451-000



Hip Pad


Multi-Adjustable



5640-1900-000



Accessories Continued

Dynamic Hip Control (Hip Belt Required)		Hip Belt (Required)	
Suspension, soft 5640-1952-000		Small Hip Size: 53-83cm 5640-6201-000	
Suspension, hard 5640-1953-000		Large Hip Size: 65-115cm 5640-6202-000	
Kit for orthoses 5640-1951-000			

Universal Handlebar	
Size 1 5641-5751-000	
Size 2 5642-5752-000	
Size 3 5643-5753-000	

Forearm Supports – 2 Axial Adjust			
Small 14 x 9.5cm 5640-5500-000		Standard 22 x 10.5cm 5640-5501-000	

Notes

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Flux_20231215_AU