

Gravity Chair




Distributed by
medifab

INDOOR SEATING

SCRIPTING AND ORDERING GUIDE

Our product codes are meaningful!


Follow the 3 steps below, to build your gravity chair product code:



3 7 7 0 - 0 0 X X - 0 X X

EXAMPLE: A large Gravity Chair with a green classic liner, would have the code 3770-0013-055.

Step 1: Select Liner

Classic Version	Wipe Down Liner
<p>Air breathable fabric on the user contact surfaces and rear back panel, with wipe-down covers in all other areas -001X-</p> 	<p>Fully wipe-down, including internal contact surfaces where user is seated -002X-</p> 


Step 2: Select Size

Measurements	Small	Medium	Large	X-Large	XL-Wide	Notes
Order Code	-00X1-	-00X2-	-00X3-	-00X4-	-00X5-	
A Seat Depth	24cm	30cm	38cm	43cm	43cm	
B Seat Width at Back	21cm	23cm	27cm	30cm	41cm	
C Seat Width at Front	28cm	30cm	37cm	39cm	50cm	
D Back Height	60cm	69cm	74cm	81cm	81cm	
E Seat Front Height	20cm	23cm	26cm	29cm	29cm	
F Arm Support Height	29cm	35cm	39cm	43cm	43cm	
G Total Width of Chair	46cm	49cm	57cm	60cm	71cm	
H Total Height of Chair	55cm	61cm	67cm	71cm	72cm	
I Total Depth of Chair	60cm	71cm	82cm	91cm	91cm	
J Tilt Wedge Length	47cm	57cm	67cm	74cm	74cm	
K Tilt Wedge Height	10cm	12cm	14cm	17cm	17cm	
Weight Limit	100kg					





Gravity Chair Package includes:

- Foam chair with cover
- Coloured liner
- Standard head support
- Tilt wedge
- Velcro side pocket

Fully wipe-down, including internal contact surfaces where user is seated.




Step 3: Select Liner Colour

Red Liner	Green Liner	Blue Liner	Black Liner
-022 	-055 	-066 	-099 


Head Supports

Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
Merino/Spacer Fabric 2010-4110-099		Small (17cm internal width)		Small (17cm internal width)	
		Merino/Spacer Fabric 2010-4210-099		Merino/Spacer Fabric 2010-4310-099	
		Wipe Down Cover 2010-4210-209		Wipe Down Cover 2010-4310-209	
Wipe Down Cover 2010-4110-209		Large (21cm internal width)		Large (21cm internal width)	
		Merino/Spacer Fabric 2010-4211-099		Merino/Spacer Fabric 2010-4311-099	
		Wipe Down Cover 2010-4211-209		Wipe Down Cover 2010-4311-209	

Mobile Base

For ease of maneuverability in the home or classroom				
Small 3780-8701-000	Medium 3780-8702-000	Large 3780-8703-000		
XL 3780-8704-000	XL-Wide 3780-8705-000			



Tray Table

Encourages hand-to-eye coordination			
To fit Size Small and Medium 3770-7201-000	To fit size Large and XL 3770-7202-000	To fit size XL-Wide 3770-7203-000	

Foot Support


Small		Large	
(Suitable for Small or Medium Gravity chairs) 3770-1031-000		(Suitable for Large, XL and XLW Gravity Chairs) 3770-1032-000	

Abductor Pommel

Small		Large	
<p>(Suitable for Small or Medium Gravity chairs)</p> <p>3770-1021-000</p>		<p>(Suitable for Large, XL and XLW Gravity Chairs)</p> <p>3770-1022-000</p>	


Side Filler Pads

For reducing width of Gravity Chair by up to 6cm, adjustable in 1cm increments

Small 3770-1011-000	Medium 3770-1012-000	Large 3770-1013-000	XL/XLW 3770-1014-000	
-------------------------------	--------------------------------	-------------------------------	--------------------------------	-------------------------------------------------------------------------------------

Depth Adjustment Liner

For reducing the depth of the Gravity Chair

Small 3770-0321-000	Medium 3770-0322-000	Large 3770-0323-000	
XL 3770-0324-000	XLW 3770-0325-000		

Notes

This image shows a single sheet of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Gravity_20231220