



# Lying and Positioning



A 24-hour postural management program is essential to ensuring a user's posture is maintained at its optimum level.

Because everyone is different and has different needs for positioning, Medifab's lying and positioning range includes from basic pillows to highly customised sleeping systems.

We offer a large array of components that can even be customised for any user needs and body shapes.

---

Dreama	4
Simple Stuff Works	7
Care Wave	13
Kid-In-Form	19
Poz-In-Form	23

---

# Resources to help you configure sleep systems to accurately meet your needs and those of the client!




Explore our extensive online resources available for downloading at [www.medifab.com.au/knowledge-base](http://www.medifab.com.au/knowledge-base)

# Dreama Positioning

- Combined pressure reducing mattress and postural support system, suitable for children and adults.
- Supporting at risk limbs and reducing tonal spasms, Dreama contributes towards more recuperative periods of rest, which in turn has a positive contribution towards minimising the risk of muscular skeletal deformity.
- Dreama can be used as a preventative system or to reduce the risk of further postural deterioration, and is appropriate for any age.
- Users who may benefit from Dreama include those with the following conditions: Cerebral Palsy, Muscular Dystrophy, Spinal Muscular Atrophy, Arthrogryposis, Rheumatoid Arthritis and associated conditions, Post-trauma, Post-surgery, Spina Bifida and Scoliosis.




- Contour the flexible aluminium base for a custom lying surface. Suits profiling beds, rolls, wedges.



- Quiet with no noise or disturbance when adjusting positions due to Glide-lock supports.





- Quick, simple assembly, and exceptional flexibility with numerous accessories.



# Dreama Base & Mattress


- Includes flexible aluminium base, full set of mattress cells, full set of neutral covers, and pair of anchor straps.
- Pressure-reducing 15 cm deep modular mattress removes the need for overlays.
- Individual covers mean any that are soiled can be removed and replaced without the need to remove the user from the bed, minimising sleep disruption and the need for assisted transfers.
- Spaces between mattress cushions allow airflow, helping to maintain a comfortable temperature. Covers are upholstered in fire-retardant, water resistant, 2-way stretch, antibacterial fabric for infection control.
- Add locking feet (#7120-6100-000) to mark the position where an accessory should be placed. Locking feet keep accessories in place for users with high movement.



	Cot	Single bed
<b>Code</b>	7120-0010-000	7120-0020-000
<b>Age range</b>	9mths - 4yrs	9mths - adult
<b>Overall size</b>	56 x 120cm	90 x 190cm
<b>Overall height</b>	17cm	17cm
<b>Single mattress cushion</b>	56 x 9.5 x 15 cm	120 x 9.5 x 15 cm
<b>Double mattress cushion</b>	56 x 19.5 x 15 cm	120 x 19.5 x 15 cm
<b>Number of mattress units</b>	8 single, 2 double	15 single, 2 double
<b>Max user weight</b>	80kg	110kg

# Dreama Positioning Pads


- Made from comfortable and hygienic PU. Positioning Pad Combos include angle-adjustable positioning pad, memory foam pad, and towelling cover.
- Use to position the trunk, arms or legs.
- Slimline Support option is an alternative slimmer support than the standard positioning pads.
- Removable memory foam pad (shown in green) for extra comfort and pressure relief. Ideal for support of areas with bony prominences.
- Additional covers available for all positioning pads in white towelling.



	Pad Combo	Spare Cover	L	W	H
<b>Slimline</b>	7120-3250-500	7120-3280-000	-	108 mm	159 mm
<b>Small</b>	7120-3251-500	7120-3281-000	102 mm	83 mm	95 mm
<b>Medium</b>	7120-3252-500	7120-3282-000	102 mm	83 mm	140 mm
<b>Large</b>	7120-3253-500	7120-3283-000	203 mm	89 mm	165 mm
<b>X large</b>	7120-3254-500	7120-3284-000	203 mm	89 mm	267mm

# Dreama Abduction Block

- Can be used individually or in pairs to ensure the hips are abducted and the body is straight.
- Ideal for users with increased abductor tone to help with hip preservation.



	Size 1	Size 2
<b>Code</b>	7120-1601-000	7120-1602-000
<b>Cover</b>	7120-1661-000	7120-1662-000
<b>Width</b>	130 mm	240 mm
<b>Depth</b>	200 mm	300 mm
<b>Height</b>	260 mm	310 mm

# Dreama Knee Blocks

- Comfortably holds the knees in a specified position.


	Small	Large
<b>Code</b>	7120-1751-000	7120-1753-000
<b>Cover</b>	7120-1781-000	7120-1783-000
<b>Width</b>	77 mm	108 mm
<b>Depth</b>	102 mm	102 mm
<b>Height</b>	51 mm	102 mm



# Dreama Side Lying System

- Complete system includes combinations of positioning pads, straps and Side Lying (abduction) foam pad. The straps attach to standard Dreama positioning pads.
- Side Lying Support Strap can be used to secure the upper leg when side lying. Includes fittings to fit to two standard positioning pads.
- Side Lying foam pad is available in different thicknesses to encourage a neutral leg separation when side lying.


	Size 1	Size 2	Size 3
<b>Lying Pad</b>	7120-1851-000	7120-1852-000	7120-1853-000
<b>Pad Cover</b>	7120-1861-000	7120-1862-000	7120-1863-000
<b>Pad Width</b>	345 mm	445 mm	540 mm
<b>Pad Depth</b>	120 mm	140 mm	159 mm
<b>Pad Height</b>	35 mm	51 mm	76 mm
<b>Strap Code</b>	7120-1872-000		
<b>Strap Cover</b>	7120-1882-000		
<b>Strap Dims</b>	140 mm wide x 388 mm long x 57 mm high		



# Dreama Positioning Hoop

- Functional item which keeps bed clothes away from the body.
- Allows movement within a limited range for a comfortable night's sleep.


	Small	Medium	Large
<b>Code</b>	7120-1951-000	7120-1952-000	7120-1953-000
<b>Cover</b>	7120-1961-000	7120-1962-000	7120-1963-000
<b>Width</b>	70 mm	70 mm	70 mm
<b>Length</b>	432 mm	560 mm	700 mm
<b>Height</b>	15 mm	15 mm	15 mm




# Dreama T-Roll

- Versatile, height adjustable accessory. Supports the knees, reduces abdominal pressure and prevents scissoring and wind-sweeping.
- Sides are height adjustable and the central pommel can be removed.

	One Size
<b>Code</b>	7120-1771-000
<b>Cover</b>	7120-1761-000
<b>Width</b>	527 mm
<b>Length</b>	108 mm
<b>Height</b>	121 mm





From our family to yours

The Simple Stuff Works way provides high quality education and therapeutic positioning equipment with an aim to eradicate the preventable problems associated with body shape distortion.

Simple Stuff Works advisors and colleagues are all highly experienced, having written some of the leading research into Postural Care and its effects. Medifab educators have been upskilled by these individuals and are ambassadors to promote protection of body shape.



Therapeutic positioning requires a comprehensive understanding from everyone involved if it is to be implemented successfully. Medifab with Simple Stuff Works is equipped to offer a range of training to support families, healthcare professionals and indeed anyone with an interest in therapeutic positioning. Our courses are accredited and license delegates to join a worldwide community of Postural Care Leaders with access to all our training materials and support.

Simple Stuff Works take pride in the fact that their equipment is silent, hygienic and simple to use. Training programs are available in Australia and New Zealand and will allow you to implement this powerful form of therapy as safely and as humanely as possible.



# Simple Stuff Works Packages

Packages are available in three sizes with common items

	Infant Package	Child Package	Adult Package
<b>Code</b>	7150-7001-C00	7150-7002-S00	7150-7003-S00
<b>User Age</b>	For infants from 0 to 5 years of age	For children from 5 to 12 years of age	For adults from 12+ years of age
<b>Package Includes</b>	<ul style="list-style-type: none"> <li>• Cot bed mesh</li> <li>• 2x cot toppers</li> <li>• 2x cot cotton sheets</li> <li>• 2x small lateral supports</li> <li>• 2x small wedges</li> <li>• Small roll pillow</li> <li>• Small horseshoe pillow</li> <li>• Travel bag</li> </ul>	<ul style="list-style-type: none"> <li>• Single bed mesh</li> <li>• 2x single toppers</li> <li>• 2x single cotton sheets</li> <li>• 2x medium lateral supports</li> <li>• 2x medium wide wedges</li> <li>• Standard horseshoe pillow</li> <li>• Size 3 supine stander</li> <li>• Travel bag</li> </ul>	<ul style="list-style-type: none"> <li>• Single bed mesh</li> <li>• 2x single toppers</li> <li>• 2x single cotton sheets</li> <li>• 2x medium lateral supports</li> <li>• 2x medium wide wedges</li> <li>• Standard horseshoe pillow</li> <li>• Size 4 supine stander</li> <li>• Travel bag</li> </ul>



## Simple Stuff Works Stabilising Mesh

- The base of the sleep system starts with a set of two layers of non-slip mesh. Supports are placed between the two layers around the individual to create the ideal position for them.
- One of the most hygienic bases for a positioning system being anti-microbial, complete with SteriTouch running through it, anti-mite, anti-fungal, entirely washable and conforms to Oeko-Tex Standard 100.
- Firmly holds any supports and is completely silent to use, making changes in position throughout the night as gentle and calm as possible.
- Can be used to also hold other pieces of equipment in place as well as holding the individual which means that the overall positioning system is very secure whilst allowing for easy changes in position.

	Cot Bed Size	Single Bed Size
<b>Code</b>	7150-7011-C00	7150-7012-S00
<b>Length</b>	<b>Base layer:</b> 1500mm <b>Top layer:</b> 1800mm	<b>Base layer:</b> 1800mm <b>Top layer:</b> 1800mm
<b>Width</b>	<b>Base layer:</b> 700mm <b>Top layer:</b> 1300mm	<b>Base layer:</b> 1000mm <b>Top layer:</b> 1500mm



## Simple Stuff Works Temperate Topper

- Pair of soft mattress toppers which is non-dynamic, neither cooling nor warming the individual so as not to interfere with their unique temperature regulatory system; an essential safety feature for those who find it difficult to control their own temperature.
- Used underneath the sheet on top of the layers of Stabilising Mesh.
- Made with a state-of-the-art material which is machine washable to 90°
- Wicks moisture away from the skin and allows air to circulate around an individual, helping to maintain a dry, temperate environment for them to rest in.
- Single size can also be used on larger beds.

	Cot Bed Size	Single Bed Size
<b>Code</b>	7150-7051-C00	7150-7052-S00
<b>Length</b>	1000mm	1000mm
<b>Width</b>	590mm	1000mm





# Simple Stuff Works Fitted Sheets

- Made from a beautiful material which is hard wearing but soft to ensure maximum comfort.
- Extra big and stretchy so that supports can be placed around an individual and the sheet will still conform to their body.
- Double bed sizes available on request.




	Cot Bed Size	Single Bed Size
Code	7150-7121-C00	7150-7122-S00
Length	1400mm	1900mm
Width	700mm	900mm

# Simple Stuff Works Wedges

- Made with a specialist fibre and covered in a soft yet durable polyester-mix material.
- Fibre is highly breathable and wicks moisture away from the skin allowing the person to remain comfortable whilst reducing the risk of tissue damage.
- Tapered at one end to facilitate use in a variety of positions.
- Also available in option for air mattresses designed to grip directly to the mattress surface. See scripting and ordering guide for more information.



	Small Narrow	Small Wide	Medium Narrow	Medium Wide	Large
Code	7150-7311-000	7150-7312-000	7150-7313-000	7150-7314-000	7150-7316-000
A	150mm	220mm	360mm	310mm	490mm
B	80mm	120mm	90mm	130mm	130mm
C	80mm	100mm	80mm	110mm	110mm
D	80mm	100mm	200mm	200mm	330mm




# Simple Stuff Works Padded Lateral Supports

- Perfect for providing firmer support for an individual using our therapeutic positioning equipment.
- L-shaped with soft, washable, removable pads which cover them.
- Positioned between the two layers of stabilising mesh which holds them in place securely.
- Can also be used to support other support items such as legs supports or pillows.
- Also available in option for air mattresses designed to grip directly to the mattress surface. See scripting and ordering guide for more information.



	Small	Medium	Large
Code	7150-7201-000	7150-7202-000	7150-7203-000
A	180mm	220mm	280mm
B	120mm	140mm	220mm
C	180mm	220mm	280mm




# Simple Stuff Works Supine Stabiliser

- Perfect for providing firmer support for an individual using our therapeutic positioning equipment.
- L-shaped with soft, washable, removable pads which cover them.
- Positioned between the two layers of stabilising mesh which holds them in place securely.
- Can also be used to support other support items such as legs supports or pillows.
- Available with either a standard terry toweling outer, or combination of wipeable inner and a terry toweling outer.



	Size 1	Size 2	Size 3	Size 4	Size 5
Standard	7150-7371-000	7150-7372-000	7150-7373-000	7150-7374-000	7150-7375-000
Wipeable inner/terry outer	7150-7381-000	7150-7382-000	7150-7383-000	7150-7384-000	7150-7385-000
A	350mm	420mm	490mm	580mm	720mm
B	240mm	290mm	300mm	320mm	380mm
C	100mm	130mm	180mm	220mm	220mm
D	130mm	170mm	190mm	210mm	210mm
E	90mm	90mm	100mm	110mm	110mm
F	40mm	70mm	70mm	70mm	70mm




# Simple Stuff Works Side Lying Leg Support

- Designed to protect the hips and pelvis when a person is lying on their side.
- The leg channel creates a perfect place for the top leg to lie and the bottom leg wraps gently around the outside.
- When combined with lots of support along the person's back, a supportive side lying position can be achieved.
- Perfect for individuals who need to lie on their side for safety reasons.



	Size 1	Size 2	Size 3
Standard	7150-7971-000	7150-7972-000	7150-7973-000
Wipeable inner/terry outer	7150-7981-000	7150-7982-000	7150-7983-000
A	270mm	400mm	400mm
B	380mm	560mm	610mm
C	280mm	340mm	400mm
D	80mm	100mm	110mm
E	130mm	200mm	150mm
F	80mm	90mm	70mm




# Simple Stuff Works Foot Support

- Designed to hold the feet in a neutral position for long periods of time.
- Made from the same fibre as the Temperate Toppers so they help to keep the feet feeling fresh whilst they are protected from falling forwards or rotating.
- Terry toweling outer makes for a comfortable and soft feel on the feet.
- Easy clip system makes it easy to put on and take off.




	Size 1	Size 2	Size 3
Code	7150-5351-000	7150-5352-000	7150-5353-000
A	210mm (max)	250mm (max)	300mm (max)
B	200mm	220mm	300mm



# Simple Stuff Works Neck Support Pillow

- Great option for supporting an individual's neck.
- Made from a soft fibre and covered in a comfortable yet durable polyester-mix material.


	Neck Support Pillow
Code	7150-4050-000
A	350mm
B	300mm



# Simple Stuff Works Horseshoe Pillow

- Made from a soft fibre and covered in a comfortable yet durable polyester-mix material.
- Can be used for a variety of different functions and is normally used to support the user's legs.
- Can be used between the two layers of stabilising mesh, or on top of the sheets

	Small	Large
Breathable inner/terry outer	7150-7711-000	7150-7712-000
Wipeable inner/terry outer	7150-7721-000	7150-7722-000
A	460mm	500mm
B	190mm	200mm



# Simple Stuff Works Head Positioning Pillow

- Universal supportive positioning pillow to provide the extra support where required.
- Can be used between the two layers of stabilising mesh, or on top of the sheets.
- Covered by a soft and comfortable, yet durable terry toweling cover.

	Small	Large
Code	7150-7501-000	7150-7502-000
Length	480mm	680mm
Width	350mm	450mm




## CareWave Cylindrical



- Can be used between legs, under heels and ankles in order to reduce pressure.
- Prevents the user from falling forward from a decubitus position.
- Provides sufficient support to achieve a corrected position for many users.
- Use for corrected side lying positioning for sufficient angles to keep the airway open and to ensure a neutral position in the involved joints.


	Small	Large
<b>Code</b>	7007-1502-045	7007-1501-066
<b>Length</b>	45 cm	70 cm
<b>Diameter</b>	14 cm	17 cm
<b>Cover Code</b>	7007-1502-270	7007-1501-270



## CareWave Abduction

- Use when the patient is lying down, placed in crotch area.
- Can also be placed under the arm.

	Small	Large	XX Large
<b>Code</b>	7007-1552-045	7007-1551-066	7007-1554-066
<b>Width</b>	12 cm	15 cm	90 cm
<b>Depth</b>	23 cm	30 cm	95 cm
<b>Height</b>	15 cm	20 cm	30 cm
<b>Cover Code</b>	7007-1552-270	7007-1551-270	



## CareWave Semi-Fowler

- Gives support and holds the user comfortably in a semi-fowler position.
- Place the cushion under the legs, the abduction module forward between the knees.


	Small	Large
<b>Code</b>	7007-1852-045	7007-1851-066
<b>Width</b>	50 cm	70 cm
<b>Depth</b>	35 cm	50 cm
<b>Height</b>	12 cm	14 cm
<b>Cover Code</b>	7007-1852-270	7007-1851-270



# CareWave Decubitus

- Supports shoulders, back, hips and knees of the user in a decubitus position. Reduces rubbing.
- Provides excellent trunk stability while the hip and knee angle can be varied.
- Lower part of the “S” can be drawn up between the knees so that the right median and good body symmetry can be achieved.
- Can also be shaped into a U-shaped cushion by adjusting the fill and rotating one end to thus provide a neutral position for the knee and hip joints, while the pad is pushed under the seat to fix the pelvis.
- The ends also provide support for the elbows.

	Small	Large	XX Large
<b>Code</b>	7007-3002-045	7007-3001-066	7007-3004-066
<b>Width</b>	25 cm	30 cm	36 cm
<b>Length</b>	140 cm	180 cm	200 cm
<b>Cover Code</b>	7007-3002-270	7007-3001-270	



# CareWave Decubitus Slim

- Designed for positioning of lower limbs to prevent them from crossing. Also prevents lateral sliding.


	Small	Large
<b>Code</b>	7007-3012-045	7007-3011-066
<b>Width</b>	9 cm	12 cm
<b>Length</b>	224 cm	300 cm
<b>Cover Code</b>	7007-3012-270	7007-3011-270



# CareWave Back Cushion

- Maintains the back in a semi-sitting position when user is sitting in bed or in a lateral decubitus position. Helps reduce pressure on ischium and heels.
- Can also be used in a sitting position in a chair.

	Small	Large
<b>Code</b>	7007-2452-045	7007-2451-066
<b>Length</b>	200 cm	200 cm
<b>Width</b>	40 cm	60 cm
<b>Cover Code</b>	7007-2452-270	7007-2451-270



## CareWave Body Bolster

- New addition to the CareWave range, which allows full support of the body from head to feet.
- An integrated pillow facilitates installation and user comfort and successfully maintains user position.
- When used in lateral position, the Body Bolster protects against pressures and avoids friction and pressure points on arms, knees and calves.

	Small	Large
<b>Code</b>	7007-3502-045	7007-3501-066
<b>Length</b>	135 cm	210 cm
<b>Width</b>	20 cm	35 cm
<b>Cover Code</b>	7007-3502-270	7007-3501-270



## CareWave Half-Moon

- It gives support and holds the patient comfortably in a lying or sitting position, and can be used to stabilize the upper body. The filling can be moved within the cushion so you can arrange it to adjust the head and arm height.
- If one wishes to counteract unwanted lateral movement of the body, the ends of the Half-Moon can be pushed under the hip area. Body weight will stabilize the pelvis and prevent unwanted movement. The idea is to eliminate movement that will take the body out of the desired symmetrical position.

	Small	Large
<b>Code</b>	7007-4052-045	7007-4051-066
<b>Length</b>	140 cm	180 cm
<b>Width</b>	25 cm	30 cm
<b>Cover Code</b>	7007-4052-270	7007-4051-270



## CareWave Ring Cushion

- Comfortably supports the user's head. If support is required especially on one side, the special filling can be shaken to one side of the cushion and positioned slightly skewed to one side if needed.
- Can be positioned with the opening towards the patient's visual direction which allows for a stable position for the head and neck, backwards and sideways.
- If the patient has the ability to make a lot of head movement the pillow can move out of place so the patient must be able to move the pillow or get help with this.




	Small	Large
<b>Code</b>	7007-4082-045	7007-4081-066
<b>Width</b>	35 cm	40 cm
<b>Cover Code</b>	7007-4082-270	7007-4081-270

## CareWave Heel Cushion

- Fixes to the mattress at the heel area with an adjustable strap.
- Supports ankles and relieves pressure on the heel zone.

	One Size
Code	7007-5201-066
Length	90 cm
Width	18 cm



## CareWave Heel Bootie

- Designed to protect and prevent the feet/heels from injury and bedsores and to provide comfort whatever the patient's position.


	Small	Large
Code	7007-5302-045	7007-5301-067
Length of foot surface	25 cm	30 cm



## Bow Cushion

- Can be used between the knees as an abduction cushion.
- Use under heels to reduce pressure.

	Small	Large
Code	7007-1602-045	7007-1601-066
Length	28 cm	34 cm
Width	19 cm	24 cm
Height	6 cm	10 cm
Cover Code	7007-1602-270	7007-1601-270





# Stabilegs

- Stablegs uses the same filling as the CareWave cushions.
- Positioning device has removable and adjustable modules.
- Abduction module (removable/adjustable) to limit friction and prevent leg-crossing in cases of spasticity.
- Adduction module (removable/adjustable) ensures alignment of lower limbs and pelvis stability.
- Heel relief and prevention equinus for prevention of bedsores in the heel area. End modules (removable/adjustable) are necessary to ensure pressure relief between the bed footboard and Stablegs.



	One Size
Code	7007-1831-066
Length	58 cm
Width	55 cm
Height	12 cm



# Fitleg Non-Slip

- Positioning device for proper orthopaedic alignment of lower limbs and full or partial relief on the heels.
- Prevents knee hyper-extension.
- Full relief on heels with the Non-Slip Module.
- High resiliency foam for stable positioning & optimized product life.
- Two distinct modules allow for customised patient positioning.



	Fitleg Non-Slip	Fitleg Module
Code	7007-1821-066	
Length	55 cm	50 cm
Width	58 cm	21 cm
Height	12 cm	10.5 cm



# Fitleg

- Positioning device for proper orthopaedic alignment of lower limbs which prevents front sliding, friction and shear. Anatomic shape provides maximum comfort and stable positioning.
- Fitleg Evo has abduction module to limit friction around knee articulations and prevent leg-crossing in case of spasticity.
- Prevention of knee hyper-extension due to the anatomic shape.
- Viscoelastic foam combined with HR foam base to provide maximum comfort.
- Reduces pressure on heels area with support on calves and thigh.



	Fitleg Small	Fitleg Large	Fitleg Evo
Code	7007-1802-066	7007-1801-066	7007-1811-066
Length	58 cm	58 cm	58 cm
Width	40 cm	55 cm	55 cm
Height	12 cm	12 cm	12 cm

# Ear Protection Cushion

- Comfortably supports the patient's head and relieves contact on the ear area in order to prevent bedsores.
- Comes standard with two cotton covers.

	One Size
<b>Code</b>	7007-4751-066
<b>Length</b>	30 cm
<b>Width</b>	30 cm



# Floor Pod

- Floor positioner which can be used for supine and prone positioning and postural drainage.
- This positioning cushion adapts to the user's body shape and comfortably keeps them in a half-seated position.
- Filling is not free-flow like a bean bag; unique CareWave filling supports the user in a stable position.

	Small	Large
<b>Code</b>	7007-0402-066	7007-0401-066
<b>Diameter</b>	110 cm	130 cm
<b>Cover Code</b>	7007-0402-270	7007-0401-270



# Hemi-Care Cushion

- Supports and holds the upper limbs in a stable position. Both parts of the Hemi Arm are fixed together with Velcro strips.
- Maintains the arm in a position which cannot slide and reduces pressure on elbow area so that there is no sensation of tear strength in the shoulder. For users with kyphosis, the Hemi-Care cushion helps maintain the back straight.

	Small	Large
<b>Code</b>	7007-5522-045	7007-5521-066
<b>Length</b>	45 cm	48 cm long
<b>Width</b>	30 cm	35 cm



## Kid'In'Form Starter Kit

- Includes 2 medium cylindrical rolls, large cylindrical roll, medium belt strap, 2 extra-small cylindrical rolls, backrest and comfort mat.
- Suitable for lesser supportive positioning requirements.

	Starter Kit
Code	7008-0002-279



## Kid'In'Form Ultimate Kit

- Includes 2 small cylindrical rolls, 2 medium cylindrical rolls, large cylindrical roll, 2 extra-small cylindrical rolls, pierced cushion, medium belt strap, large belt strap, backrest and comfort mat.
- Wide variety of positioning accessories allows flexibility of positioning options.
- If using in multi-user environments, will suit more users with increased sizes and options.

	Ultimate Kit
Code	7008-0003-279



## Kid'In'Form Positioning System

- Introductory lying and positioning system for early intervention. Portable and adaptable, with simple supports to facilitate active sitting, supine and prone positioning.
- Use the complete kit or components of the system to achieve the desired supported position.
- A versatile option for occupational therapists and physiotherapists who need to position the child securely in order to practice exercises or for stimulation/ therapy.
- Can also be used by parents and caregivers at home to offer positioning options to the child alternative to rigid positioning systems. Parents can position their child in a wide range of positions, to stimulate them and encourage perception of their surroundings.
- Recommended in the child's discovery and exploration phases and for directed practical work due to the following benefits:
  - Fabric covers feel pleasant to children and meet requirements of sensory tactile stimulation.
  - Different designs allow set-up of psychomotor exercises, to securely place the user and to offer a safe and relaxed environment.
  - Flexibility of the system encourages creativity and enables positions to be set for the child (to encourage protraction of the shoulders, to get down on hands and knees, to sit up, and more).




	Dimensions	Component
Extra-small cylindrical rolls	Ø 8 x 10 cm	
Small Cylindrical Roll	Ø 8 x 20 cm	
Medium Cylindrical Roll	Ø 8 x 70 cm	
Large Cylindrical Roll	Ø 11 x 50 cm	
Medium Belt Strap	60 x 14 cm	
Large Belt Strap	80 x 14 cm	
Backrest	76 x 24 x 17 cm	
Pierced Cushion	25 x 25 x 5 cm	
Comfort Mat	140 x 76 cm	

## Kid'In'Form Comfort Mat

- Covers 1m<sup>2</sup> area and allows combination of several positions of the supports.
- Top surface is jersey brushed soft fleece which is hook sensitive for attachment of the cushions.
- High density polyurethane foam ensures a sustainable cushioning and firm support.
- The underside is coated polyester which is waterproof and anti-slip.
- Machine washable at 40°. Do not bleach, tumble dry, iron or dry clean.
- Portable – can be carried using the handle.

	One Size
Length	140 cm
Width	76 cm



## Kid'In'Form Cylindrical Rolls

- Soft, waterproof, removable covers which can be machine washed.
- Rolls can be fixed to the mat or backrest using the self-fastening strips.
- Filling: machine wash at 40°. Do not bleach, tumble dry, iron or dry clean.
- Cover: machine wash at 95°, can be tumble dried. Do not bleach, iron or dry clean.

### Extra-Small Cylindrical

- Little cushion used for positioning, especially for support of the head. Can be very useful for complementing positioning.

### Small Cylindrical

- Used for support and pressure relief on bearing zones.

### Medium Cylindrical

- Long cushion is used for wrapping, for support and pressure relief on bearing zones.

### Large Cylindrical

- Long cushion with large diameter is used for support and pressure relief on bearing zones for larger users.



	Extra-Small	Small
Specifications	Ø 8 x 10 cm	Ø 8 x 20 cm
Colour	Purple	Blue

	Medium	Large
Specifications	Ø 8 x 70 cm	Ø 11 x 50 cm
Colour	Yellow-Green	Purple

## Kid'In'Form Pierced Cushion

- Cushion filled with memory foam fibres with central pressure relief.
- Centre hole assists with movement
- Soft, waterproof and removable cover which can be machine washed.
- Filling: machine wash at 40°. Do not bleach, tumble dry, iron or dry clean.
- Cover: machine wash at 95°, can be tumble dried. Do not bleach, iron or dry clean.

	One Size
Length	25 cm
Width	25 cm
Depth	5 cm
Hole	Ø 12 cm



## Kid'In'Form Belt Strap

- Medium belt strap helps to maintain smaller children on the mat or backrest.
- Large belt strap is used for older children.
- Belt strap can be used together with the cushions.
- Machine wash at 40°. Do not bleach, tumble dry, iron or dry clean.

	Medium	Large
Length	60 cm	80 cm
Width	14 cm	14 cm



## Kid'In'Form Backrest Cushion

- Positionable backrest is made of a structure in flexible foam that can be fixed to the mat with patches.
- Adjustable depending on the size of the user and positioning required.
- Can be used as a prone wedge so that the inclination assists stomach lying to enable the hands to be free for other activities.
- Cover allows fixture of other cushions to it and can be removed.
- Foam: Clean with sponge & soapy water. Replace cover only when foam is completely dry.
- Cover: Machine wash at 40°. Do not bleach, tumble dry, iron or dry clean.



	One-Size
Length	76 cm
Height	17 cm
Width	24 cm



## Kid'In'Form Notes

**KID'IN'FORM**

- This system is a medical device which must be used under permanent adult supervision.
- These belts, cushions and backrests are insufficient to keep the child in a secure position. This device must be used on the ground.
- It must not be used in case of any visible damage, delamination of fabric and foam and deterioration of seams. This system is designed to be used by children whose size is between 40 cm and 1m.
- Do not carry the child with the mat or the cushions



# Poz'In'Form Kit

- Transport bag
- Universal Cushion 55cm x 40cm
- Half moon cushion
- Semi-lateral decubitus cushion
- Cylindrical cushion
- Ring cushion
- Heel cushion with strap
- Occiput head cushion
- Pharmatex covers or Lenzing covers

### Effective Positioning

Poz'In'Form® cushions are designed to perform various support and comfort functions for long-term bedridden people. The microfibers prevent from creep effects that are regularly observed in classical positioning and from early deterioration of the support. The various forms, sizes and models of Poz'In'Form® allow for a wide range of positioning possibilities.


	With Pharmatex Cover	With Lenzing FR® Cover
<b>Code</b>	7009-0001-269	7009-0001-279



# Poz'In'Form Delta

- This cushion adapts to all morphologies. Placed under knees, it supports the knees' joint, prevents the patient from sliding downwards and reduces friction and shearing effects on the sacrum's area.
- It can be used in a semi-fowler position and helps to maintain a perfect position while relieving pressure in sensitive areas (sacrum, ischial tuberosities, heels) and shifting it to low-risk zones (calves, thighs and back).


	With Pharmatex Cover	With Lenzing FR® Cover
<b>Code</b>	7009-0551-269	7009-0551-279
<b>Length</b>	56 cm	56 cm
<b>Height</b>	28 cm	28 cm
<b>Depth</b>	32 cm	32 cm



# Poz'In'Form Universal

- Thanks to the three sizes, the universal cushion can be used in many different ways. In lying down the 55 x 40 cm size can be used for the knee flexion and a moderate hip abduction.
- The smaller cushions help to hold the patient in optimal positioning thus reducing bedsores risks and abnormal postures.

	With Pharmatex Cover				With Lenzing FR® Cover		
<b>Code</b>	7009-1251-269	7009-1252-269	7009-1253-269	<b>Code</b>	7009-1251-279	7009-1252-279	7009-1253-279
<b>Length</b>	35 cm	55 cm	40 cm	<b>Length</b>	35 cm	55 cm	40 cm
<b>Depth</b>	25 cm	40 cm	15 cm	<b>Depth</b>	25 cm	40 cm	15 cm



# Poz'In'Form Cylindrical

- Placed between legs, it reduces friction and shearing effects on thighs and knees' area.
- It can also be used under ankles to relieve pressure.


	With Pharmatex Cover	With Lenzing FR® Cover
<b>Code</b>	7009-1501-269	7009-1501-279
<b>Length</b>	75 cm	75 cm
<b>Diameter</b>	18 cm	18 cm



# Poz'In'Form Abduction

- It is used in crotch area when the patient is lying down to hold hips abducted and legs apart to reduce friction and shearing effects between legs, on the inner sides of knees and calves.


	With Pharmatex Cover	With Lenzing FR® Cover
<b>Code</b>	7009-1551-269	7009-1551-279
<b>Length</b>	35 cm	35 cm
<b>Height</b>	15 cm	15 cm
<b>Small Width</b>	10 cm	10 cm
<b>Large Width</b>	30 cm	30 cm



# Poz'In'Form Half-Moon

- In supine or semi-fowler position it supports and holds the back in proper position. It can also be used in (wheel) chairs for a comfortable and upright sitting position.

	With Pharmatex Cover	With Lenzing FR® Cover
<b>Code</b>	7009-1851-269	7009-1851-279
<b>Length</b>	200 cm	200 cm
<b>Width</b>	30 cm	30 cm






# Poz'In'Form Semi-Lateral Decubitus

- As multi-purpose cushion, it comfortably supports shoulders, back, hips and knees when the patient is lying in lateral position. Moreover, it reduces friction and shearing effects.


	With Pharmatex Cover	With Lenzing FR® Cover
<b>Code</b>	7009-3001-269	7009-3001-279
<b>Length</b>	185 cm	185 cm
<b>Width</b>	35 cm	35 cm



# Poz'In'Form Ring & Half Ring

- In combination with the cylindrical cushion, the half ring cushion reduces friction and shearing effects between legs, on the inner sides of knees and thighs and holds hips abducted. It can also be used under the neck.
- This cushion can be used in addition to universal, triangle or cylindrical cushions to support the posture correction. Put around the neck, this cushion comfortably supports the patient's head.

	With Pharmatex Cover	With Lenzing FR® Cover
<b>Half Ring Code</b>	7009-4051-269	7009-4051-279
<b>Ring Code</b>	7009-4081-269	7009-4081-279
<b>Length</b>	105 cm	135 cm
<b>Width</b>	22 cm	20 cm



# Poz'In'Form Occiput

- The hollow ensures an optimal positioning and a pressure relief on the occipital area. The combination of memory foam and microfiber allows a regular relief on the contact zones.

	With Pharmatex Cover	With Lenzing FR® Cover
<b>Code</b>	7009-4751-269	7009-4751-279
<b>Length</b>	50 cm	50 cm
<b>Height</b>	9 cm	9 cm
<b>Width</b>	40 cm	40 cm



## Poz'In'Form Heel Cushion

- It is placed at the bottom of the bed and fixed with an adjustable strap.
- The pressure on the heels is completely discharged.
- The heel cushion holds the ankles in proper position.

	With Pharmatex Cover	With Lenzing FR <sup>®</sup> Cover
<b>Code</b>	7009-5201-269	7009-5201-279
<b>Length</b>	85 cm	85 cm
<b>Diameter</b>	25 cm	25 cm



## Poz'In'Form Heel Wedge

- This cushion contributes to shift the pressure points from the high-risk heel area to low-risk zones (calves and knees).
- The heel is completely off the mattress and the weight of the lower leg is distributed evenly across the cushion. It holds the ankles in proper position and leads to a better blood circulation.
- Its size allows a pressure relief on the Achilles' heel.

	With Pharmatex Cover	With Lenzing FR <sup>®</sup> Cover
<b>Code</b>	7009-5401-269	7009-5401-279
<b>Length</b>	48 cm	48 cm
<b>Height</b>	12 cm	12 cm
<b>Width</b>	30 cm	30 cm



## Poz'In'Form Hand Wedge

- This positioning cushion for the hand reduces the pressure on the skin on the zones with risk of pressure sore. It is advisable after surgical operations.
- It reduces the number of pressure sensitive bearing points and allows a discharge of the elbows area. Its size allows a pressure relief on the wrist.

	With Pharmatex Cover	With Lenzing FR <sup>®</sup> Cover
<b>Code</b>	7009-5411-269	7009-5411-279
<b>Length</b>	39 cm	39 cm
<b>Height</b>	9 cm	9 cm
<b>Width</b>	22 cm	22 cm

