



medifab

Shaping Better Lives

SHARKY  
Specialised  
Tilt-in-Space  
Commode

Multifunctional  
commode for  
bath, shower  
and toilet use



## Versatile, practical, and functional!

These aspects have been the primary focus throughout the development and engineering of the Sharky commode system to make it easy and simple to use for parents and caregivers.

The Sharky features a 45° tilt-in-space mechanism with locking pins, guaranteeing safe tilting of the user whilst showering.

Height adjustment is easy to operate and the Sharky's two stainless steel gas pressure springs make it comfortable to control.

The Sharky has multifunctional applications, being able to be adapted for use over a toilet and in a bath. The bath frame adaptor along with the slim design of the removable seat makes it very easy to use in the bathtub, and with the seat unit so close to the tub floor, water consumption is reduced by not needing to fill the tub with as much water.

Unique amongst shower and bath commodes, the Sharky is easy to fold and reduces to a very compact folding size with a height of just 28 cm, all while having the seat still attached!

Available in sizes 1 and 2, with two back heights to choose from.



# Functional Benefits



## Removable Padding

Easy care removable padding for cleaning.



## Compact Folding

Fold down to compact height of just 28cm, even with seat attached, makes transport and storage easy.



## Back Support Recline

Simple to use system for easy and secure handling.



### **45° Tilt-in-Space**

The special construction of the tilt mechanism allows a generous 45° tilt of seat to ensure comfort and high safety especially when showering.



### **Height Adjustment**

Works with two stainless steel gas pressure springs which makes it easy and comfortable to use.



### **Hygiene Opening**

Removable section in rear area to assist with hygiene care.

# Accessories



## Seat Width Reducing Hip Supports

Size 1:

24cm: 6085-1110-000

27cm: 6085-1120-000

30cm: 6085-1130-000

Size 2:

29cm: 6085-1210-000

32cm: 6085-1220-000

35cm: 6085-1230-000



## Lateral Pads

11 x 6 cm: 6085-3310-000

14 x 8cm: 6085-3320-000



## Abduction Block

Size 1: 6085-1500-000

Size 2: 6085-1520-000



## Neoprene Ankle Huggers

Size 1: 6085-6101-000

Size 2: 6085-6102-000



## Neoprene Shoulder Harness

Size 1: 6085-6610-000

Size 2: 6085-6620-000



## Hip Belt with Neoprene Pads

Size 1: 6085-1910-000

Size 2: 6085-1920-000



## Neoprene Calf Strap

Size 1: 6085-6160-000

Size 2: 6085-6260-000



## Foot Support System

Size 1: 6085-1810-000

Size 2: 6085-1820-000



**Head Support**

6085-4750-000




**Push Handles**

6085-8600-000



**Seat Infill Pad**

6085-0401-000



**Toilet Adapter**

Size 1: 6085-0600-000

Size 2: 6085-0620-000



**Depth Adjustable  
Tray Table**

Size 1: 6085-7100-000

Size 2: 6085-7200-000




**Commode Pan  
with Retainer**

6085-8210-000



**Removable  
Splash Guard**

6085-7675-000



**Bath Subframe with  
Seat Angle Adjustment**

Size 1: 6085-0500-000

Size 2: 6085-0520-000



**Height-Adjustable  
Arm Support**

6085-5500-000

# Specifications

Back Support Height	Sharky Size 1	Sharky Size 2
35cm Back	6081-0010-100	N / A
45cm Back	6081-0010-000	6082-0015-100
55cm Back	N / A	6082-0015-000



**Sharky Size 1**



**Sharky Size 2**







A young child with short brown hair is sitting in a specialized medical chair. The chair is white with blue accents and has a blue safety harness with black straps. The child is looking off to the side with a neutral expression. The background is a wall of light-colored, vertically-oriented tiles.

medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)

DISTRIBUTED BY: