

**WonderSeat**  
FOR WHEELCHAIR

+ Zippie  
**IRIS**




## ULTIMATE TILT-IN-SPACE MANUAL WHEELCHAIR

AUSTRALIAN  
SCRIPTING AND ORDERING GUIDE

# Wheelchair Seating Package Options



Space




**Size 0**  
1420-1020-011

**Size 1**  
1420-1021-011

Wonderseat Constructa Cushion, 10" x 12"	Wonderseat Constructa Cushion, 12" x 16"
Growth-Adjust Spex Classic Back Support, 11-15" Growth	Growth-Adjust Spex Classic Back Support, 12-17" Growth
Drop Mount Hardware - Paediatric	
10"x12" Spex Seat Base with Quick Release Hardware	12"x16" Spex Seat Base with Quick Release Hardware
Spex Paediatric Standard Lateral Head Support Pad	
Wonderseat Quick-Release Head Support Hardware	
Narrow Straight Lateral Trunk Support Pad, Size 0, 8.5 x 11.5cm	Narrow Curved Lateral Trunk Support Pad, Size 1, 9.5 x 14.5cm
Paediatric Swing-Away Lateral Trunk Support Hardware	
Wonderseat Arm & Elbow Support Pad	Wonderseat Arm & Elbow Support Pad
Wonderseat Snug Thigh Support Pad, 19 x 8.5cm	Wonderseat Thigh Support Pad, 19 x 8.5cm
Wonderseat Side-Release Arm & Thigh Support Hardware	Wonderseat Complete-Release Arm & Thigh Support Hardware
Wonderseat Adjustable Medial Thigh Support, Pad height 7cm	
Wonderseat Padded Centrepoint Shoulder Harness	Wonderseat Padded Centrepoint Shoulder Harness
Wonderseat 4-Point Padded Hip Belt, Pad width 24cm	
Wonderseat Tray Table	
Wonderseat Tray Table Padding 1311-7410-013-C	
Spare Cushion Outer Cover 1102-1012-WSC-C	Spare Cushion Outer Cover 1102-1216-WSC-C

# Wheelchair Seating Scriptable Option

## Colour



## Cushion

### Wonderseat Constructa Cushion



Size 0, 10" x 12"  
1102-1012-WSC  
Size 1, 12" x 16"  
1102-1216-WSC

Cushion Spare Cover

## Back Support

### Growth-Adjust



Size 0, 12" x 11-15"  
1252-1211-OGA  
Size 0, 12" x 14-18"  
1252-1214-OGA  
Size 1, 14" x 12-17"  
1252-1412-OGA  
Size 1, 14" x 16- 21"  
1252-1416-OGA

### Growth-Adjust Spex Classic (Adaptable Contour)



Straight T Shape

Size 0, 12" x 11-15"  
1255-1211-OGA  
Size 0, 12" x 14-18"  
1255-1214-OGA  
Size 1, 14" x 12-17"  
1252-1412-OGA  
Size 1, 14" x 16-21"  
1255-1416-OGA

## Seat Base

### Quick Release Hardware



Size 0, 10" x 12"  
1008-1012-QRB  
Size 1, 12" x 16"  
1008-1216-QRB

### Fixed Seat Hardware



Size 0, 10" x 12"  
1008-1012-FXB  
Size 1, 12" x 16"  
1008-1216-FXB

## Back Support Hardware



Drop Mount Hardware - Paediatric  
1259-1400-00P



Single Mount Hardware  
1259-2100-100



Quad Mount Hardware  
1259-2400-100

## Head Support Pad



Circle 150  
1275-0011-000



Comfort 250  
1275-0031-000



Contour 250  
1275-0021-000



Paediatric Standard Lateral  
1275-0041-000



Paediatric Extended Lateral  
1275-0042-000

## Head Support Hardware



Wonderseat Quick-Release\*  
1274-3232-100



Wonderseat Lite\*  
1274-0232-100

## Lateral Trunk Support Hardware



Paediatric Axial Lateral  
1266-2613-100-[L/R]

Pair



Paediatric Swing-Away Lateral  
1266-2623-100-[L/R]

Pair

\*These products are designed to work as part of a complete Wonderseat Seating System, and are not compatible with other components of the Spex Seating range



## Lateral Trunk Support Pads

**Teardrop**



Pair

**Size 0, 14 x 8cm**  
1266-305[1/2]-400

**Broad**



Pair

**Size 1, 10.5 x 13 cm**  
1266-3031-400

**Size 2, 12.5 x 14 cm**  
1266-3032-400

**Narrow Straight**



Pair

**Size 0, 8.5 x 11.5cm**  
1266-3010-400

**Size 1, 9.5 x 14.5cm**  
1266-3011-400

**Size 2, 9.5 x 18cm**  
1266-3012-400

**Narrow Curved**



Pair

**Size 1, 9.5 x 14.5cm**  
1266-3021-400

**Size 2, 9.5 x 18cm**  
1266-3022-400

## Arm Support Pads

**Wonderseat Arm Support Pad\***



Pair

**Size 0, 9 x 20cm**  
1156-6410-310

**Size 1, 9 x 25cm**  
1156-6411-310

**Wonderseat Arm & Elbow Support Pad\***



Pair

**Size 0, 9 x 20cm**  
1156-6420-310

**Size 1, 9 x 25cm**  
1156-6421-310

## Arm & Thigh Support Hardware



Pair

**Wonderseat Side-Release\***  
1156-0403-130-[L/R]



Pair

**Wonderseat Complete-Release\***  
1156-0403-110-[L/R]

## Lateral Thigh Support Pads



Pair

**Wonderseat Thigh Support Pad, 5" deep, 11.5 x 8.5cm\***  
1156-1105-400



Pair

**Wonderseat Snug Thigh Support Pad, 8" deep, 19 x 8.5cm\***  
1156-1208-400-[L/R]



Pair

**Wonderseat Thigh Support Pad, 8" deep, 19 x 8.5cm\***  
1156-1108-400



Pair

**10" Thigh Support Pad, 26.5 x 9cm**  
1156-1010-400

## Medial Thigh Support



**Wonderseat Adjustable Medial Thigh Support, Pad height 7cm\***  
1162-1401-000



**Medial Thigh/Knee Block with Axial-Adjust Pad, Size 1, Pad height 9.5cm**  
1162-1601-000

## Tray Table



**Wonderseat Tray Table, 52cm x 42cm\***  
1311-7410-013




**Wonderseat Tray Table Padding\***  
1311-7410-013-C

## Anterior Pelvic Support



**Spex 2 Point Padded Hip Belt, Pad width 30cm**  
1405-6220-027



**Wonderseat 4-Point Padded Hip Belt, Pad width 24cm**  
1406-5420-027

## Anterior Trunk Support

**Wonderseat Padded Centrepoint Shoulder Harness\***



**Size 0, Pad height 22cm**  
1409-5660-017

**Size 1, Pad height 28cm**  
1409-5661-017

## Bit-Box



**Bit-Box\***  
Space Bounce Zoomy

\*These products are designed to work as part of a complete Wonderseat Seating System, and are not compatible with other components of the Spex Seating range

# Wheelchair Frame Package Options

## Standard Package:

1562-2012-001

Black frame only, supplied in a 12" x 14" seat width and depth configuration. Full growth range listed below.

### Includes:

- True-fit growth system
- Tie downs
- Width adjusting system
- 5" x 1.5" semi pneumatic front castors
- 16" Rear mag wheels; pneumatic with Airless Insert
- 90° legrest hangers with angle adjustable footplates
- Offset hanger mounts
- Height adjustable push handle
- Anti-tips
- Black Opal frame colour
- Foot operated tilt
- Foot operated brakes
- Medium seat depth range: 12-16"
- 0-55 degrees of rotation in space



# Wheel Lock

Push to Lock as additional to footlock

Push to Lock 6" Extensions

Drum Brakes 16" Pneumatic

1568-8455-000

1579-84-161930

1568-8452-16P/16S

# Footrest Hangers

Change to 80° Front Mount Hanger (pair)

(when purchased with chair)  
1568-5080-000



Front Mount Footrest Extension Tube

2" Pair  
1568-111505

4" Pair  
1568-111506

6" Pair  
1568-111507



# Footplates

## Heel Loops for Footplates

1403-6181-017



## Footplate Padding

1415-5312-300



## Spex Impact Guards

80 & 90  
Degree Hangers  
1415-5361-300



## Spex Calf Strap

12"  
1402-6171-017



# Tray

## Padded Travel Tray

**Small**  
1311-7441-209

(Hard tray included as standard with Wonderseat Space package).



# Canopy Options

## Bingo Canopy

Fits wheelchairs from 10" wide to 16" wide. Includes free rain cover suitable for smaller users.

1800-7500-000



# Notes

## Further Information



Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script: Zippie Iris\_20240513\_AU