

# Smirthwaite Bath Corner Chair



## BATHROOM AIDS

### SCRIPTING AND ORDERING GUIDE

Distributed by  
**medifab**

# Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

## Size Options

	Size 1	Size 2
Code	6101-4501-000	6101-4502-000
<b>Seat Unit Measurements</b>		
Seat height (mm)	75	75
Seat width (mm)	340	340
Seat depth (mm)	260	310
Seat depth with pads (mm)	245	295
Back support height (mm)	350	460
Load capacity (kg)	30	35
Product weight (kg)	3	4

### Base package includes:

- Bath corner chair complete with removable padding
- Medial thigh support
- Pelvic support belt





# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)

DISTRIBUTED BY: