

# Juni Chair



## INDOOR SEATING

### SCRIPTING AND ORDERING GUIDE

Distributed by  
**medifab**

# Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

## Size & Colour Options

	Size 1	Size 2	Size 3	Size 4
Black	3209-1001-099	3209-1002-099	3209-1003-099	3209-1004-099
Blue	3209-1001-060	3209-1002-060	3209-1003-060	3209-1004-060
Red	3209-1001-022	3209-1002-022	3209-1003-022	3209-1004-022
Lime	3209-1001-054	3209-1002-054	3209-1003-054	3209-1004-054

Seat Unit Measurements				
Seat height (cm)	16-22	20-28	24-34	30-40
Seat width (cm)	24	28	32.5	36.5
Seat depth (cm)	18-24	22-30	26-36	32-42
Back support height (cm)	26	29	32.5	40
Load capacity (kg)	30	40	50	60
Product weight (kg)	4	6	8	10

### Base package includes:

- Seat Unit
- Hip Belt



# Postural Support Accessories


Side Cushions - For Size 1 & 2 Juni	Side Cushions - For Size 3 & 4 Juni
-------------------------------------	-------------------------------------

Reduce seat width using side support (supplied as a pair).	
To reduce width by 5cm 3209-0001-025	To reduce width by 5cm 3209-0003-025
To reduce width by 10cm 3209-0001-050	To reduce width by 10cm 3209-0003-050
To reduce width by 15cm 3209-0011-050	To reduce width by 15cm 3209-0013-050



Lateral Supports - For Size 1 & 2 Juni	Lateral Supports - For Size 3 & 4 Juni
--	--

3209-3001-000	3209-3003-000
---------------	---------------



## Medial Thigh Support

Padded support for increased medial thigh support and seat contouring.	
<b>For size 1, 2 &amp; 3</b> 3209-1501-099	
<b>For size 4</b> 3209-1504-099	



## Anterior Chest Support

Size 1	Size 2	Size 3	Size 4
<b>Red</b> 3201-6401-022	<b>Red</b> 3201-6402-022	<b>Red</b> 3201-6403-022	<b>Red</b> 3201-6404-022
<b>Lime</b> 3201-6401-054	<b>Lime</b> 3201-6402-054	<b>Lime</b> 3201-6403-054	<b>Lime</b> 3201-6404-054
<b>Blue</b> 3201-6401-060	<b>Blue</b> 3201-6402-060	<b>Blue</b> 3201-6403-060	<b>Blue</b> 3201-6404-060
<b>Black</b> 3209-6351-099	<b>Black</b> 3209-6352-099	<b>Black</b> 3209-6353-099	<b>Black</b> 3209-6354-099



## Footsupport


Adjustable depth tray table, providing a work surface and an upper limb support.
<b>Size 1</b> 3209-5001-000
<b>Size 2</b> 3209-5002-000
<b>Size 3</b> 3209-5003-000
<b>Size 4</b> 3209-5004-000



# Base Options

Standard Skis		Standard Skis with Footboard	
Makes moving Juni chair in a straight line easier (under tables).		Standard skis with built in foot support, great for users who will push and could cause their chair to tip.	
<b>Size 1</b> 3209-5021-000		<b>Size 1</b> 3209-5061-000	
<b>Size 2</b> 3209-5022-000		<b>Size 2</b> 3209-5062-000	
<b>Size 3</b> 3209-5023-000		<b>Size 3</b> 3209-5063-000	
<b>Size 4</b> 3209-5024-000		<b>Size 4</b> 3209-5064-000	


Mobile Skis	
Makes the Juni chair easier to move, ideal for pushing under tables!	
<b>Size 1</b> 3209-5051-000	
<b>Size 2</b> 3209-5052-000	
<b>Size 3</b> 3209-5053-000	
<b>Size 4</b> 3209-5054-000	

Sliding Footboard (for use with mobile skis)	Height Adjustable Push Handle
Add foot support to mobile skis to support users legs.	Great for use with standard or mobile skis.
<b>Size 1</b> 3209-5071-000	
<b>Size 2</b> 3209-5072-000	
<b>Size 3</b> 3209-5073-000	
<b>Size 4</b> 3209-5074-000	
	<b>Size 1</b> 3209-8611-000
	<b>Size 2</b> 3209-8612-000
	<b>Size 3</b> 3209-8613-000
	<b>Size 4</b> 3209-8614-000

# Growth Accessories

## Backrest Extension (Increases back height as below)

Size 1: 11cm	Size 2: 11cm	Size 3: 13cm	Size 4: 16cm
<b>Red</b> 3209-2101-022	<b>Red</b> 3209-2102-022	<b>Red</b> 3209-2103-022	<b>Red</b> 3209-2104-022
<b>Lime</b> 3209-2101-054	<b>Lime</b> 3209-2102-054	<b>Lime</b> 3209-2103-054	<b>Lime</b> 3209-2104-054
<b>Blue</b> 3209-2101-060	<b>Blue</b> 3209-2102-060	<b>Blue</b> 3209-2103-060	<b>Blue</b> 3209-2104-060
<b>Black</b> 3209-2101-099	<b>Black</b> 3209-2102-099	<b>Black</b> 3209-2103-099	<b>Black</b> 3209-2104-099



## Outriggers

Increases seat to floor height by 5cm.

<b>Size 1</b> 3209-2811-000	<b>Size 2</b> 3209-2812-000	<b>Size 3</b> 3209-2813-000
--------------------------------	--------------------------------	--------------------------------



# Table Options

## Adjustable & Foldable Table      Quick Release Tray

Adjustable height work surface for children which can be easily folded for storage, a great companion to the Junii chair. This table comes in a beautiful 'natural' finish that complements the colour of the legs. (We don't currently offer the green finish shown in the image).

**Size 1**  
Table height 31-45cm  
3209-7191-000



**Size 2**  
Table height 40-60cm  
3209-7192-000



Adjustable depth tray table, providing a work surface and an upper limb support.

**Size 1**  
3201-7020-000

**Size 2**  
3201-7021-000

**Size 3**  
3201-7022-000

**Size 4**  
3201-7023-000





# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY: