

# Multiroller



**SPECIAL  
PURPOSE  
STROLLER**


SCRIPTING AND  
ORDERING GUIDE

Distributed by  
**medifab**

# Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

## Standard Package



<b>Multiroller Stroller Package</b>	
<p>Stroller frame with easy to fold function, with upholstered seat, height adjustable flip up footplate, interchangeable front end stroller frame wheelset, drum brakes and canopy.</p> <p>2610-1000-000</p>	

## Specifications



Push Handle Height (cm)	86 - 120
Overall Width (cm)	68
Wheel Diameter, Front/Buggy (cm)	18
Wheel Diameter, Rear (cm)	41
Folded size (with rear wheels on) (cm)	50 (H) x 60 (W) x 84 (L)
Folded size (with rear wheels removed) (cm)	34 (H) x 56 (W) x 84 (L)
Load Capacity Stroller (kg)	75
Weight of Stroller (kg)	9
Seat Width (cm)	30 - 38
Seat Depth (cm)	33
Lower Leg Length (cm)	20 - 30
Foot Support (cm)	38 (W) x 23 (L)
Back Support Width (cm)	30 - 40
Back Support Height (cm)	29 - 48

# Accessory Options

Running Wheel Add On	5 Point Harness
----------------------	-----------------

2610-8416-000		2610-6500-000	
---------------	---	---------------	---


Seat Depth Extension	Infill for Smaller Users
----------------------	--------------------------

2610-0281-000		2610-0393-000	
---------------	---	---------------	---

2 Point Hip Belt – Centre Release – Small	2 Point Hip Belt – Pen Key Release – Small
---	--

1405-6211-028		1405-6211-025	
---------------	---	---------------	---

Anti-Tip
----------

2610-8420-000	
---------------	---

## Notes

---

---

---

---



# medifab

Distributed by Medifab Ltd

[solutions@medifab.com](mailto:solutions@medifab.com)

[medifab.com](http://medifab.com)



DISTRIBUTED BY:

Script Multiroller\_20241016