

Jenx Atom



MODULAR SEATING SYSTEM

AUSTRALIAN
SCRIPTING AND ORDERING GUIDE

Distributed by
medifab

Client Details

Client/User Name:

Date:

NDIS/Funding:

Measurements	User Measurements
A Hip width	
B Seat depth	
C Seat to top of shoulder	
D Seat to top of head	
E Chest width	
F Seat to axilla	
G Leg length	

Product Specifications


Average age	6 Months-5 Years
Weight capacity	25kg
Seat width	18cm-27cm
Seat depth	9cm-34cm
Chest support width	20cm-28cm
Seat to top of backpad height	27cm-38cm
Foot support angle	-20° to 15°
Foot support height	11cm-36cm
Back support angle	-75° to 105°
Tilt range	30° posterior, 10° anterior

Get it Faster

9999-ASMB-002

Medifab Assembly Service for Easy Delivery

Medifab offers our popular assembly service on this product, meaning the product arrives ready to be delivered, straight from the box, saving time and speeding up the delivery process.



Therapists Normally Choose

Seat Unit/Colour Choice


Base Options


Red	Pink	Y-Base
<p>Includes:</p> <ul style="list-style-type: none"> Padded seat with memory foam Shaped back pad Arm supports Nursery tray Adjustable flip away hip support pads Push handle Upper back pad Flip up angle adjustable footplate <p>3431-1001-022</p>	<p>Includes:</p> <ul style="list-style-type: none"> Padded seat with memory foam Shaped back pad Adjustable flip away hip support pads Push handle Upper back pad Arm supports Nursery tray Shoulder protractors High rise flexi thoracic pads Thoracic pad covers 2 Link multigrip head support Flip up angle adjustable footplate <p>3431-0010-028</p>	<p>Black Y-base with lockable swivel castors Floor to seat height: 30-75cm Tilt range: 30° posterior 10° anterior 3615-8740-000</p> <p></p>

Head Supports

1 Link Multigrip	2 Link Multigrip	
Single link wings for some degree of lateral support. 33cm (W) x 5cm (D) x 7.5cm (H) For Red Atom: 3431-4601-022 For Pink Atom: 3431-4601-028	2 Links in each wing for extra support and positioning.  47.5cm (W) x 5cm (D) x 12cm (H) For Red Atom: 3431-4602-022 For Pink Atom: 3431-4602-028	
Wipeable Cover (Optional) 5356-4671-000	Wipeable Cover (Optional) 5356-4672-000	

Back Pad Options





Lumbar Pad	Back Pad Infill
Provides additional lumbar support. 3431-2501-000	
	A 50mm back pad to reduce the seat depth for smaller children. 3431-2370-000
	

Washable Cover	
For Seat and Back Support. Comes in Red. 3431-2470-022	

Shoulder & Trunk Supports

Shoulder Protractors	High-Rise Short Thoracic Pad	High-Rise Flexible Thoracic Pads
Flip away when not required. Independently adjustable for width and angle. 13cm (H) x 15.5cm (D) 3431-2750-022	Short, slim trunk support pads with single-handed flip away action. No strap. 3430-3151-000	Slimline and super flexible trunk support pads with straps and single-handed flip away action.  8cm (H) x 18cm (D) 3431-3251-022
		
	Cover (Optional) 3615-3171-022	Cover (Required) 3431-3271-022

Harnesses & Hip Belt

Waistcoat Harness	Wonderseat Harness	3 Point Lap Strap	4-Point Hip Belt
<p>Provides additional trunk support to children with low tone. 3615-6701-000</p>	<p>Size 0, Pad height 22cm Padded Centrepoint Shoulder Harness. 1409-5660-017</p>	<p>Dual lapstrap with additional safety of a groin harness. 3431-6251-000</p>	<p>Wonderseat Padded Hip Belt, Pad width 24cm. 1406-5420-027</p>
			

Leg Supports

Leg Positioning Pads	Long Leg Sitting Board	Abduction Block	Abduction Block
<p>Ideal for controlling windsweeping or creating an abducted seating position. Required fixed angle foot support. 3431-1751-000</p>	<p>For very young children or those needing hamstring stretch. Used in place of foot support. Comes in red. 3431-0281-022</p>	<p>Size 0 Length: 9cm Height above Seat: 8.5cm Overall height: 14.5cm Width Front: 9cm Width Rear: 7cm 3615-1650-000</p>	<p>Size 1 Length: 18cm Height above Seat: 7.5cm Overall height: 15.5cm Width Front: 10cm Width Rear: 8.5cm 3615-1651-000</p>
			

Foot Supports

Foot Support Pad
<p>Attaches to foot support for softer surface. Can only be used on the flip-up foot support. 3cm thickness. 3431-5701-000</p>


Foot Supports (cont.)

Sandals	Sandal Raising Blocks	Ankle Snugs
<p>Ideal for positioning feet where a firmer hold on the foot is required. 8cm (W) x 17cm (L) 5356-6081-000</p>	<p>Easily accommodate leg length discrepancies by giving sandals extra height where needed. 4 x 1.5cm Blocks 5356-6091-000</p>	<p>Features unique magnetic buckle fastening mechanism. Sandals required. Inner Circumference: 12cm-22.5cm 5356-6131-000</p>



Arm Supports

Covers for Arm Supports	
<p>(Red) Provide soft top for arm supports when tray is not used. 3431-5510-022</p>	

Tray & Accessories

Nursery Tray Infill	
<p>Designed to work with the Nursery Tray to reduce the cut-out for smaller children. 3431-7182-000</p>	

Nursery Tray Padding	Elbow Blocks
<p>Half Padding 3431-7442-000</p>	<p>Comes as a pair. Mounts to Nursery Tray corners for elbow support. 9.5cm (H) x 10cm (W) 3430-7350-000</p>
<p>Full Padding 3431-7452-000</p>	

Activity Frame	Grab Bar
<p>Great for encouraging head control and hand/eye co-ordination. Requires tray. 54cm (W) x 16cm (D) x 35cm (H) 3615-7475-000</p>	<p>Nursery Tray required. Bar dimensions: 54cm (W) x 9.5cm (D) x 7cm (H) 3615-7200-000</p>
	

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Jenx Atom_20241030_AU