

Jenx Junior+



**SEATING SYSTEM
FOR MILD, MODERATE
AND COMPLEX NEEDS**

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

Client Details


Client/User Name:	
Date:	
NDIS/Funding:	

User Measurements	
A Hip width	
B Seat depth	
C Seat to top of shoulder	
D Seat to top of head	
E Chest width	
F Seat to axilla	
G Leg length	

Get it Faster 9999-ASMB-002

Medifab Assembly Service for Easy Delivery

Medifab offers our popular assembly service on this product, meaning the product arrives ready to be delivered, straight from the box, saving time and speeding up the delivery process.



Notes

Specifications

	Size 1	Size 2	Size 3
Average age	18 Months–5 Years	3–8 Years	5–12 Years
Weight capacity (kg)	50	50	50
Seat width (cm)	15*–37.5	15*–37.5	15*–37.5
Seat depth (cm)	21–30.5	24–35.5	26–44
Chest support width (cm)	18.5–28	18.5–28	21–30
Seat-top of backpad height (cm)	34.5–51**	34.5–51**	41–57.5**
Footplate angle (°)	-35–5	-35–5	-35–5
Footplate height (cm)	19.5–27	23.5–35	25–40.5
Backrest recline (°)	-15–15	-15–15	-15–15
Tilt range (°)	30 (Posterior)–10 (Anterior)	30 (Posterior)–10 (Anterior)	30 (Posterior)–10 (Anterior)

* 15cm achieved only when using the hip reducing pads on Page 3.
 ** 51cm back height in Size 1 + 2 and 57.5cm in Size 3 only achieved with back post extension on Page 4.


Colour Options

Red	Lime	Blue
		
Fabric Cover -022	Fabric Cover -054	Fabric Cover -066
Wipeable Cover -W22	Wipeable Cover -W54	Wipeable Cover -W66




Medifab recommends that items outlined in blue are selected as standard package.

 **Therapists Normally Choose**

Base


Y-Base	
Height adjustment of 34 - 76cm, lockable swivel castors, tilt-in-space; 30° posterior, 10° anterior. 3615-8740-000	

Seat Unit


Size 1	Size 2	Size 3
Age range: 18 months–4 years. Max user weight: 50kg. 3610-1001-000	Age range: 3–8 years. Max user weight: 50kg. 3610-1002-000	Age range: 6–12 years. Max user weight: 50kg. 3610-1003-000
		
Contoured PU seat with removable abduction pad and memory foam seat pad, adjustable and flip-away PU hip pads, angle-adjustable footplate, 4-point hip belt, height-adjustable PU back pad, an independently adjustable PU lumbar/sacral pad, arm supports, tray table and push handle.		


Covers Required:		
Seat Cover Fabric 3612-0421-XXX	Seat Cover Fabric 3612-0423-XXX	Seat Cover Fabric 3612-0423-XXX
Back Cover Fabric 3612-2421-XXX	Back Cover Fabric 3612-2421-XXX	Back Cover Fabric 3612-2423-XXX
Hip Cover Fabric 3612-1072-XXX	Hip Cover Fabric 3612-1072-XXX	Hip Cover Fabric 3612-1072-XXX
Cover set (wipe clean only) 3612-0521-WXX	Cover set (wipe clean only) 3612-0522-WXX	Cover set (wipe clean only) 3612-0523-WXX

Hip Width Reducing Pads (Pair)


2.5cm	
Reduces Hip Width by 2.5cm on each side. 3612-1040-000	

Head Supports


Multigrip - 1 Link	Multigrip - 2 Link	Multigrip - 3 Link	
Single link wings for some degree of lateral support. 33cm (W) x 5cm (D) x 7.5cm (H) 5356-4601-000	2 Links in each wing for extra support and positioning. 47.5cm (W) x 5cm (D) x 12cm (H) 5356-4602-000 <div style="text-align: right; font-size: 24px; color: #00AEEF; font-weight: bold;">✓</div>	3 adjustable links in each wing for ultimate lateral headrest control and comfort. 62cm (W) x 5cm (D) x 12cm (H) 5356-4603-000	
Wipeable Cover (Optional) 5356-4671-000	Wipeable Cover (Optional) 5356-4672-000	Wipeable Cover (Optional) 5356-4673-000	

Oval Head Support	
19cm (W) x 10cm (H) 5356-4611-XXX	

Shoulder Protractors (Pair)

Size 0	Size 1	Size 2	
Can be used on all three size seat units. Pad size 13cm (H) x 13.5cm (D) 3612-2750-000	Can be used on all three size seat units. Pad size 21cm (H) x 14.5cm (D) 3612-2751-000	Can be used on all three size seat units. Pad size 23.5cm (H) x 17.5cm (D) 3612-2752-000	
Covers (required) 3612-2760-XXX	Covers (required) 3612-2761-XXX	Covers (required) 3612-2762-XXX	

Trunk Supports (Pair)


Flexi Thoracics			
Size 1 18cm (L) x 8cm (H) 3615-3251-000	Size 1 Long 22cm (L) x 8cm (H) <div style="text-align: right; font-size: 24px; color: #00AEEF; font-weight: bold;">✓</div> 3615-3252-000	Size 2 22cm (L) x 11cm (H) 3615-3253-000	
Covers (required) 3615-3271-XXX	Covers (required) 3615-3272-XXX	Covers (required) 3615-3273-XXX	

Waistcoat Harness

Size 1	Size 2	Size 3
A 5cm B 28cm C 12cm D 13.5cm 3615-6701-000	A 5cm B 32cm C 14cm D 15.5cm 3615-6702-000	A 6.5cm B 37cm C 17.5cm D 18cm 3615-6703-000




Back Post Extension

7.5cm	
7.5cm Back Post Extension for extra backrest height 3610-2480-000	


Abduction Block


Size 0	Size 1	
Length: 9cm Height above Seat: 8.5cm Overall height: 14.5cm 3615-1650-000	Width Front: 9cm Width Rear: 7cm Length: 18cm Height above Seat: 7.5cm Overall height: 15.5cm 3615-1651-000	Width Front: 10cm Width Rear: 8.5cm

Leg Positioning


Leg Control Pads	
Control leg position independently of pelvic position, flip completely out of the way to facilitate smooth and easy transfers. 12.5cm (H) x 7.5cm (D) 3615-1861-000	
Covers (required) 3615-3171-XXX	

Sandals


Size 1	Size 2	Size 3	
17cm (L) x 8cm (W) 5356-6081-000	22cm (L) x 9.5cm (W) 5356-6082-000	27.5cm (L) x 10.5cm (W) 5356-6083-000	

Sandal Raiser Blocks			
Size 1 Kit of 4x1.5cm blocks. 5356-6091-000	Size 2 Kit of 4x1.5cm blocks. 5356-6092-000	Size 3 Kit of 4x1.5cm blocks. 5356-6093-000	


Ankle Snugs

Size 1	Size 2 & 3	
Inner Circumference: 12cm - 22.5cm 5356-6131-000	Inner Circumference: 16.5cm - 29.5cm 5356-6132-000	

Foot Plate Pad

Foot Plate Pad	
Footplate pad not compatible with sandals. 4cm thick. 3615-5301-000	

Accessories for Arm Support & Tray

Tray Infill		
Great for attaching suction toys.		
Size 1 3610-7180-000	Size 2 3610-7182-000	

Half Tray Padding	Full Tray Padding
Comfort for forearms and elbows, whilst still enjoying a hard tray surface.	Added comfort for hands.
Size 1 3612-7441-000	Size 1 3612-7451-000
Size 2 3612-7442-000	Size 2 3612-7452-000

Grab Bar	Activity Frame
Used for children who need to fix one hand to increase functionality. 54cm (W) x 9.5cm (D) x 7cm (H) 3615-7200-000	Great for encouraging head control and hand to eye co-ordination. 54cm (W) x 16cm (D) x 35cm (H) 3615-7475-000

Elbow Blocks	Arm Support Covers
Used to encourage arms into a mid-line position and/or to block excessive retraction. 10cm (W) x 14cm (H) 3615-7350-000	3615-5671-XXX

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: JenxJunior+_20241024