

Sharky



SPECIALISED TILT-IN-SPACE COMMODE

Distributed by
medifab

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE

January 2025

Client Details

Client/User Name:		Notes
Date:		
DOB:		

Product Details

Options Available	Size 1	Size 2
35cm Back Support Height	6081-0010-100	N / A
45cm Back Support Height	6081-0010-000	6082-0015-100
55cm Back Support Height	N / A	6082-0015-000
Seat Width (cm)	24, 27, 30	29, 32, 35
Seat Depth (cm)	27 – 36	27 – 40
Maximum Traversable Toilet Height (cm)	52	
Tilt Adjustment Angle (°)	-10 to +45	
Back Support Recline Angle (°)	72 – 120	
Minimum Distance Between Lateral Supports (cm)	19	24
Leg Elevation Angle (°)	78 – 138	
Foot Support Angle (°)	80 – 100	
Arm Support Height (cm)	14.5 – 20.5	
Wheel Diameter (mm)	75	
Load Capacity (kg)	60	
Frame Length (cm)	70	
Frame Depth (cm)	56	



Sharky Size 1



Sharky Size 2

Upgrade Kit to Sharky Size 2

Upgrade Sharky mobility base from size 1 to size 2 using this kit.
6085-8081-086



Mobility Base (required)

(Sizes are specific for each size Sharky)

Sharky mobility base with height adjustable compact folding base, -10 to 45° tilt in space and goes over toilet heights from 42 - 53cm.	
For Sharky Size 1 6081-0015-000	For Sharky Size 2 6082-0010-000



Head Support

Suitable for use on either size Sharky

Sharky headrest, multi adjustable. 6085-4750-000



Lateral Supports

Both sizes suitable for use on either size Sharky

Sharky Lateral Pads. 11 x 6cm 6085-3310-000
Sharky Lateral Pads. 14 x 8cm 6085-3320-000



Hip Support (required)

Options for Size 1 Sharky

24cm		27cm		30cm	
Sharky seat width reducing hip supports for seat width of 24cm. 6085-1110-000		Sharky seat width reducing hip supports for seat width of 27cm. 6085-1120-000		Sharky seat width reducing hip supports for seat width of 30cm. 6085-1130-000	

Options for Size 2 Sharky

29cm		32cm		35cm	
Sharky seat width reducing hip supports for seat width of 29cm. 6085-1210-000		Sharky seat width reducing hip supports for seat width of 32cm. 6085-1220-000		Sharky seat width reducing hip supports for seat width of 35cm. 6085-1230-000	

Seat Options

Seat Infill Pad		Armrests	
Sharky removable seat infill pad for hygiene opening 6085-0401-000		Sharky removable and height adjustable armrests 6085-5500-000	
Medial Thigh Support (Both sizes suitable for use on either size Sharky)			
Size 1 8cm wide x 12cm deep 6085-1500-000	Size 2 10cm wide x 14cm deep 6085-1520-000		
Splash Guard		Commode Pan	
Sharky removable splash guard 6085-7675-000		Sharky removable commode pan with retainer 6085-8210-000	

Belt Options (recommended)

Pelvi-Loc (for use on either size Sharky)		Standard Hip Belt (for use on either size Sharky)	
Pelvi-Loc 3D for maximum pelvic security		Sharky hip belt with neoprene pads	
X-Small 1430-6340-024		Size 1 Pad Width 12.5cm 6085-1910-000	
Small 1430-6341-024		Size 2 Pad Width 16.5cm 6085-1920-000	

Foot Positioning

Foot Support System			
Sharky removable foot support system (These options are specific for each size Sharky).	Size 1 For Size 1 Sharky 6085-1810-000	Size 2 For Size 2 Sharky 6085-1820-000	
Calf Strap		Ankle Huggers (for use on either size Sharky)	
For Size 1 Sharky 6085-6160-000		Size 1 Ankle circumference 20-37cm approx 6085-6101-000	
For Size 2 Sharky 6085-6260-000		Size 2 Ankle circumference 28-47cm 6085-6102-000	

Harnesses


Shoulder Harness (for use on either size Sharky)		
Sharky neoprene shoulder harness.		
Size 1 24.5cm high x 19cm wide 6085-6610-000	Size 2 31.5cm high x 21.5cm wide 6085-6620-000	

Handles

Sharky Push Handles	
6085-8600-000	

Accessories

Toilet Adapter		
Size 1 (for size 1 Sharky) 6085-0600-000	Size 2 (for size 2 Sharky) 6085-0620-000	

Bath Subframe		
Sharky bath subframe with seat angle adjustment. Minimum width of bath: 36cm		
Size 1 (for size 1 Sharky) 6085-0500-000	Size 2 (for size 2 Sharky) 6085-0520-000	

Tray Table			
Sharky depth adjustable, removable tray table. Note: Arm Supports required for use with Tray Table	Size 1 (for size 1 Sharky) 6085-7100-000	Size 2 (for size 2 Sharky) 6085-7200-000	



Distributed by Medifab Ltd

solutions@medifab.com

medifab.com

A large, abstract image of a Sharky product, likely a shark-shaped container or bin, in blue and white. The word "SHARKY" is printed in blue, bold, uppercase letters on a white diagonal band across the product. The background is a light gray with a blue triangle in the top right corner.

SHARKY

DISTRIBUTED BY:

Script: Sharky_20250120_NZ