

BINGO Evo



PAEDIATRIC ADAPTIVE STROLLER

Distributed by
medifab

AUSTRALIAN
SCRIPTING AND ORDERING GUIDE

June 2025

Client Details

Client/User Name:	
Date:	
DOB:	

Get it Faster 9999-ASMB-002

Medifab Assembly Service for Easy Delivery

Medifab offers our popular assembly service on this product, meaning the product arrives ready to be delivered, straight from the box, saving time and speeding up the delivery process.

	Size 1	Size 2	Size 3
BINGO Stroller Package	2881-0008-000	2882-0008-000	2883-0008-000
BINGO Switch Stroller Package	2884-0001-500	2884-0002-500	N/A

Mobility Base Specifications			
Push handle height (cm)	75–110	75–110	
Overall width (cm)	63	68	
Wheel diameter, front/swivel (cm)	19	19	
Wheel diameter, rear (cm)	29	29	
Tilt in space (Forward and rear facing)	45°	45° (Size 3 – Forward facing only)	
Folded size, mobility base (cm) (with rear wheels on)	80 x 63 x 38	80 x 68 x 38	
Folded size, mobility base (cm) (with rear wheels removed)	80 x 63 x 31	80 x 68 x 31	
Weight, mobility base (kg) (swivelling front wheels)	8.6	8.7	
Load capacity, mobility base (kg)	60	60	

Seat Unit Specifications			
Seat width (cm)	18–30	23–35	23–35
Seat depth (cm)	16–34	22–40	30–45
Lower leg length (cm)	19–40	19–40	19–40
Back support height (cm)	45–63	54–75	62–83
Height of slots for shoulder straps (cm)	28–46	34–55	43–64
Seat to floor height (cm) (dependant on tilt position of seat)	46	46	46
Back support recline	90–160°	90–160°	90–135°
Weight of seat (kg)	6.4	7.4	8.4
Load capacity of seat unit (kg)	50	50	50
Weight of seat and mobility base combined (kg)	15	16.1	17.1

Stroller packages include:

- Matte black stroller frame with black wheels (BINGO Switch frame features interchangeable front end capabilities for mounting the Swivelling Front Wheel Attachment – page 6)
- Removable seat unit
- Tie-downs for transport
- Removable hip supports (width & depth adjustable)

- Fabric bag (max. 3kg)
 - Adjustable suspension set
 - Receiver for grab rail or therapy tray
 - Elevating foot support hanger (angle adjustable 90-160°)
- Note:** Upholstery not included in package

Upholstery Colour Options

Berry Soft		Turquoise Soft		Black	
Upholstery (required)					
Size 1 2881-8801-347		Size 1 2881-8801-349		Size 1 2881-8801-346	
Size 2 2882-8802-347		Size 2 2882-8802-349		Size 2 2882-8802-346	
Size 3 2883-8803-347		Size 3 2883-8803-349		Size 3 2883-8803-346	

Hip Support Upholstery (required)					
Size 1 2881-1075-347		Size 1 2881-1075-349		Size 1 2881-1075-346	
Size 2 2882-1076-347		Size 2 2882-1076-349		Size 2 2882-1076-346	
Size 3 2883-1077-347		Size 3 2883-1077-349		Size 3 2883-1077-346	

Tubie Fun Covers

Cover Sets

Cover set includes seat liner with waterproof layer along with covers for hip support, shoulder wing and spacer fabric head support pad.

(Note: Spacer fabric head support pads required. Must be used in conjunction with standard covers selected on page 2).

Llamas Cover Set - For Size 1 BINGO
2880-9561-121


Llamas Cover Set - For Size 2 BINGO
2880-9562-121

Llamas 5-Point Harness Cover (One size)
2880-6511-121

Planets Cover Set - For Size 1 BINGO
2880-9561-036

Planets Cover Set - For Size 2 BINGO
2880-9562-036

Planets 5-Point Harness Cover (One size)
2880-6511-036





Head Support Options

Spacer Fabric Head Support Pads		Head Support Speakers	
Size 1 11cm (D) x 16cm (H) 2880-4471-346		Requires Spacer Fabric Head Support Pads. 2880-8330-000	
Size 2 13cm (D) x 18cm (H) 2880-4472-346			

Occipital Head Support		Standard Lateral Head Support (12cm Deep)		Extended Lateral Head Support (13-15cm Deep)	
Merino/Spacer Fabric 2010-4110-099		Small (17cm Internal Width)		Small (17cm Internal Width)	
		Merino/Spacer Fabric 2010-4210-099		Merino/Spacer Fabric 2010-4310-099	
Wipe Down Cover 2010-4110-209		Wipe Down Cover 2010-4210-209		Wipe Down Cover 2010-4310-209	
		Large (21cm Internal Width)		Large (21cm Internal Width)	
Wipe Down Cover 2010-4110-209		Merino/Spacer Fabric 2010-4211-099		Merino/Spacer Fabric 2010-4311-099	
		Wipe Down Cover 2010-4211-209		Wipe Down Cover 2010-4311-209	

Lateral Support Options

Standard Lateral Supports			
Height and angle adjustable. All sizes are compatible with all BINGO seat units			
Size 1 8cm (H) x 10cm (L) 2880-3010-346	Size 2 10cm (H) x 12cm (L) 2880-3020-346	Size 3 12cm (H) x 13cm (L) 2880-3030-346	

Narrow Straight Lateral Supports		Broad Lateral Supports		Teardrop Lateral Supports	
Size 0 8.5cm (H) x 11.5cm (D) 2880-3008-100		Size 1 10.5cm (H) x 13cm (D) 2880-3008-310		Size 0 15cm (H) x 10.5cm (D) 2880-3008-500	

Pelvis and Lower Limb Support Options

Sacral Support (including spacer)	
22cm (W) x 9cm (H) x 1cm (D) 2880-2510-060	

Pelvis and Lower Limb Support Options

Thigh Wedge

22.5cm (W) x 10.5cm (D)
Height: 2cm tapering to 1cm.
(Recommended for
size 1 BINGO Evo only).
2006-0505-000



Trochanter Wedges

(Recommended for
size 1 BINGO Evo only).
2006-0660-000



Swing-Down Abductor Pommel (Fits all sizes)

Size 1
9cm (H) x 6cm (D)
2880-1651-346

Size 2
10cm (H) x 8cm (D)
2880-1652-346

Size 3
12.5cm (H) x 9cm (D)
2880-1653-346



Footplate Spacer

Spacer block to reduce lower leg length.
4cm thickness.

Size 1
2880-5321-000

Size 2 & 3
2880-5322-000



Ankle Footfasts

Requires technician to install. Purchase only.
Inner circumference listed below:

Size 1
16-25cm
2880-6101-000

Size 2
21-29cm
2880-6102-000



Hip Belts & Harnesses

Spex 2-Point Hip Belt

X-small
2 point slim line, dual pull
with neoprene sleeve padding.
1405-6220-027



Wonderseat 4-Point Hip Belt

Stabilises and controls the
pelvis. Can be quickly removed
for washing.
1406-5420-027



Groin Strap

To fit Size 1
2880-6361-346

To fit Size 2 & 3
2880-6362-346



Chest Shoulder Harness

Size 1
2880-6751-346

Size 2
2880-6752-346

Size 3
2880-6753-346



Wonderseat Padded Centrepont Shoulder Harness

Size 0
Pad height 22cm
1409-5660-017

Size 1
Pad height 28cm
1409-5661-017



Hip Belts & Harnesses

5-Point Harness

Includes front pull adjusters and neck padding.
2880-6511-346



Spex 5-Point Harness

Includes tee pad for babies and very small children. (Recommended for size 1 BINGO Evo only).
2006-6551-000



Tray & Grab Rail Options

Height Adjustable Tray

Features upgraded mechanism. (Cannot be used with Grab Rail).

To fit Size 1
2880-7121-000

To fit Size 2 & 3
2880-7122-000



Therapy Tray Infill (Fits all sizes)

For smaller trunk cut out.
2880-7180-000



Grab Rail

Features upgraded mechanism. (Cannot be used with therapy tray).

To fit Size 1
2880-7271-346

To fit Size 2 & 3
2880-7272-346



Critical Care Options

IV Pole (Fits all sizes)

With mobility bracket for attachment to mobility base.
2880-8275-000

Extra Bracket
For QR mounting of IV pole to second base
2880-8275-050



Oxygen Bottle Holder (Fits all sizes)

With mounting bracket for attachment to mobility base. 120mm diameter for 'C Type' bottle.
2880-8253-000



Chassis Mount Tray

Size 1
2880-8058-000

Size 2 & 3
2880-8059-000



BINGO Switch Accessories

Bingo Switch Swivelling Front Wheel Attachment		
Suitable for use with Bingo Switch frame only.		
Swivel Lock Suitable for both sizes. 2884-9580-000	Size 1 2884-9571-000	Size 2 2884-9572-000



Sibling Seat Mount	
(Suitable for use with Bingo Switch frame only).	
Mounting frame for toddler seat or capsule to enable transport of 2 children, compatible with Phil and Teds Voyager - adapters available for other brands (purchase separately)	
For Size 1 - 2884-4101-000 For Size 2 - 2884-4103-000	



(Sibling seat sold seperately).

Base Frame Options


Separate Seat Unit		
Separate seat unit for use with the Cobra Base if a combined stroller and seat unit is not purchased. This is also a great option if two seat units are preferred to avoid swapping the seat unit between the stroller and Cobra bases. Must be purchased with upholstery.		
Size 1 2881-0011-000	Size 2 2882-0012-000	Size 3 2883-0013-000



Cobra Indoor Base		
Cobra Hi-Low Indoor Frame with 75mm Twin Castors	Size 1 2880-8711-000	Size 2 & 3 2880-8721-000



Twin Bingo Stroller	
This twin stroller connection is only possible between strollers of the same size, or for connection of the mini and Size 1 together, it is not available for the size 3.	
Please fill out an additional script form to order the second stroller unit with a unique configuration and accessories.	
Note that this can only be fitted at the time of ordering both stroller units. (Medifab cannot retrofit this attachment). 2880-8630-000	
Attendant Brake Recommended for purchase with Twin Stroller configuration due to the increased combined weight and inertia when both strollers are connected. 2880-8452-000	



Accessories


Hand Strap

Wrist safety tether strap.
2005-6930-000



Swivel Lock (Fits all sizes)

To hold castors straight when brisk walking.
2880-8440-000



Canopy

Includes mounts and rain cover.
2880-7500-346



Notes

This image shows a single page of white paper with horizontal blue lines, resembling notebook paper. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: BINGO_20250630_AU