

BINGO Mini



EARLY INTERVENTION ADAPTIVE STROLLER

Distributed by
medifab

AUSTRALIAN
SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name:

Date:

NDIS/Funding:

Get it Faster

9999-ASMB-002

Medifab Assembly Service for Easy Delivery

Medifab offers our popular assembly service on this product, meaning the product arrives ready to be delivered, straight from the box, saving time and speeding up the delivery process.

Product Details

Bingo Mini Package Options

Standard Package	2880-0008-000
Switch Package (Standard package with Bingo Switch stroller base)	2884-0000-500
Conversion Kit from Mini to Size 1 (Stroller and Cobra) Converts BINGO Mini frames to accept the BINGO Evo size 1 seat unit, extending the growth of the solution (size 1 seat unit must be purchased separately).	2880-8080-086

Seat Unit		Mobility Base with swivelling wheels	
Seat width	13-25cm	Push handle height	75-110cm
Seat depth	12-27cm	Overall width	63cm
Lower leg length	16-36cm	Wheel diameter front (swivel)	19cm
Back support height	37-51cm	Wheel diameter rear	29cm
Height of slots for shoulder straps	25.5-39.5cm	Tilt in space (forward and rear facing)	45°
Seat to floor height (dependant on tilt position of seat)	50cm	Folded size mobility base (incl. rear wheels)	80 x 63 x 38cm
Backrest recline	80-180°	Folded size mobility base (excl. rear wheels)	80 x 63 x 31cm
Weight of seat unit*	7.8kg	Weight mobility base (swivelling front wheels)	9.8kg
Load capacity of seat unit	25kg	Load capacity mobility base	60kg
Weight of seat unit and mobility base combined*	17.6kg		

*Including upholstery and standard accessories

Stroller package includes:

- Matte black foldable stroller frame with black wheels
- Removable seat unit
- Tie-downs for transport
- Removable hip supports (width & depth adjustable)
- Fabric bag (max. 3kg)
- Adjustable suspension set
- Receiver for grab rail or therapy tray
- Elevating foot support hanger (angle adjustable 85-180°)
- **Note: Upholstery not included in package.**



Items outlined in blue are required as standard package

Upholstery Colour Options

Berry Soft		Turquoise Soft		Black	
Seat Upholstery (required)					
2880-8800-347		2880-8800-349		2880-8800-346	

Berry Soft		Turquoise Soft		Black	
Hip Support Upholstery (required)					
2880-1074-347		2880-1074-349		2880-1074-346	

Upholstery Upgrades

Upholstery Pressure Relief Foam Upgrades


Seat
2880-0391-000

Back
2880-2392-000

Head Supports

Head Support Pads




2880-4470-346




Head Supports (continued)

Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
Merino/Spacer Fabric 2010-4110-099		Small (17cm Internal width)		Small (17cm Internal width)	
		Merino/Spacer Fabric 2010-4210-099		Merino/Spacer Fabric 2010-4310-099	
		Wipe Down Cover 2010-4210-209		Wipe Down Cover 2010-4310-209	
Wipe Down Cover 2010-4110-209		Large (21cm Internal width)		Large (21cm Internal width)	
		Merino/Spacer Fabric 2010-4211-099		Merino/Spacer Fabric 2010-4311-099	
		Wipe Down Cover 2010-4211-209		Wipe Down Cover 2010-4311-209	

Lateral Supports

Standard Lateral Supports		Narrow Straight Lateral Supports	
Size 0 80mm (H) x 100mm (D) 2880-3009-346		Size 0 85mm (H) x 115mm (D) 2880-3007-100	
Broad Lateral Supports		Teardrop Lateral Supports	
Size 1 105mm (H) x 130mm (D) 2880-3007-310		Size 0 150mm (H) x 105mm (D) 2880-3007-500	

Sacral Support

Sacral foam insert including Spacer	
On velcro attachment. Also includes foam spacer which adheres behind Sacral insert where required. 10cm (H) x 18cm (W) 2880-2510-000	

Seat Options

Thigh Wedge	Trochanter Wedges
Velcro attached under upholstery. 22.5cm (W) x 10.5cm (D). Height: 2cm tapering down to 1cm 2006-0505-000	Option to contour standard flat cushion. Velcro attached under upholstery. 2006-0660-000

Abductor Pommel

Swing down abductor pommel

Size 1
 For increased abduction of the legs.
 9cm (D) x 6.5cm (W) x 9cm (H)
 2880-1651-346

Hip Belts & Harnesses

Spex 5-Point Harness		Wonderseat Harness		Wonderseat Hip Belt	
2006-6551-000		Wonderseat Padded Centrepoint Shoulder Harness - Size 0. Pad height 22cm 1409-5660-017		Wonderseat 4-Point Padded Hip Belt. Pad width 24cm. 1406-5420-027	

Foot Positioning

Footplate Cover		Footplate Spacer		Ankle Footfasts	
Neoprene fabric protects feet from the hard footplate surface. 2880-5305-000		Spacer block to reduce lower leg length by 4cm. 2880-5321-000		Strap inner circumference 15-23cm. 2880-6101-000	

Tray & Grab Rail Options

Height Adjustable Tray		Therapy Tray Infill	
Features upgraded mechanism. (Cannot be used with Grab Rail). 2880-7120-000		For smaller trunk cut out. 2880-7180-000	
Grab Rail			
Features upgraded mechanism. (Cannot be used with therapy tray). 2880-7270-346			

Critical Care Options

IV Pole		Oxygen Bottle Holder	
With mounting bracket for attachment to mobility base. Fits all sizes. 2880-8275-000		With mounting bracket for attachment to mobility base. 120mm diameter for "C Type" bottle. Fits all sizes. 2880-8253-000	
Extra Bracket For QR mounting of IV pole to second base. (Cobra or other indoor base). 2880-8275-050			

Chassis Mount Tray	
Size 1 2880-8058-000	

Swivel Lock

To hold castors straight when brisk walking	
2880-8440-000	

Hand Strap

Wrist safety tether strap	
2005-6930-000	

Canopy

Includes mini mounts and rain cover	
2880-7501-346	

BINGO Switch Accessories

Bingo Switch Swivelling Front Wheel Attachment

(Suitable for use with Bingo Switch frame only).

Size 1
2884-9571-000

Swivel Lock
2884-9580-000



Sibling Seat Mount

(Suitable for use with Bingo Switch frame only).

Mounting frame for toddler seat or capsule to enable transport of 2 children, compatible with Phil and Teds Voyager - adapters available for other brands (purchase separately)

2884-4101-000



(Sibling seat sold seperately).

Base Frame Options

Cobra Hi-Low Indoor Frame with 75mm Twin Castors

The seat unit easily transfers out of the stroller base and into the Cobra indoor base for feeding and activity times. Lower down the floor for playing and activities, and raise up to a table for meal time! Tilt settings can still be adjusted in the same way as the stroller base.

Size 0
2880-8720-000



Separate Seat Unit

Separate seat unit for use with the Cobra Base if a combined stroller and seat unit is not purchased. This is also a great option if two seat units are preferred to avoid swapping the seat unit between the stroller and Cobra bases. Must be purchased with upholstery.

One Size
2880-0010-000



Twin Bingo Stroller

Connection mount that allows two strollers to be connected one in front of the other for an additional sibling. The strollers are able to be easily detached and used separately where required, providing flexibility when travelling with or without siblings.

This is not available for use with a size 3 Bingo Evo. This can be used to connect either: bases of the same size (and corresponding sized Bingo Spex strollers), or the Bingo Mini and Evo Size 1. It is not suitable for use between any other size combinations.

Please fill out an additional script form to order the 2nd stroller unit with a unique configuration and accessories.

Please note that this can only be fitted at the time of ordering both stroller units (Medifab cannot retrofit this attachment).

2880-8630-000



Attendant Brake
Recommended for purchase with Twin Stroller configuration due to the increased combined weight and inertia when both strollers are connected.
2880-8452-000

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: BINGO Mini_20250630_AU