

spex for BINGO



FOR HOGGI STROLLERS

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE



Client Details

Client/User Name	
Date	
NDIS/Funding	

Notes

Select Size

		Size 1		Size 2	
BINGO Spex Stroller Package		2885-0031-032		2885-0032-032	
BINGO Spex Switch Stroller Package		2884-1001-000		2884-1002-000	
Seating specs	A Hip width	17-28cm (6.5-11")		22-33cm (8.5-13")	
	B Seat depth	18-30cm (7-12")		23-38cm (9-15")	
	C Back height	36-44cm (14-17")		40-48cm (16-19")	
	D Lower leg length*	25-40cm (10-16")		25-45cm (10-18")	
	E Trunk lateral width	12-28cm (5-11")		15-30cm (6-12")	
	F Trunk lateral height	15-34cm (6-13.5")		17-40cm (6.5-16")	
	G Arm support height	15-22cm (6-8.5")		15-22cm (6-8.5")	
	H Internal arm support width	20-32cm (8-13")		28-38cm (11-15")	
	I Head support height	45-65cm (18-25.5")		50-70cm (20-28")	
	J Cushion total width	33cm (13") Size 1.1	33cm (13") Size 1.2	38cm (15") Size 2.1	38cm (15") Size 2.2
	K Cushion total depth	33cm		38cm	
	L Internal cushion width	23cm (9") Size 1.1	28cm (11") Size 1.2	28cm (11") Size 2.1	33cm (13") Size 2.2
	M Centre of thigh contours	15cm (6") Size 1.1	17cm (6.5") Size 1.2	18cm (7") Size 2.1	20cm (8") Size 2.2
	N Thigh contour depth	10cm (4") Size 1.1	12cm (5") Size 1.2	12cm (5") Size 2.1	15cm (6") Size 2.2
	O Cushion side height	14cm (5.5")		14cm (5.5")	
	P Thigh contour height	9cm (3.5")		9cm (3.5")	
Q Abductor contour height	5cm (2")		5cm (2")		
R Back cushion width	33cm (13")		38cm (15")		
S Trunk cushion width	26cm (10")		30cm (12")		
Stroller specs	Frame seat depth capacity	18-30cm (7-12")		20-38cm (8-15")	
	Frame back height capacity	36-46cm (14-18")		41-55cm (16-22")	
	Backrest recline angle	90-160°			
	Height of push bar	75-110cm (29.5-43")			
	Overall frame width	63cm (25")		68cm (27")	
	Max load of seat unit	50kg (110lb) (body weight including other load)			
	Base frame capacity	60kg (132lb)			



Stroller packages include:

- Matte black frame with tie downs (BINGO Switch frame features interchangeable front end capabilities for mounting the Swivelling Front Wheel Attachment – page 6)
- Fabric basket (max load: 3kg)
- Tray table
- Sun canopy

Items outlined in red are required as standard package

Upholstery Colour Options


Jet	Royal	Chilli	Mango	Bumblebee
		 Express Delivery		
-099	-077	-022	-030	-044
Ocean	Leaf	Calypso	Granite	Ink/Wipe Down
 Express Delivery		 Express Delivery	 Express Delivery	
-060	-054	-028	-090	-209

(Express delivery on colours stocked in our Devonport Warehouse).


Cushion Options

Size 1.1	Size 1.2	Size 2.1	Size 2.2	
Fits BINGO Size 1		Fits BINGO Size 2		
33cm (13") total width 23cm (9") internal width (Recommended for 18-23cm seat depth) 2885-0610-XXX	33cm (13") total width 28cm (11") internal width (Recommended for 23-30cm seat depth) 2885-0611-XXX	38cm (15") total width 28cm (11") internal width (Recommended for 23-33cm seat depth) 2885-0612-XXX	38cm (15") total width 33cm (13") internal width (Recommended for 28-38cm seat depth) 2885-0613-XXX	
Spare Cover: 2885-0660-XXX	Spare Cover: 2885-0661-XXX	Spare Cover: 2885-0662-XXX	Spare Cover: 2885-0663-XXX	


Back Support Options

Spex Contouring Back System		
For Size 1 BINGO 2885-2012-XXX	For Size 2 BINGO 2885-2014-XXX	
Spare Cover: 2885-2062-XXX	Spare Cover: 2885-2064-XXX	



Head Supports

Standard Lateral		Extended Lateral		Mount and Stem Hardware	
With Stylo hardware & BINGO attachment					
For moderate head positioning needs. 2885-4241-XXX		For complex head positioning needs. 2885-4341-XXX		For use with other Spex head pads 2885-4510-000	
Spare Cover: 2885-4261-XXX		Spare Cover: 2885-4361-XXX			

Arm & Pelvic Supports

Arm Supports		Pelvic Supports	
With height-adjustable, quick-release hardware and neoprene cover.		Mounts onto arm support bracket (pair) behind cushion trough (arm supports must be selected).	
To fit Size 1 2885-5621-291		To fit Size 1 2885-1111-000	
To fit Size 2 2885-5622-291		To fit Size 2 2885-1112-000	

Lateral Supports

Narrow lateral supports		Broad lateral supports	
Soft padding/neoprene cover (pair) with Axial Technology. 8.5cm high x 11.5cm deep 2885-3410-291		Soft padding/neoprene cover (pair) with Axial Technology. 10.5cm high x 13cm deep 2885-3431-291	
			

Shoulder Harness

X-Small Centrepont	Small Centrepont	
Pad height 26cm, total pad width 22cm 1409-6660-017	Pad height 36cm, total pad width 32cm 1409-6661-017	

Hip Belt

2 Point Hip Belt		4 Point Hip Belt	
2-point, padded dual centre pull, side release buckle.		4-point, padded dual centre pull, side release buckle.	
X-Small (Pad width 30cm) 1405-6220-027		X-Small (Pad width 30cm) 1406-6420-027	
Small (Pad width 36cm) 1405-6221-027		Small (Pad width 36cm) 1406-6421-027	

Feet Positioning

Footplate Spacer Block*		Foot positioning sandals	Height kit / raiser blocks		
Size 1 (10/11") 2885-5321-000		May require height kit/s. 2885-6081-000		For foot positioning sandals (kit of 4). 2885-6091-000	
Size 2 (12/13") 2885-5322-000					

*Can not be used with foot positioning sandals

Critical Care Options

IV Pole (Fits all sizes)		Oxygen Bottle Holder	
Includes mobility bracket for attachment to mobility base. 2880-8275-000		With mounting bracket for attachment to mobility base. 120mm diameter for 'C Type' bottle. Fits both Size 1 and Size 2. 2880-8253-000	
Extra Bracket For QR mounting of IV pole to second base. 2880-8275-050			

Chassis Mount Tray	
For high dependency equipment - suction/oxygen etc.	
Size 1 (for BINGO Size 1) 2880-8058-000	
Size 2 (for BINGO Size 2) 2880-8059-000	

Accessories

Swivel Lock (Fits all sizes)		Hand Strap	
To hold castors straight when brisk walking. 2880-8440-000		2005-6930-000	
Therapy Tray Infill			
For smaller trunk cut out. 2880-7180-000			

Bingo Switch Accessories

Bingo Switch Swivelling Front Wheel Attachment		
Suitable for use with Bingo Switch frame only.	Size 1 2884-9571-000	Size 2 2884-9572-000
Swivel Lock Suitable for both sizes 2884-9580-000		



Sibling Seat Mount	
(Suitable for use with Bingo Switch frame only).	
Mounting frame for toddler seat or capsule to enable transport of 2 children, compatible with Phil and Teds Voyager - adapters available for other brands (purchase separately)	
For Size 1 - 2884-4101-000	
For Size 2 - 2884-4103-000	



(Sibling seat sold separately).

Base Frame Options

Separate Seat Unit		
Separate seat unit for use with the Cobra Base if a combined stroller and seat unit is not purchased. This is also a great option if two seat units are preferred to avoid swapping the seat unit between the stroller and Cobra bases. Must be purchased with upholstery.		
Size 1 2885-0011-032	Size 2 2885-0012-032	

Cobra high-low indoor frame with 75mm castors		
Size 1 (for BINGO Size 1) 2880-8711-000		
Size 2 (for BINGO Size 2) 2880-8721-000		



Manufactured by Spex Ltd

solutions@spexseating.com

spexseating.com



DISTRIBUTED BY:

Script: BINGOSpex_20250630_NZ