

Carrot 3000 Car Seat

Various extension accessories and multi-adjustable features allow for growth and adaptation as a child's needs change over time. The Carrot 3000 caters for safety and the importance of keeping the child in a supported posture for travel.

Our most comprehensive special needs car seat



STANDARD PACKAGE INCLUDES

- Carrot seat unit
- 5-point anterior chest support
- Headrest wings
- Tee-pad
- Free recline back support




BEHAVIOURAL PACKAGES:

Small Package with Anti-Escape Buckle:
4116-0010-026



Large package with Anti-Escape Buckle:
4116-0000-026


Small Package with Pen-Key Buckle:
4116-0010-025



Large package with Pen-Key Buckle:
4116-0000-025

Huge growth for pre-schoolers to older children

When requesting quote, please specify the user's shoulder height and seat depth.



The Carrot 3000 conforms to EU regulation ECE R44/04 and FMVSS 213 standard.

If user is too large for the Carrot 3000 Car Seat, the Carrot XL Seat is available for larger users – see page 56.

Size	Code	Age Range	User Weight Range	Seat Depth	Seat Width	Shoulder Height	Foot Support Height	Seat Weight
Small	4116-0010-259	3-15 years	13.6-49kg	25-40cm	31cm	32-57cm 45-67cm (when shoulder protector wings are added)	Small: 25-30cm Medium: 30-40cm Large: 40-50cm	8.7kg
Large*	4116-0000-259							

*Includes larger chest pad and deeper headrest wings.



Air Travel Friendly

The Carrot 3000 Car Seat has been approved for use in aircraft of some major airlines including the Federal Aviation Administration in the USA. This is NOT a guarantee that your chosen airline will allow the Carrot 3000 Car Seat to be used in their aircraft, and you must contact the airline prior to flying.

Accessories

Swivel Base Attachment*

4116-0180-104

Swivel Base Protection Plate

4116-0190-000

*Standard adjustment of mandatory footprop is 34cm-40cm

*Requires ISOFIX



Foot Support

Sml lower leg (D): 25-30cm

4116-5241-250

Med lower leg (D): 30-40cm

4116-5242-250

Lrg lower leg (D): 40-50cm

4116-5243-250



Swivel Base Footplate Adaptor

4116-5120-000




Seat Belt Holder

4116-6790-017




Removable Incontinence Cover

4116-0475-209



Rodded 50mm Seat Depth Pad

4116-0281-259



100mm Rodded Waterfall Seat Base Extension

4116-0283-259




50mm Seat Depth Growth Pad

4116-0285-259




100mm Seat Depth Growth Pad (only with Waterfall)

4116-0289-259



Abduction Pad

4116-1510-000



50mm Back Support Height Extension Infill Pad

4116-2481-259



100mm Back Support Height Extension Infill Pad

4116-2482-259



Shoulder Protector Wings

4116-2753-259



Shoulder Strap Retainer

Small: 4116-6721-000

Medium: 4116-6722-000

Large: 4116-6723-000



Cross Strap Anterior Chest Support

Anti-Escape 4604-3004-026

Pen Release 4604-3004-025

Side Release 4604-3112-019



Tether Strap

4116-6953-017



ISOFix Lower Straps (Optional)

(Seat belt must be used in conjunction)

4116-6910-000



Hinge Lock Belt

(for easy carrying)

4116-8170-000



Tray / Forearm Support

4116-5534-259



Tilt Wedge

4116-0291-259



Firm Contoured Lateral Trunk Support Wedges

4116-3250-090



Air Mesh Cushion Pad

4116-8880-099



See the Carrot 3000 script for our full range of accessories.