

Carrot 3000



CAR SEAT

Distributed by
medifab

AUSTRALIAN
SCRIPTING AND ORDERING GUIDE

Client Details


Client/User Name:	
Date:	
NDIS/Funding:	

User Measurements	
A Hip width	
B Seat depth	
C Seat to top of shoulder	
D Seat to top of head	
E Chest width	
F Seat to axilla	
G Leg length	

Get it Faster 9999-ASMB-002

Medifab Assembly Service for Easy Delivery

Medifab offers our popular assembly service on this product, meaning the product arrives ready to be delivered, straight from the box, saving time and speeding up the delivery process.



Notes

--

Seat Specifications

Age Range	User Weight Range	Seat Depth	Shoulder Height	Seat Weight
3-15 years	13.6-49kg	25-40cm	32-47cm 45-60cm (when shoulder protector wings are added)	8.7kg

Packages

Basic Packages:	Small		Large	
- Carrot seat unit 5-point harness - Headrest wings - Tee-pad - Free recline backrest	4116-0010-259		Larger chest pad and deep headrest wings. 4116-0000-259	
Behavioural Packages:	Small		Large - larger tee pad and deep headrest wings	
- Carrot seat unit 5-point harness - Headrest wings - Tee-pad - Free recline backrest - Supplied with emergency webbing slicer kit	Small package with anti-escape buckle. 4116-0010-026		Large package with anti-escape buckle. 4116-0000-026	
	Small package with pen-key buckle. 4116-0010-025		Large package with pen-key buckle. 4116-0000-025	

Swivel Base Attachment		Swivel Base Protection Plate	
ISOfix essential. (Standard adjustment of mandatory footprop is 34cm-40cm, if this is too short please contact Medifab Customer Service for Quote). 4116-0180-104 Select for footwell depth of 40-55cm only 4116-0186-104		Foam plate to protect car seat cover, mounts under swivel base. 4116-0190-000	

Seat Depth Growth


For growth, please select one of the following 3 rodded options, and associated extensions below, if required:

(Rodded components cannot be used in combination with each other).

Option 1

Rodded 50mm Seat Depth Extension		50mm Growth Pad	
4116-0281-259	+	4116-0285-259	
			

Option 2

Rodded Waterfall Seat Depth Extension		50mm Growth Pad*		100mm Growth Pad*	
4116-0283-259	+	4116-0285-259		4116-0289-259	
					

Option 3

Rodded Foot Support		50mm Growth Pad*		100mm Growth Pad*	
Small Lower leg (D): 25-30 cm 4116-5241-250	+	4116-0285-259		4116-0289-259	
Medium Lower leg (D): 30-40 cm 4116-5242-250					
Large Lower leg (D): 40-50 cm 4116-5243-250					
					
		Adaptor			
		For use of foot support with swivel base. 4116-5120-000			

*50mm and 100mm growth pads cannot be used in combination with each other.

Back Support Height Growth

50mm Height Extension	100mm Height Extension	Shoulder Protector Wings
50mm back support height extension infill pad. 4116-2481-259	100mm back support height extension infill pad. 4116-2482-259	(For children over 45cm shoulder height). 4116-2753-259

*50mm and 100mm infill pads cannot be used in combination with each other

Installation Accessories

Tether Strap (Optional)	ISOfix Lower Straps (Optional)	Airplane Fixing Strap
No installation option removes the need for the vehicle safety belt to be used, this must be used in every situation. 4116-6953-017	No installation option removes the need for the vehicle safety belt to be used, this must be used in every situation. 4116-6910-000	Required for securing to aircraft seat. 4116-6960-017

Cover Options

Wipe-down			
Wipe-down seat base 4116-0476-209		Tee-pad, Standard 4116-6553-209	
Wipe-down backrest (Main unit only). 4116-2476-209		Tee-pad, Large 4116-6558-209	
Removable Incontinence Cover		Air Mesh Cushion Pad	
Attaches using Velcro. Unable to be attached to wipe down covers. 4116-0475-209		Attaches using Velcro. Unable to be attached to wipe down covers. 4116-8880-099	

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Carrot_20250710_AU