

Carrot 3000



CAR SEAT

Distributed by
medifab

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

User Measurements	
A Hip width	
B Seat depth	
C Seat to top of shoulder	
D Seat to top of head	
E Chest width	
F Seat to axilla	
G Leg length	

Seat Specifications

Age Range	User Weight Range	Seat Depth	Shoulder Height	Seat Weight
3-15 years	13.6-49kg	25-40cm	32-47cm 45-60cm (when shoulder protector wings are added)	8.7kg

Packages

Basic Package:		
- Carrot seat unit 5-point harness - Headrest wings - Tee-pad - Free recline backrest	4116-0010-259	

Behavioural Packages:		
- Carrot seat unit 5-point harness - Headrest wings - Tee-pad - Free recline backrest - Supplied with emergency webbing slicer kit	Small Package with Anti-Escape Buckle 4116-0010-026	
	Small Package with Pen-Key Buckle 4116-0010-025	

Swivel Base Attachment		Swivel Base Protection Plate	
ISOfix essential. (Standard adjustment of mandatory footprop is 34-40cm). 4116-0180-104		Foam plate to protect car seat cover, mounts under swivel base. 4116-0190-000	
Select for footwell depth of 40-55cm only 4116-0186-104			

Seat Depth Growth

For growth, please select one of the following 3 rodded options, and associated extensions below, if required:

(Rodded components cannot be used in combination with each other).

Option 1

Rodded 50mm Seat Depth Extension

4116-0281-259



+

50mm Growth Pad

4116-0285-259



Option 2

Rodded Waterfall Seat Depth Extension	50mm Growth Pad*	100mm Growth Pad*
4116-0283-259 	+	
	4116-0285-259 	4116-0289-259 

Option 3

Rodded Foot Support		50mm Growth Pad*		100mm Growth Pad*	
Small Lower leg (D): 25-30 cm 4116-5241-250	+	4116-0285-259		4116-0289-259	
Medium Lower leg (D): 30-40 cm 4116-5242-250					
Large Lower leg (D): 40-50 cm 4116-5243-250					
		Adaptor For use of foot support with swivel base. 4116-5120-000			

*50mm and 100mm growth pads cannot be used in combination with each other.

Back Support Height Growth


50mm Height Extension	100mm Height Extension	Shoulder Protector Wings
50mm back support height extension infill pad. 4116-2481-259 	100mm back support height extension infill pad. 4116-2482-259 	(For children over 45cm shoulder height). 4116-2753-259 

*50mm and 100mm infill pads cannot be used in combination with each other

Installation Accessories

Tether Strap (Optional)	ISOfix Lower Straps (Optional)	Airplane Fixing Strap
No installation option removes the need for the vehicle safety belt to be used, this must be used in every situation. 4116-6953-017 	No installation option removes the need for the vehicle safety belt to be used, this must be used in every situation. 4116-6910-000 	Required for securing to aircraft seat. 4116-6960-017 

Cover Options

Wipe-down			
Wipe-down seat base 4116-0476-209		Tee-pad, Standard 4116-6553-209	
Wipe-down backrest (Main unit only). 4116-2476-209		Tee-pad, Large 4116-6558-209	
Removable Incontinence Cover		Air Mesh Cushion Pad	
Attaches using Velcro. Unable to be attached to wipe down covers. 4116-0475-209		Attaches using Velcro. Unable to be attached to wipe down covers. 4116-8880-099	

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Carrot_20250710_NZ