

Bed Comparison Chart



SPECIFICATIONS	GENIE BED	SAFE SURROUND 135	SAFE SURROUND 170	SAFESPACE	SAFESPACE HI LOW	COSYFIT	COSYFIT HI LOW	COSYFIT HI LOW PROFILING	SIESTA	SIESTA HI LOW	VOYAGER**
EXTERNAL HEIGHT - FLOOR TO TOP OF MATTRESS (CM)	42-102	55-95	55-95	28	59-99	15	40-80	40-80	34	59-101	20 (Approx) ***
EXTERNAL HEIGHT - BED BASE (NO MATTRESS) TO TOP OF SIDES (CM)	98	135	170	204	136	107, 122, 147	107, 122, 147	107, 122, 147*	140	142	161
EXTERNAL HEIGHT - FLOOR TO TOP OF SIDES (CM)	135-195	184	219	204	180-220	122, 137, 162	147/162/187 - 187/202/227	147/162/187 - 187/202/227	160	184-226	205 (Support post height)
EXTERNAL DIMENSIONS W X L (CM)	101 x 223	117 x 215	117 x 215	109 x 201 (Minimum), 320 x 320 (Maximum)*	107 x 214	109 x 198 or 109 x 213	110 x 234	110 x 219	116 x 225	116 x 225	149 x 223
INTERNAL HEIGHT - MATTRESS TO TOP OF SIDES (CM)	93 (Standard mattress)	91-131	126-166	169	113	122, 137, 162	92, 107, 132	92, 107, 132	125	125	161
INTERNAL MATTRESS DIMENSIONS W X L (CM)	90 x 200 (Before padding added)	100 x 200 (Before padding added)	100 x 200 (Before padding added)	91x183 or 122 x 183	91 x 198	91 x 183/198, 122 x 183/198, 183 x 183, 198 x 198	91 x 198	91 x 183	90 x 200	90 x 220	72.5 x 185
HIGH LOW BASE Y/N	✓	✓	✓	✗	✓	✗	✓	✓	✗	✓	✗
MOBILE BASE Y/N	✓	✓ (Optional accessory)	✓ (Optional accessory)	✗	✓	✗	✓	✓	✗	✓	✗
ELECTRONIC HEAD RAISE/KNEE BREAK Y/N	✓	✓	✓	✗	✗ (Manual option only)	✗ (Manual option only)	✗ (Manual option only)	✓	✗ (Manual option only)	✗ (Manual option only)	✗
PADDING Y/N	✓	✓	✓	✓ (Padded floor and corners with flexible sides)	✗ (Flexible sides)	✓ (Padded head and foot board with flexible sides)	✓ (Padded head and foot board with flexible sides)	✓ (Padded head and foot board with flexible sides)	✓ (Padded corners with flexible sides)	✓ (Padded corners with flexible sides)	✗
CANOPY Y/N	✗	✗	✗	✓	✓	✗	✗	✗	✗	✗	✓
COLOUR CHOICE Y/N	✗	✗	✗	✓	✓	✓	✓	✓	✓	✓	✗
OPENING MECHANISM	Hinged door - both sides	Hinged door - single side	Hinged door - single side	Zip - single or multiple sides	Zip - single or multiple sides	Zip - single or multiple sides	Zip - single or multiple sides	Zip - single or multiple sides	Zip - single or multiple sides	Zip - single or multiple sides	Zip - single side
WEIGHT LIMIT (KG)	150	150	150	120	100 (standard) 250 (heavy-duty)	120	120	120	120	120	120 (standard) 150 (heavy-duty)

* Risk assessment required for 147cm high sides on Cosyfit Hi Low Profiling
** Specifications refer to the standard Voyager model – please contact us for specifications for compact Voyager model
*** Depends on inflation levels for inflatable mattress