

Standing Frame Comparison Chart



SPECIFICATIONS	JENX MONKEY	DANDY			JENX MULTISTANDER			JENX STANDZ		ELLA STANDER				JENX SUPINE		
		SIZE 1	SIZE 2	SIZE 3	SIZE 1	SIZE 1 (ABDUCTION)	SIZE 2	SIZE 1	SIZE 2	SIZE 1	SIZE 2	SIZE 3	SIZE 4	SIZE 1	SIZE 2	SIZE 3
USAGE	Prone	Prone			Supine/Prone			Supine/Prone		Supine/Prone				Supine		
APPROX AGE*	9mths-4yrs	1-4yrs	4-8yrs	6-15yrs	9mths-6yrs	9mths-6yrs	3yrs-13yrs	1-9yrs	5-14yrs	1-3yrs	2-8yrs	6-12yrs	10-18yrs	9mths-4yrs	3-11yrs	10-Adult
MAX WEIGHT (kg)	25	35	45	55	35	35	60	45	70	30	45	55	75	25	50	100
MAX HEIGHT OF HEAD SUPPORT (cm)	N/A	85	125	145	116	116	155	138	160	95	130	155	180	106	146	190
FOOTPLATE TO TOP OF CHEST SUPPORT (cm)	51-75	39-57	60-86	76-110	49-85	50-89	75-120	58-100	82 - 116	40-60	56-88	77-106	89-113	The chest supports can be fitted at any available slot on the supine stander.		
CHEST SUPPORT WIDTH (cm)	18.5-28.5	17-20	25-30	25-35	14-27.5	14-27.5	20-32	17.5-30.5	22-36	14-37				13-31	13-31	22-44
HIP SUPPORT WIDTH (cm)	18.5-28.5	17-20	25-30	25-35	14-27.5	14-27.5	20-32	17.5-30.5	25-37	14-37				13-31	13-31	22-44
TILT RANGE (90° AS UPRIGHT)	20-90	19-90	14-90	11-90	10-90	10-90	0-85 (Electric) 10-85 (Gas)	5-90	0-90	0-90				10-90	10-90	5-90
TILT OPERATION	Gas	Gas			Gas			Gas		Electric				Gas		Electric
POSTURAL SUPPORT LEVEL	Mild/Mod	Mild			Mild-Complex			Mild-Complex		Mod-Complex				Mild-Complex		
FRAME FOOTPRINT (cm)	58x88.5	63x67	68x80	72x87	61.5x76	61.5x76	67.5x107	64.5x89	76x107	58x75	58x75	70x100	75x108	57.5x66	65x99	75x107
LEG ABDUCTION RANGE (degrees)**	0-approx. 10	10° each leg			0-approx. 10	0-30 each leg	0-approx. 10	0-30		0-33/45***	0-22/34***	0-21/30***	0-18/27***	0-approx. 10		
KNEE FLEXION CONTRACTURE ACCOMMODATION	✗	✓			Supine: ✓ Prone: ✗			Supine: ✓ Prone: ✗		✓				✓		

* Indicator only, for accurate product selection client measurements are the best guide.
**Indicates abduction range per leg, measured from neutral hip position (0 degrees)
***Maximum achievable abduction angle is relative to stander set up and height of user