

Stroller Comparison Chart



SPECIFICATIONS	BINGO MINI	BINGO EVO			BINGO SPEX		DURO		MULTIROLLER
		Size 1	Size 2	Size 3	Size 1	Size 2	Size 1	Size 2	
DESCRIPTION	A smart stroller base with ultimate early-intervention seating.	A versatile and modular stroller system for users with either basic or critical care needs.			A smart-looking stroller with superior postural support that enables parents to take their child outside and explore.		Configurable stroller for everyday challenges, with a great range of support and positioning accessories.		Incredibly light, incredibly tough!
APPROX AGE OF USER*	4mths-4yrs	1-5yrs	3-12yrs	8-15yrs	2-5yrs	3-7yrs	3-12yrs	8-15yrs	8-18
MAX LOAD (kg)	25	50			50		75		75
SEAT WIDTH (cm)	13-25	18-30	23-35	23-35	17-28	22-33	23-34	28-39	30-38
SEAT DEPTH (cm)	12-27	16-34	22-40	30-45	18-30	23-38	26-40	31-45	33
SHOULDER HEIGHT (cm)	25.5-39.5	28-46	34-55	42-63	36-44	40-48	30-48	42-62	29, 38, 48
FOOT SUPPORT HEIGHT RANGE (cm)	16-34	19-40	19-40	19-40	25-40	25-45	20-44	20-44	20 – 30
FRONT/REAR FACING	Both	Both			Both		Front		Front
INDOOR FRAME OPTIONS	✓	✓			✓		✗		✗
APPROX PRODUCT WEIGHT (kg)	17.5	15	16.1	18	29	30	19.2	21	9
TILT IN SPACE	45°	45°			45°		Static 16° or 20°		✗
RECLINE	80°-180° in 10° steps	90°-160°		90°-135°	90°-160°		78-120°		✗
FOOT SUPPORT ELEVATION	✓	✓			✓		✓		✗
GROWTH ADJUST SEAT	✓	✓			✓		✓		✓
SWIVEL FRONT WHEELS	✓	✓			✓		✓		✓
ACTIVITY TRAY	✓	✓			✓		✓		✗
SUN CANOPY	✓	✓			✓		✓		✓
RAIN SHIELD	✓	✓			✓		✓		✗
BASKET STORAGE	✓	✓			✓		✓		✗
POSTURAL SUPPORT LEVEL	Moderate/Complex	Moderate/Complex			Complex		Moderate		Mild
TYRE SOLID / PNEUMATIC	Solid	Solid			Solid		Solid		Pneumatic
TRANSPORT APPROVED	✓	✓			✓		✓		✗

* Indicator only, for accurate product selection client measurements are the best guide.