

Walking Aids Comparison Chart



SPECIFICATIONS	HOGGI CRUISER		FLUX WALKER				GRILLO			
	Size 1 (Pastel Blue)	Size 2 (Pastel Green)	Size 1	Size 2	Size 3	Size 4	Mini	Small	Medium	Large
DESCRIPTION	Delivering perfect support for users from 18-months to 16 years' old, the Hoggi Cruiser Walker really is the ultimate anterior walking frame.		Lightweight yet sturdy walking aid with unique design that promotes natural walking patterns.				Versatile aid that can be either posterior or anterior. Interchangeable supports make it suitable for wide range of users.			
APPROX AGE OF USER*	18-months to 16 years		1-6yrs	3-8yrs	8-18yrs	Adults	1-6yrs	3-8yrs	8-18yrs	Adults
MAX USER WEIGHT (kg)	50	65	60	75	75	100	35	40	90	110
USER HEIGHT RANGE (cm)	N/A		65-115	95-130	110-170	Adult	75-90	90-135	110-165	150-190
FOOT - AXILLA (cm)	N/A		N/A				52-77	70-100	76-125	81-143
FOOT - HIP (cm)	N/A		N/A				42-61	49-84	54-100	63-117
PRODUCT WEIGHT (kg)	5.6	5.9	6.5	6.8	7.1	9.1	13.5	Posterior: 15.7 Anterior: 16.2	Posterior: 17.7 Anterior: 19.5	Posterior: 22 Anterior: 23
ANTI-REVERSE MECHANISM	✗		✓				✓			
CHEST SUPPORT	✗		✗				✓			
HIP POSITIONER	✗		✓				✓			
FOREARM SUPPORTS	✗		✓				✓			
LEG DIVIDER/SEPARATOR	✗		✗				✓			
FOOTPRINT (W x L) (cm)	51x55	65x58	64x65	69x81	74x88	78x91	57x72	65x75	76x91	80x111
FOOT TO HAND GRIP HEIGHT (cm)	32-59	43-83	48-61 (Std...) 43-68 (uni.)	57-81 (Std...) 51-88 (uni.)	69-87 (Std...) 62-94 (uni.)	82-102 (Std...) 75-102 (uni.)	Variable			
SADDLE SEAT	✗		✗				✓			
SLING SEAT	✓		✓				✓			
ATTENDANT HANDLE	✗		✗				✓			
FRICTION/SPEED BREAK	✓		✓				✓			
FOLDING SEAT OPTION	✓		✓				Posterior Only			

* Indicator only, for accurate product selection client measurements are the best guide.