

Corzino Stroller



**FULLY FOLDABLE
STROLLER**

Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Size Options & Specifications

	Corzino 30	Corzino 34	Corzino 38	Corzino 42
Code	2510-0001-030	2510-0001-034	2510-0001-038	2510-0001-042
Approx. Age of User	4-8 years	6-10 years	8-12 years	10-14 years
Max. Load (kg)	45	55	75	75
Seat Width (cm)	30	34	38	42
Seat Depth (cm)	28-34	29-35	32-39	33-41
Seat to Floor Height (cm)	58	58	58	58
Back Support Height (cm)	55	67	65 (85 with extension – supplied as Std.)	69 (88 with extension – supplied as Std.)
Footplate Height (cm)	27-41	29-43	30-47	32-48
Approx. Stroller Weight (kg)	13	13.1	15	16
Folded Dimensions (cm)	117 (L) x 32 (W) x 34 (H)	122 (L) x 32 (W) x 36 (H)	122 (L) x 32 (W) x 36 (H)	126 (L) x 37 (W) x 36 (H)
Backrest Angle	Fixed	Fixed	Fixed	Fixed

Standard Package Includes:

- Adjustable seat depth
- Soft and padded upholstery
- Removable upholstery
- Comfortable 5-point support harness
- Comfort side protection
- Back support recline adjustment (2 positions)
- Height-adjustable and foldable foot supports
- Solid tyres
- Removable wheels
- Heel strap



Standard Package
(Storage Basket not included).

Showing Back Support Extension
(Included as Std. on Corzino 38 and 44).

Accessories

Extended Side Protection

Small (Corzino 30) 2510-6071-000	
Medium (Corzino 34) 2510-6072-000	
Large (Corzino 38/42) 2510-6073-000	

Upper Limb Support

2510-6140-000	
---------------	---

Sun Canopy

2510-7510-000	
---------------	---

Rain Cover - Must be used with Sun Canopy

For Corzino 30/34 2510-7522-000
For Corzino 38/42 2510-7523-000

Storage Basket

Small (Corzino 30/34) 2510-8010-100	
Large (Corzino 38/42) 2510-8010-200	

Pen Key Release Hip Belt

X-Small (30cm Pad) 1405-6210-025	
Small (36cm Pad) 1405-6211-025	
Medium (46cm Pad) 1405-6212-025	

Pen Key Release 5 Point Harness

2880-6511-025

Notes

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY: