

Duro Stroller



PAEDIATRIC ADAPTIVE STROLLER

Distributed by
medifab

AUSTRALIAN
SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name:	
Date:	
NDIS/Funding:	

Get it Faster 9999-ASMB-002

Medifab Assembly Service for Easy Delivery

Medifab offers our popular assembly service on this product, meaning the product arrives ready to be delivered, straight from the box, saving time and speeding up the delivery process.

Package Options

Specifications	Size 1	Size 2
Order Code	2780-0001-052	2780-0002-052
Seat Specifications		
Seat width (cm)	23-34	28-39
Seat depth (cm)	26-40	31-45
Seat angle (cm)	16 or 20°	16 or 20°
Back height (cm)	45-63	61-83
Back angle	78-120°	78-120°
Knee angle	85-160°	85-160°
Lower leg length (cm)	20-44	20-44
Folded size (L x W x H) (cm)	84 x 63 x 32	84 x 68 x 38
Base Specifications		
Max weight (kg)	75	75
Front wheel diameter (cm)	20	20
Rear wheel diameter (cm)	30	30

Stroller Package Includes:

- Swivelling Front Wheels
- Foot Support Bracket, Angle Adjustable, Flip Up
- Foot Support With Heel Stop
- Accessories Bag (Max Weight: 3kg)
- Tie Down Kit
- Retainer for Accessories
- Hip Guides (Low Version)
- **Note: Upholstery Not Included In Package**

Upholstery Colour Options


Turquoise	Black
Seat Upholstery (required)	
Size 1 2780-8801-349	Size 1 2780-8801-346
Size 2 2780-8802-349	Size 2 2780-8802-346

Head Supports


Head Support Pads	
Size 1 11cm (D) x 16cm (H) 2880-4471-346	
Size 2 13cm (D) x 18cm (H) 2880-4472-346	

Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
Merino/Spacer Fabric 2010-4110-099		Small (17cm Internal width)		Small (17cm Internal width)	
		Merino/Spacer Fabric 2010-4210-099		Merino/Spacer Fabric 2010-4310-099	
Wipe Down Cover 2010-4110-209		Large (21cm Internal width)		Large (21cm Internal width)	
		Merino/Spacer Fabric 2010-4211-099		Merino/Spacer Fabric 2010-4311-099	
Wipe Down Cover 2010-4110-209		Large (21cm Internal width)		Large (21cm Internal width)	
		Wipe Down Cover 2010-4211-209		Wipe Down Cover 2010-4311-209	

Lateral Supports

Standard Lateral Supports			
Height and angle adjustable. All sizes are compatible with both Duro sizes.			
Size 1 8cm (H) x 10cm (L) 2880-3010-346	Size 2 10cm (H) x 12cm (L) 2880-3020-346	Size 3 12cm (H) x 13cm (L) 2880-3030-346	
Teardrop Lateral Supports		Contoured Lateral Supports	
13.5cm (H) x 7.5cm (L) 2880-3008-500		Spex Discovery Axial Design 9cm (H) x 13cm (L) 2006-3030-099	

Harness Options

Spex Centrepont Harness		5-Point Harness	
X-Small Recommended for size 1 Duro Total pad dimensions: 26cm (H) x 22cm (W) 1409-6660-017		With front pull adjusters and neck padding 2880-6511-346	
Small Recommended for size 1 & 2 Duro Total pad dimensions: 36cm (H) x 32cm (W) 1409-6661-017			
Medium Recommended for size 2 Duro Total pad dimensions: 36cm (H) x 32cm (W) 1409-6662-017		Pen-Key Release Option 2880-6511-025	

Hip Belt Options

Spex Push Button		Pen Popper Buckle	
X-Small Pad width 30cm 1405-6210-028		Not supplied with Pen Popper tool – If required separately please order using product code: 1400-6010-025	
Small Pad width 36cm 1405-6211-028		X-Small Pad width 30cm 1405-6210-025	
Medium Pad width 46cm 1405-6212-028		Small Pad width 36cm 1405-6211-025	
		Medium Pad width 46cm 1405-6212-025	

Tray & Grab Rail Options

Therapy Tray		Therapy Tray Infill	
Cannot be used with Grab Rail			For smaller trunk cut out. Fits both tray sizes. 2880-7180-000 
To fit Size 1 2880-7111-000	To fit Size 2 2880-7112-000		

Grab Rail	
Cannot be used with therapy tray	
To fit Size 1 2880-7251-346	To fit Size 2 2880-7252-346

| | |

Hand Strap

Wrist safety tether strap

2005-6930-000



Swivel Lock

To hold castors straight when brisk walking

2880-8440-000



Canopy

Includes mounts and rain cover

Size 1
2780-7501-346

Size 2
2780-7502-346





Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY: