

Duro Stroller



PAEDIATRIC ADAPTIVE STROLLER

Distributed by
medifab

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE

October 2025

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Package Options

Specifications	Size 1	Size 2
Turquoise	2780-0011-349	2780-0012-349
Black	2780-0011-346	2780-0012-346
Seat Specifications		
Seat width (cm)	23-34	28-39
Seat depth (cm)	26-40	31-45
Seat angle (cm)	16 or 20°	16 or 20°
Back height (cm)	45-63	61-83
Back angle	78-120°	78-120°
Knee angle	85-160°	85-160°
Lower leg length (cm)	20-44	20-44
Folded size (L x W x H) (cm)	84 x 63 x 32	84 x 68 x 38
Base Specifications		
Max weight (kg)	75	75
Front wheel diameter (cm)	20	20
Rear wheel diameter (cm)	30	30

- Stroller package includes:**
- Hoggi Duro Rehab-Buggy
 - Hip Support, Aluminium
 - Seat and Backrest Upholsteries
 - Small Spex Hip Belt - Side-Release Buckle
 - Canopy
 - PU Wheels with Black Rim
 - Swivelling Front Wheels
 - Supporting Brackets
 - Foot Support Bracket, Angle Adjustable, Flip Up
 - Foot Support with Heel Stop
 - Accessories Bag (Max weight: 3kg)
 - Tie Down Kit (According to ISO 7176-19)
 - Retainer for Accessories
 - User Manual



Head Supports

Head Support Pads	
Size 1 11cm (D) x 16cm (H) 2880-4471-346	
Size 2 13cm (D) x 18cm (H) 2880-4472-346	


Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
Merino/Spacer Fabric 2010-4110-099		Small (17cm Internal width)		Small (17cm Internal width)	
		Merino/Spacer Fabric 2010-4210-099		Merino/Spacer Fabric 2010-4310-099	
Wipe Down Cover 2010-4110-209		Large (21cm Internal width)		Large (21cm Internal width)	
		Merino/Spacer Fabric 2010-4211-099		Merino/Spacer Fabric 2010-4311-099	
Wipe Down Cover 2010-4110-209		Large (21cm Internal width)		Large (21cm Internal width)	
		Wipe Down Cover 2010-4211-209		Wipe Down Cover 2010-4311-209	

Lateral Supports

Standard Lateral Supports			
Height and angle adjustable. All sizes are compatible with both Duro sizes.			
Size 1 8cm (H) x 10cm (L) 2880-3010-346	Size 2 10cm (H) x 12cm (L) 2880-3020-346	Size 3 12cm (H) x 13cm (L) 2880-3030-346	

Teardrop Lateral Supports		Contoured Lateral Supports	
Size 0 13.5cm (H) x 7.5cm (L) 2880-3007-500	Size 1 15.5cm (H) x 10.5cm (L) 2880-3008-500		Spex Discovery Axial design 9cm (H) x 13cm (L) 2006-3030-099 


Harness Options

Spex Centrepont Harness		5-Point Harness	
X-Small Recommended for size 1 Duro Total pad dimensions: 26cm (H) x 22cm (W) 1409-6660-017		With front pull adjusters and neck padding 2880-6511-346	
Small Recommended for size 1 & 2 Duro Total pad dimensions: 36cm (H) x 32cm (W) 1409-6661-017			
Medium Recommended for size 2 Duro Total pad dimensions: 36cm (H) x 32cm (W) 1409-6662-017			

Hip Belt Options

Spex Push Button		Pen Popper Buckle	
X-Small Pad width 30cm 1405-6210-028		Not supplied with Pen Popper tool – If required separately please order using product code: 1400-6010-025	
Small Pad width 36cm 1405-6211-028		X-Small Pad width 30cm 1405-6210-025	
Medium Pad width 46cm 1405-6212-028		Small Pad width 36cm 1405-6211-025	
		Medium Pad width 46cm 1405-6212-025	

Tray & Grab Rail Options

Therapy Tray		Therapy Tray Infill	
Cannot be used with Grab Rail			For smaller trunk cut out. Fits both tray sizes. 2880-7180-000 
To fit Size 1 2880-7111-000	To fit Size 2 2880-7112-000		
Grab Rail			
Cannot be used with therapy tray			
To fit Size 1 2880-7251-346	To fit Size 2 2880-7252-346		

Hand Strap

Wrist safety tether strap

2005-6930-000



Swivel Lock

To hold castors straight when brisk walking

2880-8440-000





Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY: