

Ella Stander



**STANDARD & COMPLEX
STANDING AID**


Distributed by
medifab

SCRIPTING AND ORDERING GUIDE

September 2025

Client details

Client/User name:	
DOB:	
Therapist:	
Date:	
User Measurements	
A Hip width	
B Axilla height	
C Chest width	
D Overall user height	
E Heel to knee	
F Hip height	
User weight	



Specifications

	Size 1	Size 2	Size 3	Size 4
Approx Age	0 – 4 years	3 – 7 years	6 – 12 years	11 – 18 years
User Height	60 – 95 cm	90 – 130 cm	120 – 155 cm	130 – 180 cm
Max User Weight	30 kg	40 kg	55 kg	80 kg
Trunk Pad Size	20 x 15 cm	25 x 20 cm	28 x 24 cm	34 x 28 cm
Trunk Width - Inset Laterals	12 – 20 cm	17 – 25 cm	19 – 30 cm	25 – 41 cm
Trunk Width - Straight Laterals	16 – 24 cm	21 – 29 cm	25 – 36 cm	31 – 47 cm
Hip Abduction Angle (Minimum)	0 – 18° per side	0 – 18° per side	0 – 18° per side	0 – 18° per side
Tilt Angle Adjustment	0 – 90°	0 – 90°	0 – 90°	0 – 90°

Package Options

Standard Configuration	Complex Configuration
	 Material cover is the same as the Standard Configuration
Package Includes: 1. Mobile base with powered tilt function, abduction set, four lockable castors, upright column and adjustment brackets for the feet and knees. 2. Sandals 3. Knee Cups with Buckle Strap 4. Pelvic Support 5. Trunk Support 6. Belt 7. Straight Laterals 8. Tray 9. Head Support with Hardware 10. Measuring Tape on Upright Extrusion 11. Tilt Angle Indicator Note: Standard package includes wider brackets at the knees and feet (not pictured).	Package Includes: Complex package includes all standard options with the following changes: 1. Multi Axis, rotational support boards for supporting flexion contractures or complex postures 2. Flexi laterals for trunk support 3. Rotational knee support upgrade 4. Shoulder Supports with Harness fitting kit

Product Codes		Size 1	Size 2	Size 3	Size 4
Standard Package	Manual:	5220-1021-000	5220-1022-000	5220-1023-000	N/A
	Powered:	5221-0020-000	5222-0020-000	5223-0020-000	5224-0020-000
Complex Package	Manual:	N/A	5220-1031-000	5220-1032-000	N/A
	Powered:		5222-0030-000	5223-0030-000	5224-0030-000

Support Sizes for Package:

Sandals	Small	Medium	Large	Large
Knee Cups	Small	Medium	Large	Large
Pelvic Lateral Supports	Small	Medium	Large	Large
Trunk Lateral Supports	Small	Medium	Large	Large
Support Belts	Size 1	Size 2	Size 3	Size 4
Shoulder Support (Complex Only)	N/A	Small	Large	Large
Head Support	10" Plush	10"	14"	14"

Build Your Own Stander

The options below allow you to build your own stander from scratch, specifically to the requirements of the user. Where components would be considered an "upgrade" to a package (they are not already within the package in some form) these are highlighted in blue.

Scan QR code
to configure stander
precisely for your client:



Base:				
Product Codes	Size 1	Size 2	Size 3	Size 4
Manual	5220-1181-000	5220-1182-000	5220-1183-000	N/A
Power	5220-0181-000	5220-0182-000	5220-0183-000	5220-0124-000

Foot Supports (Pair)		Knee Supports (Pair)	
Includes sandals with straps.		Includes knee cups, pads and straps.	
Small 13.4 x 7.6cm 5220-6071-000		Small Diameter: 8cm 5220-1761-000	
Medium 18.4 x 8.3cm 5220-6072-000		Medium Diameter: 10cm 5220-1762-000	
Large 23.5 x 11.4cm 5220-6073-000		Large Diameter: 12cm 5220-1763-000	

Knee Support Rotation - Module	
Enables the rotation of the Knee Support according to the rotation of the footplate. 5220-7200-000	

Pelvic Support Boards:				
Product Codes	Size 1 (21cm W x 16cm H)	Size 2 (26cm W x 21cm H)	Size 3 (29cm W x 25.5cm H)	Size 4 (35.5cm W x 29cm H)
Standard	5220-1911-000	5220-1912-000	5220-1913-000	5220-1914-000
Rotational (to support complex postures)	5220-1941-000	5220-1942-000	5220-1943-000	5220-1944-000

Trunk Support Boards:				
Product Codes	Size 1 (21cm W x 21cm H)	Size 2 (26cm W x 24cm H)	Size 3 (29cm W x 29.5cm H)	Size 4 (35.5cm W x 32.5cm H)
Standard	5220-2911-000	5220-2912-000	5220-2913-000	5220-2914-000
Rotational (to support complex postures)	5220-2961-000	5220-2962-000	5220-2963-000	5220-2964-000

Lumbar Support Pad

(Provides additional support for the lumbar spine when used between the pelvic and trunk support boards)

Small (26cm wide)

5220-2601-000

Medium (29cm wide)

5220-2602-000

Large (35.5cm wide)

5220-2603-000



Support Belts

(Support Belt required for standard swing-away laterals – optional for Flexi Pad swing away)

Size 1 (7.5cm W x 40-73cm L)

5220-6821-000

Qty:

Size 2 (10cm W x 45-33cm L)

5220-6822-000

Qty:

Size 3 (10cm W x 51-92cm L)

5220-6823-000

Qty:

Size 4 (12.5cm W x 61-112cm L)

5220-6824-000

Qty:



Lateral Supports - includes hardware and pads - for pelvic and trunk support (pair)

**Straight
Swing-Away**

Medium Swing-Away Straight

11.5cm D x 11.5cm H

5220-3321-000

Qty:

Large Swing-Away Inset 27mm Per Side

15.5cm D x 14cm H

5220-3330-000

Qty:

Large Swing-Away Straight

15.5cm D x 14cm H

5220-3331-000

Qty:

**Flexi-Pad
Swing-Away**

Medium Flexi-Pad Swing-Away Straight

21.5cm D x 12.5cm H

5220-3421-000

Qty:

Large Flexi-Pad Swing-Away Inset

27mm Per Side. 27cm D x 14cm H

5220-3430-000

Qty:

Large - Flexi-Pad Swing-Away Straight

27cm D x 14cm H

5220-3431-000

Qty:



Elbow Blocks - Compatible with swing-away laterals only

Small

(11.5cm W x 11.5cm H)


5220-7351-000


Large


(15.5cm W x 14cm H)

5220-7352-000

Standard Tray	Tray with Bowl
Tray size must match the size of the support board selected.	
Size 1 5220-7121-000	Size 1 5220-7181-000
Size 2 5220-7122-000	Size 2 5220-7182-000
Size 3 & 4 5220-7123-000	Size 3 & 4 5220-7183-000

Centrepoint Harness	
X-Small 1: 23cm 2: 26cm 3: 22cm 1409-6660-017	
Small 1: 29cm 2: 36cm 3: 32cm 1409-6661-017	
Medium 1: 34cm 2: 40cm 3: 36cm 1409-6662-017	
Large 1: 40cm 2: 45cm 3: 42cm 1409-6663-017	

Shoulder Protraction Pads (pair)	Harness Fitting Kit	
Includes pads and Harness Fitting Kit.		
Small 3cm W x 13cm H 5220-2751-000		
Medium 15.5cm W x 14cm H 5220-2752-000		

Adjust-A-Plush Head Support Options	Other Head Support Options	
Head Support Hardware Required when selecting an Adjust-A-Plush option.	For Ella Size 1 & 2 5220-4672-000	10" Plush (10" W x 5" H) Includes hardware. 5220-4610-000
	For Ella Size 3 & 4 5220-4674-000	
Adjust-A-Plush Standard 8" - 18" (W) x 5" (H) 1480-4330-281		14" Plush (14" W x 7.5" H) Includes hardware. 5220-4614-000
Adjust-A-Plush Medium 7" - 16" (W) x 4" (H) 1480-4320-281		
Adjust-A-Plush Narrow 6" - 14" (W) x 3" (H) 1480-4310-281		

Tool Free Adjustment Kit
Converts all postural adjustment points to hand operated. 5220-8160-000

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Ella Stander_20250904