

Empulse M90



Distributed by
medifab

Empulse m90

SCRIPTING AND ORDERING GUIDE

September 2025

Client Details

Client/User Name:	
Date:	
NDIS/Funding:	

Max. Range	15 km
Speed	6 kph
Batteries	Lithium Ion
Max. Safe Slope	6 degrees
Electronics	Bluetooth/Wireless
Min. Weight	0.6 kg (battery)
Weight of Heaviest Part	12 kg
Total Weight	5.3 kg (each wheel, without battery)
Max User Weight	100 kg

Empulse M90 Wheels Kit

Basic Package Includes:

- 2 x wheels
 - Marathon Plus Tyres
 - Black Anodised Handrims
 - Quick Release Axles
- 2 x wheel batteries
- 1 x battery dock charger
- 1 x Joystick Controller
- 1 x Controller Charger
- 1x Tall Fixed Controller Mount

1593-1000-000

Add: Mount Kit for Wheelchair

	Ki Mobility Catalyst 5X 1593-1071-000	Progeo Joker 1593-1065-000	Quickie 2 1593-1055-000
	Ki Mobility Focus CR 1593-1073-000	Progeo Tekna Advance 1593-1068-000	Quickie IRIS 1593-1054-000
	Ki Mobility Liberty FT 1593-1072-000	Progeo Duke 1593-1066-000	Quickie Nitrum 1593-1050-000
	Ki Mobility Catalyst 5 1593-1070-000	Progeo Ego, Ego Carbon 1593-1069-000	Ti-Lite AeroX, 2GX 1593-1064-000
Breezy Basix2 / Rubix2 1593-1060-000		Progeo ExcelleVario 1593-1067-000	Ti-Lite TR 1593-1062-000
Karma Flexx 1593-1075-000			Ti-Lite ZR 1593-1063-000
Panthera S3/U3 1593-1074-000			Ti-Lite ZRA, AERO Z, TRA, AERO T 1593-1061-000

Add: Controller Mount

	Leg Mount - Includes thigh strap and hardware 1593-1252-000
	Short Mounting Bracket M90 Clamp and Joystick Bracket - 80mm 1593-1254-000
	Attendant Mount assembly (Round tube attachment on suitable Quickie Chair models) - Clamp Included 1593-1155-000
Controller Mount Round Style	
Fixed Joystick Mount - Left 1593-1151-000	Fixed Joystick Mount - Right 1593-1150-000
Swing Away Mount - Left 1593-1153-000	Swing Away Mount - Right 1593-1152-000
Receiver (Required for mounting above) 1593-1154-000	

Optional Upgrade: Handrim

Plastic Coated		The Surge LT -Oval w/ Gription Strip	
1593-1351-000		1593-1353-000	
Natural Fit			
	No Thumb Grip 1593-1355-000	Standard Grip 1593-1354-000	Super Grip 1593-1352-000

Optional Upgrade: Tyre

Full Poly Tyre	
1593-1356-000	

Additional Accessories

Extra Single Battery		Travel Bag	
Batteries are interchangeable and can be used with either wheel. 1593-1251-000		Semi hard travel case for an entire set of M90 Wheels. 1593-1253-000	
Charger Port Cover			
1593-CCP-01-08-017			

Notes

This image shows a blank sheet of white paper designed for writing. It features horizontal blue lines spaced evenly down the page. Vertical red lines are positioned on either side, creating margins. The top edge of the paper has a faint header area. There is no handwriting or other markings on the page.A full-page view of a blank sheet of graph paper. The paper has a light gray background and is covered by a uniform grid of thin blue lines. The grid consists of small squares, typical of standard graph paper used for mathematics or engineering. There are no margins, text, or other markings on the page.

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Empulse M90_20250904