

FLUX Walker



VERSATILE POSTERIOR WALKING AID

Distributed by
medifab

NEW ZEALAND
SCRIPTING AND ORDERING GUIDE

Client Details

Client/User Name		Notes
Date		
NDIS/Funding		

Measurements

	Size 1	Size 2	Size 3	Size 4
Standard grab rail grip height	48–61cm	57–81cm	69–87cm	82–102cm
Standard grab rail grip width	41cm	46cm	51cm	51cm
Universal grip grab rail grip height	43–68cm	51–88cm	62–94cm	75–102cm
Universal grip grab rail grip width	29–41cm	34–46cm	39–51cm	39–51cm
Total frame width	64cm	69cm	74cm	78cm
Total frame length	65cm	80cm	88cm	91cm
Total load capacity	60kg	75kg	75kg	100kg
Front / rear wheel diameter	15cm / 20cm	15cm / 20cm	15cm / 20cm	20cm / 30cm

Package Options

Package includes:	Size 1	Size 2	Size 3	Size 4
- Standard Handlebar - Front Swivel Wheels - Polyurethan Tyres - Bumper Wheels - Friction Brakes - Anti-Reverse Mechanism	Luminous Red 5641-0001-000 	Sky Blue 5642-0001-000  (Shown with universal handlebar)	Gentian Blue 5643-0001-000  (Shown with universal handlebar)	Grey Aluminium 5644-0001-000 

Accessories

Anti-Tippers without Castors

Size 1-3 5640-8421-000		Size 4 5640-8424-000	
----------------------------------	---	--------------------------------	---

Anti-Tippers with Castors

Size 1-3 5640-8425-000		Size 4 5640-8428-000	
----------------------------------	---	--------------------------------	---

Fold Up Seat	Hip Pad
--------------	---------

Size 1-2 5640-0201-000		Multi-Adjustable 5640-1900-000	
Size 3-4 5640-0204-000			

Universal Handlebars (Not Compatible with Forearm Supports)

Upgrade Code Only - Must select at time of order.		
Size 1 5641-5751-000	Size 2 5642-5752-000	
Size 3 5643-5753-000	Size 4 5644-5754-000	

Driving Brake	Swivel Lock
---------------	-------------

Size 1-3 Includes wheel lock. 5640-8451-000		To hold castors straight when brisk walking. 2880-8440-000	
Size 4 Includes wheel lock. 5640-8454-000			

Accessories Continued

Dynamic Hip Control (Hip Belt Required)		Hip Belt (Dynamic Hip Control Required)	
Suspension, soft 5640-1952-000		Small Hip Size: 53-83cm 5640-6201-000	
Suspension, hard 5640-1953-000		Large Hip Size: 65-115cm 5640-6202-000	
Kit for orthoses 5640-1951-000			

Forearm Supports - 2 Axial Adjust		Carrier	
22 x 10.5cm 5640-5501-000		Size 1 5640-8061-000	
		Size 2 5640-8063-000	

Hand Disc		Flux Basket	
5640-5550-000		5640-8051-000	

Sling Seat		Sling Seat Rail	
Small Hip Size: 57-77cm 5640-1921-000		Size 1 5640-1931-000	
		Size 2 5640-1932-000	
Large Hip Size: 73-103cm 5640-1922-000		Size 3 5640-1933-000	

Notes



Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Flux_20250630_NZ