

Gravity Chair




Distributed by
medifab

INDOOR SEATING

SCRIPTING AND ORDERING GUIDE

Our product codes are meaningful!


Follow the 3 steps below, to build your gravity chair product code:



3 7 7 0 - 0 0 X X - 0 X X

EXAMPLE: A large Gravity Chair with a green classic liner, would have the code 3770-0013-055.

Step 1: Select Liner


Classic Version	Wipe Down Liner
Air breathable fabric on the user contact surfaces and rear back panel, with wipe-down covers in all other areas -001X- 	Fully wipe-down, including internal contact surfaces where user is seated -002X- 

Step 2: Select Size

Measurements	Small	Medium	Large	X-Large	XL-Wide
Order Code	-00X1-	-00X2-	-00X3-	-00X4-	-00X5-
A Seat Depth	24cm	30cm	38cm	43cm	43cm
B Seat Width at Back	21cm	23cm	27cm	30cm	41cm
C Seat Width at Front	28cm	30cm	37cm	39cm	50cm
D Back Height	60cm	69cm	74cm	81cm	81cm
E Seat Front Height	20cm	23cm	26cm	29cm	29cm
F Arm Support Height	29cm	35cm	39cm	43cm	43cm
G Total Width of Chair	46cm	49cm	57cm	60cm	71cm
H Total Height of Chair	55cm	61cm	67cm	71cm	72cm
I Total Depth of Chair	60cm	71cm	82cm	91cm	91cm
J Tilt Wedge Length	47cm	57cm	67cm	74cm	74cm
K Tilt Wedge Height	10cm	12cm	14cm	17cm	17cm
Weight Limit	100kg				

Gravity Chair Package includes:

- Foam chair with cover
- Coloured liner
- Standard head support
- Tilt wedge
- Velcro side pocket



Step 3: Select Liner Colour

Red Liner	Green Liner	Blue Liner	Black Liner
-022 	-055 	-066 	-099 

Head Supports

Occipital Head Support		Standard Lateral Head Support		Extended Lateral Head Support	
Merino/Spacer Fabric 2010-4110-099		Small (17cm internal width)		Small (17cm internal width)	
		Merino/Spacer Fabric 2010-4210-099		Merino/Spacer Fabric 2010-4310-099	
		Wipe Down Cover 2010-4210-209		Wipe Down Cover 2010-4310-209	
		Large (21cm internal width)		Large (21cm internal width)	
Wipe Down Cover 2010-4110-209		Merino/Spacer Fabric 2010-4211-099		Merino/Spacer Fabric 2010-4311-099	
		Wipe Down Cover 2010-4211-209		Wipe Down Cover 2010-4311-209	


Mobile Base

For ease of maneuverability in the home or classroom				
Small 3780-8701-000	Medium 3780-8702-000	Large 3780-8703-000		
XL 3780-8704-000	XL-Wide 3780-8705-000			

Tray Table

Encourages hand-to-eye coordination			
To fit Size Small and Medium 3770-7201-000	To fit size Large and XL 3770-7202-000	To fit size XL-Wide 3770-7203-000	

Foot Support


Small		Large	
(Suitable for Small or Medium Gravity chairs) 3770-1031-000		(Suitable for Large, XL and XLW Gravity Chairs) 3770-1032-000	

Abductor Pommel

Small		Large	
(Suitable for Small or Medium Gravity chairs) 3770-1021-000		(Suitable for Large, XL and XLW Gravity Chairs) 3770-1022-000	


Side Filler Pads

For reducing width of Gravity Chair by up to 6cm, adjustable in 1cm increments

Small 3770-1011-000	Medium 3770-1012-000	Large 3770-1013-000	XL/XLW 3770-1014-000	
-------------------------------	--------------------------------	-------------------------------	--------------------------------	---

Depth Adjustment Liner

For reducing the depth of the Gravity Chair

Small 3770-0321-000	Medium 3770-0322-000	Large 3770-0323-000	
XL 3770-0324-000	XLW 3770-0325-000		

Notes

[illegible]

medifab

Distributed by Medifab Ltd

solutions@medifab.com

medifab.com



DISTRIBUTED BY:

Script: Gravity Chair_20250121