










Grillo Dimensions & Weights

	MINI	SMALL	MEDIUM	LARGE
ANTERIOR CONFIGURATION	 TOTAL WEIGHT: 13.5 kg MAX LOAD: 35 kg	 TOTAL WEIGHT: 16.2 kg MAX LOAD: 45 kg	 TOTAL WEIGHT: 19.5 kg MAX LOAD: 90 kg	 TOTAL WEIGHT: 23 kg MAX LOAD: 110 kg
POSTERIOR CONFIGURATION	 TOTAL WEIGHT: 13.5 kg MAX LOAD: 35 kg	 TOTAL WEIGHT: 15.7 kg MAX LOAD: 45 kg	 TOTAL WEIGHT: 17.7 kg MAX LOAD: 90 kg	 TOTAL WEIGHT: 22 kg MAX LOAD: 110 kg
USER DIMENSIONS (C) Grillo mini 75 cm - 100 cm Grillo Small 90 cm - 135 cm Grillo Medium 110 cm - 165 cm Grillo Large 150 cm - 190 cm  IMPORTANT Make sure that measures A and B indicated in the figure do not exceed the suggested range indicated for the supports.				