



1- Purchasing information

P.O. number: * Date:

Company Name:*

Purchasing contact:*

Phone:*

Email:*

Contact preference:* Phone ☐ Email ☐

Marked for:

2- ATP/Clinician information

ATP/Clinician:*

Phone:*

Email:*

Contact preference:* Phone ☐ Email ☐

MA083 ☐

English Manual

MA084 ☐

German Manual

MA085 ☐

Italian Manual

MA086 ☐

Dutch Manual

MA087 ☐

Spanish Manual

MA088 ☐

French Manual

MA089

Portuguese Manual

MODEL

WC02	☐	Ultralight Helio A7 wheelchair 7000 series aluminum (120 kg/265 lb weight capacity)			
WC09	☐	HD Heavy duty kit (159 kg/ 350 lb weight capacity)			
WC10	☐	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6" camber			

Frame

SEAT WIDTH³ - cm (inches)

SW	☐	35.56 (14") ☐	38.10 (15") ☐	40.64 (16") ☐	43.18 (17") ☐	45.72 (18") ¹ ☐	48.26 (19") ¹ ☐		
	☐			50.80 (20") ¹ ☐	53.34 (21") ² ☐	55.88 (22") ² ☐			

SEAT DEPTH - cm (inches)

SD	☐	35.56 (14") □	38.10 (15") □	40.64 (16") □	43.18 (17") □	45.72 (18") □	48.26 (19") □			
	☐	50.80 (20") □								

TRANSIT TIE DOWN

AC04	☐	Transit Tie Down (<i>Unoccupied</i>)			
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Short belt 50-104 cm (20" - 41")			
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Long belt 50-116 cm (20" - 46")			
AC35	☐	WC-19 Transit Tie Down without belt ⁴			

1- Available with HD option: 2- Only available with HD option 3. Width calculator 4- May interfere with swing away armrest and Clamp on height adjustable push handle:

PACKAGES

VP01	☐	Vision Package: Orange accent added to point of contact accessories including: footrest latches, "T" and flip back armrest button cap, 6" Push to lock extension and rubber inserts on Newton composite angle adjustable footplates ¹	
------	--------------------------	--	--

1- Please select footplates, footrest, armrest, push to lock wheel lock and the extensions in the appropriate sections; To know the compatible parts, look for this logo ***VP01** ; Recommended with black chair;
See picture of the product on page 7.

A7 FRONT Seat-to-floor Height - Please select your choice

		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
Forks / Position		STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"
4 in 10,2 cm	P1	33.02cm 13"	35.56cm 14"	38.10cm 15"	34.29cm 13.5"	36.83cm 14.5"	39.37cm 15.5"												
	P2	34.29cm 13.5"	36.83cm 14.5"	39.37cm 15.5"	35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"									
5 in 12,7 cm	P1																		
	P2				35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"									
	P3	35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"						
	P4	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"	40.64cm 16"	43.18cm 17"	45.72cm 18"	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"			
7 in 17,8 cm	P1										41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"			
	P2							41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"			
	P3							43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20.5"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"
	P4	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	48.26cm 19"	50.80cm 20"	53.3cm 21"
	P5	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	48.26cm 19"	50.80cm 20.5"	53.3cm 21"	49.53cm 19.5"	52.07cm 20.5"	54.6cm 21.5"

Height available $\pm 0.65\text{cm}$ (1/4") (see stem bolt and fork combination at the end of the order form)

CASTERS				Availability		
CA01-03 to CA01-08	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6", 7", 8"		
CA10 to CA11	☐	Pneumatic casters (1 ¼")		6", 8"		
CA14 to CA16	☐	Newton UltraCaster Composite (1 ½")		4", 5", 6"		
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 ½")		4", 5", 6"		
CASTER OPTIONS						
CA12	☐	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) ¹				
CA13	☐	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)				

1- Not Available with WC-19

A7 - REAR Seat-to-floor Height

WHEEL SIZES		Height available		Specify desired height: _____
☐	20" (451) wheel	31.75cm/12.5 "	to 43.81cm/17.25"	
☐	22" (501) wheel	34.29cm/13.5 "	to 46.35cm/18.25"	
☐	24" (540) wheel	36.83cm/14.5 "	to 48.89cm/19.25"	
☐	25" (559) wheel	38.10cm/15 "	to 50.16cm/19.75"	
☐	26" (590) wheel	39.37cm/15.5 "	to 51.34cm/20.25"	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1.25 cm (½") of selected Seat to Floor height"

Footrest						
SWING-IN/SWING-OUT FOOTREST						
FR03-01	☐	60° Swing-in/swing-out footrest	34.29 - 43.18cm (13.5" - 17")	Set at*:		
FR03-02	☐	60° Swing-in/swing-out footrest	36.83 - 49.53cm (14.5" - 19.5")	Set at*:		
FR03-03	☐	60° Swing-in/swing-out footrest	39.29 - 52.07cm (15.5" - 20.5")	Set at*:		
FR04-01	☐	70° Swing-in/swing-out footrest	34.29 - 43.18cm (13.5" - 17")	Set at*:		
FR04-02	☐	70° Swing-in/swing-out footrest	36.83 - 49.53cm (14.5" - 19.5")	Set at*:		
FR04-03	☐	70° Swing-in/swing-out footrest	39.37 - 52.07cm (15.5" - 20.5")	Set at*:		
FR15-01	☐	80° Swing-in/swing-out footrest ⁵	34.29 - 40.64cm (13.5" - 16")	Set at*:		
FR15-02	☐	80° Swing-in/swing-out footrest ⁵	36.83 - 49.53cm (14.5" - 19.5")	Set at*:		
FR15-03	☐	80° Swing-in/swing-out footrest ⁵	39.37 - 52.07cm (15.5" - 20.5")	Set at*:		
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension) ¹	19.05 - 27.94cm (7.5" - 11")	Set at*:		
FR07	☐	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}	24.13 - 33.02cm (9.5" - 13")	Set at*:		
FR08	☐	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}	24.13 - 36.83cm (9.5" - 14.5")	Set at*:		
FOOTREST OPTIONS						
FR09	☐	High mount (footplate attached to the hanger with clamps) ³	12.70 - 34.29cm (5" - 13.5")	Set at*:		
FR10	☐	Elevating legrest with calf pad, length adjustable ⁴ Left ☐ Right ☐	35.56 - 55.88 cm (14"-22")	Set at:		
FR99	☐	Omit Footrest				

1- Available only with two piece adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount and the elevating legrest;

* +/- 1.27 cm (1/2")

4- The elevating legrest replace the footrest side selected; Not available with one piece footplate; 5- Not available with 8" caster

FOOTPLATES						
FP01	☐	Standard - Composite ¹				
FP02	☐	Newton adjustable angle - Composite ^{*VP01}				
FP03	☐	Newton adjustable angle - Aluminum				
FP04	☐	Newton adjustable angle - Carbon fiber with rubber inserts				
FP06	☐	One piece adjustable angle flip-up - Aluminum				
FP07	☐	One piece adjustable angle flip-up - Carbon fiber with rubber inserts				
FOOTPLATE ACCESSORIES						
FP08	☐	Length adjustable heel loop ²				
FP09	☐	Calf strap				
FP30	☐	Bodypoint Aeromesh Padded calf strap				

1- Not available with front mount; 2- not available with one piece footplate

SWING AWAY RESIDUAL LIMB SUPPORT						
FR27	☐	Residual limb support - flat 15.24x20.32cm (6"x8")	Left ☐ Right ☐			
FR28	☐	Residual limb support - flat 20.32x24.4cm (8"x10")	Left ☐ Right ☐			
FR29	☐	Residual limb support - flat 25.4x25.4cm (10"x10")	Left ☐ Right ☐			
FR30	☐	Residual limb support - flat 25.4x35.56cm (10"x14")	Left ☐ Right ☐			
FR32	☐	Residual limb support - contour 15.24x20.32cm (6"x8")	Left ☐ Right ☐			
FR33	☐	Residual limb support - contour 20.32x24.4cm (8"x10")	Left ☐ Right ☐			

Backrest									
FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE <i>(push handle included - no further selection needed)</i>									
BC01-16	☐	8° bend aluminum back cane	40.64cm (16")						
BC01-18	☐	8° bend aluminum back cane	45.72cm (18")						
BC01-20	☐	8° bend aluminum back cane	50.08cm (20")						
BC59	☐	8° bend aluminum back cane	53.34 - 60.96cm (21" - 24")	increment of 1"	Set at ² :				
FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS <i>(MUST select ONE push handle options)</i>									
BC11-912	☐	Straight aluminum back cane	22.86 - 30.48cm (9" - 12")	increment of 1/2"	Set at ² :				
BC11-1215	☐	Straight aluminum back cane	30.48 - 38.10cm (12" - 15")	increment of 1/2"	Set at ² :				
BC11-1518	☐	Straight aluminum back cane	38.10 - 45.72cm (15" - 18")	increment of 1/2"	Set at ² :				
BC11-1821	☐	Straight aluminum back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at ² :				
BC12-1315	☐	8° bend aluminum back cane	33.02 - 38.10cm (13.5" - 15")	increment of 1/2"	Set at ² :				
BC12-1518	☐	8° bend aluminum back cane	38.10 - 45.72cm (15" - 18")	increment of 1/2"	Set at ² :				
BC12-1821	☐	8° bend aluminum back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at ² :				
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS <i>(MUST select ONE push handle options)</i>									
BC13-1518	☐	Straight aluminum back cane	38.10 - 45.72cm (15" - 18")	increment of 1/2"	Set at ² :				
BC13-1821	☐	Straight aluminum back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at ² :				
BC14-1518	☐	8° bend aluminum back cane	38.10 - 45.72cm (15" - 18")	increment of 1/2"	Set at ² :				
BC14-1821	☐	8° bend aluminum back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"	Set at ² :				
		Select back cane angle:	85° ¹ ☐	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐	
PUSH HANDLES OPTIONS <i>(ONE option at a time)</i>									
PH01	☐	Push handles							
PH05	☐	Fold down push handles ^{3, 5}							
PH06-920	☐	clamp on height adjustable push handles 22.86 - 50.8cm (9" - 20") ^{3, 4, 5}							
PH06-923	☐	clamp on height adjustable push handles 22.86 - 58.4cm (9" - 23") ^{3, 4, 5}							
PH02	☐	Omit push handles							
STROLLER HANDLE AND STABILIZER BAR <i>(only available with push handle)</i>									
PH04	☐	Stroller handle <i>(not available with the HD option)</i>							
PH03	☐	Newton folding stabilizer bar <i>(not available with stroller handle)</i>							

1- Not available with 8 degree bend cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265 lb weight limit; 4- Not compatible with swing away armrests, 8° bend back canes, 21-24" back canes and amputee axle plates; 5- Not compatible with stroller handle, the stabilizer bar, headrest and attendant lock

BACK UPHOLSTERY									
BH07	☐	Slip-on Active-Mesh back upholstery							
BH08	☐	Adjustable tension Active-Mesh back upholstery							
BH99	☐	Omit Back upholstery							
SEAT CUSHION									
SC01	☐	5.08 cm (2") Cushion							
SC02	☐	7.62 cm (3") Cushion							
SC99	☐	Omit seat cushion							

Rear wheels											
ONE ARM DRIVE ¹		Only available with Newton One rear wheels permanent axles, 0° camber									
HR22	☐	One arm drive - aluminum inner handrim ²	Left ☐	Right ☐							
HR60	☐	One arm drive - plastic coated inner handrim ²	Left ☐	Right ☐							
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)									
REAR WHEELS								Available			
RW01	☐	Mag ²						20", 22" 24"			
RW02	☐	Newton One - spoke wheel ⁴						20", 22" 24"			
RW03	☐	Newton Gravity - ultralight wheel						20", 22", 24", 25", 26"			
RW04	☐	Spinergy Spox	Blk STD ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	22", 24", 25", 26"			
RW05	☐	Spinergy LX	Blk STD ☐	Pur ⁵ ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	Pink ⁵ ☐	Grn ⁵ ☐	Org ⁵ ☐
RW07	☐	Spinergy CLX ⁶	Blk STD ☐	Pur ⁵ ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	Pink ⁵ ☐	Grn ⁵ ☐	Org ⁵ ☐
TIRES								Available			
AIRLESS											
TI01	☐	Soft Urethane 3.49cm (1 3/8") - Low Tread, black ⁹						20", 22" 24"			
TI02	☐	Soft Urethane 3.49cm (1 3/8") - Medium Tread, black ^{7,9}						22", 24", 25", 26"			
TI04	☐	Soft Urethane 2.54cm (1") - No tread, black - <i>Shox</i>						22", 24", 25", 26"			
TI08	☐	Pneumatic Airless insert 3.49cm (1 3/8") - Medium Tread ⁸						22", 24", 26"			
PNEUMATIC											
TI05	☐	Pneumatic 3.49cm (1 3/8") Medium Tread ⁹						20", 22", 24", 25" 26"			
TI06	☐	Pneumatic HP 2.54cm (1") - Low Tread, Puncture resistant, black - <i>SpeedLite</i>						20", 22", 24", 25", 26"			
TI07	☐	Pneumatic HP 2.54cm (1") - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>						24", 25", 26"			
TI03	☐	Pneumatic 5.08cm (2") - Large Tread, black ¹⁰						22", 24", 25"			

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide;
2- Only standard handrim position available; 3- Not compatible with HP tires; Only available with aluminum or plastic coated handrim, standard position; 4- Not compatible with Shox tires;
5- color spokes require an additionnal 10 days lead-time; 6- Only available with 1" tires; 7- 22" is only available with Mags and Newton One rear wheels;
8- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels; 9- Not available with flip-back armrest when the 6° camber is selected;
10- Not available with armrests when the 3° or 6° camber are selected;

Handrims									
HANDRIM								Available	
HR01	☐	Aluminum anodized						20", 22", 24", 25", 26"	
HR61	☐	Aluminum black hard anodized ²						20", 22", 24", 25"	
HR02	☐	Plastic coated ¹						20", 22", 24", 25", 26"	
HR03	☐	High friction coated ^{1,2,3}						22", 24", 25", 26"	
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ²	See page 6 for the color selection and price					20", 22", 24", 25", 26"	
HR28	☐	Surge ^{1,2,3} (1 1/2" Oval handrim w/ Gripton high friction strip)						22", 24", 25", 26"	
HR05	☐	Surge LT ^{1,2,3} (1 1/8" Oval handrim w/ Gripton high friction strip)						22", 24", 25", 26"	
HR04	☐	Natural Fit ^{1,2,3}						22", 24", 25", 26"	
HR62	☐	Simi H ^{1,2,3} (1 3/8" Oval handrim with high friction strip)						22", 24", 25", 26"	
HR63	☐	Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)						22", 24", 25", 26"	
HR64	☐	Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim)						22", 24", 25", 26"	
HR99	☐	Omit handrim							

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;

Handrims option									
ASSEMBLY POSITION									
HR29	☐	Standard Handrim position	Full tab						
HR030	☐	Narrow Handrim position no cut	Full tab						
HR30	☐	Narrow Handrim position with cut ³	Short tab 2 positions						
HR31	☐	Super Narrow Handrim ³	Short tab 2 positions STD ☐	Short tab 1 position ☐					

THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS									
HR23	☐	Standard grip							
HR27	☐	Super grip							
HR099	☐	Omit thumb piece							


3- Only available with Aluminum and Air Grip handrims;

AXLE									
AX01	☐	Quick release axle ¹							
AX02	☐	Permanent axle							
AX03	☐	Quad release axle ¹							
AXLE PLATE									
AX04	☐	Amputee adjustable axle plate ²							

1- Not required with one arm drive; 2 - May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down

REAR WHEEL CAMBER									
CM00	☐	0° camber							
CM03	☐	3° camber							
CM06	☐	6° camber (not available with amputee axle plate)							

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable										Most responsive			
X =	AMPUTEE	2,5 1"	3,2 1/4"	3,8 1/2"	4,4 3/4"	5,1 2"	5,7 1/4"	6,3 1/2"	7 3/4"	7,6 3"	8,3 1/4"	8,9 1/2"	9,5 3/4"
	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 3,8 cm (1.5")

Possible conflict between a seat depth under 40.6 cm (16") and a center of gravity over 6.3 cm (2.5").

Possible conflict between a 26" rear wheel and a center of gravity over 6.3 cm (2.5").

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

1-Only available with the amputee axle plate

WHEEL LOCK

WL08	☐	Newton Wheel lock - push to lock			
WL16	☐	Aluminum push to lock with extensions			
WL02	☐	Pull to lock			
WL04	☐	Newton grade aid push to lock (works best with pneumatic or high pressure tire)			
WL13 WL14	☐	6" removable extension handle ¹ *VP01			
WL05	☐	Scissors Lock			
WL06	☐	Attendant Lock (additional set of wheel lock mounted to the back)			
WL15	☐	Unilateral wheel lock - pull to lock ²	Left ☐	Right ☐	
WL99	☐	Omit wheel locks			

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04); **Only the extension for the Newton push to lock is available with the vision package (VP01)**

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

Armrests

HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) ^{1,4} *VP01					
AR01-812	☐	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")	Locked ☐	Unlocked STD ☐
AR01-1014	☐	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")	Locked ☐	Unlocked STD ☐
AR02-812	☐	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")	Locked ☐	Unlocked STD ☐
AR02-1014	☐	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")	Locked ☐	Unlocked STD ☐
HEIGHT ADJUSTABLE "T" ARMREST ¹ *VP01					
AR13-610	☐	25.4cm (10") Desk length armpad ⁵	15.25 - 25.40cm (6" - 10")		
AR03-812	☐	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")		
AR03-1014	☐	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")		
AR14-610	☐	35.56cm (14") Full length armpad ⁵	15.25 - 25.40cm (6" - 10")		
AR04-812	☐	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")		
AR04-1014	☐	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")		
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²					
AR05	☐	Tubular armrest - Short pad	20.32 - 27.94cm (8 - 11")		
AR23	☐	Tubular armrest - Long pad	20.32 - 27.94cm (8 - 11")		
AR06	☐	25.4cm (10") Desk Length flat armpad	20.32 - 27.94cm (8 - 11")		
AR17	☐	35.56cm (14") Full Length flat armpad	20.32 - 27.94cm (8 - 11")		
ARMREST OPTIONS					
AR12	☐	35.56cm (14") Full Length Armpad - Gel Ovation (pair) ³			
AR99	☐	Omit armrest			


1- May have 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid backrests; Not available with attendant lock and fixed carbon fiber side guards;

3- Only available with full length armrest; Not available with swing away and desk length armrests; 4- Not available with attendant lock; 5- May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS (not available with Flip Back and "T" Armrest)					
SG01	☐	Plastic side guards (removable)			
SG02	☐	Carbon fiber side guards (fixed)			
SG03	☐	Carbon fiber side guards (removable)			
SG06	☐	Plastic fender side guards - straight bracket (removable)			
SG08	☐	Plastic fender side guards - offset bracket (removable) ¹			
SG04	☐	Carbon fiber fender side guards (removable)			
SG99	☐	Omit side guard			
SIDE GUARDS OPTION					
SG07	☐	Side guards clip (pair) ²			

1- only compatible with 0° camber; Increase distance between fender side guards by 1 1/8"; 2- Not available with fixed carbon fiber side guard

Fork and Stem bolt combinations



Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
	16	1.0 to 2.75
22 in	17	1.0 to 3.25
	18+	1.0 to 3.75
24 in	16	1.0 to 2.0
	17	1.0 to 2.75
	18+	1.0 to 3.75

Vision Package Accents



Notes :