



Helio C2

ORDER FORM - Medifab

1- Purchasing information

P.O. number:* Date:

Company Name:*

Purchasing contact:*

Phone:*

Email:*

Contact preference:* Phone ☐ Email ☐

Marked for:

3- Shipping information

Ship to:*

Name:

Address:

City: State:

Zip Code:

Phone:

2- ATP/Clinician information

ATP/Clinician:*

Phone:*

Email:*

Contact preference:* Phone ☐ Email ☐ * Required

MA083 ☐		MA084 ☐		MA085 ☐		MA086 ☐		MA087 ☐		MA088 ☐		MA089 ☐		
English Manual (STD)		German Manual		Italian Manual		Dutch Manual		Spanish Manual		French Manual		Portuguese Manual		
MODEL														
WC05		☐	Ultralight Folding Carbon Fiber Helio C2 Wheelchair (120 kg - 265 lb weight capacity)											
WC09		☐	HD Heavy Duty package (159 kg - 350 lb weight capacity)											
WC10		☐	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber											
Frame														
SEAT WIDTH ³ - cm (inches)														
SW		☐	35.56 (14") ☐	38.10 (15") ☐	40.64 (16") ☐	43.18 (17") ☐	45.72 (18") ¹ ☐	48.26 (19") ¹ ☐						
		☐	50.80 (20") ¹ ☐			53.34 (21") ² ☐	55.88 (22") ² ☐							
SEAT DEPTH														
SD		☐	35.56 (14") ☐	38.10 (15") ☐	40.64 (16") ☐	43.18 (17") ☐	45.72 (18") ☐	48.26 (19") ☐						
		☐	50.80 (20") ☐											
TRANSIT TIE DOWN														
AC04		☐	Transit Tie Down (Unoccupied)											
AC05-S		☐	WC-19 Transit Tie Down w/Q'straint belt ⁴				Short belt 50-104 cm (20" - 41")							
AC05-L		☐	WC-19 Transit Tie Down w/Q'straint belt ⁴				Long belt 50-116 cm (20" - 46")							
AC35		☐	WC-19 Transit Tie Down without belt ⁴											

1- Available with HD option; 2- Only available with HD option [3. Width calculator](#) 4- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

PACKAGES								
PK01	☐	Newton Ultralight Package (include: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels)						
VP01	☐	Vision Package: Orange accent added to point of contact accessories including: footrest latches, "T" and flip back armrest button caps, 6" Push to lock extension and rubber inserts on Newton composite angle adjustable footplates ¹						

1- Please select footplates, footrest, armrest, push to lock wheel lock and the extensions in the appropriate sections; To know the compatible parts, look for this logo ***VP01**; Recommended with black chair;
See picture of the product on page 6.

C2 FRONT Seat-to-floor Height - Please select your choice

Forks / Position		7.62cm (3") caster			10.16cm (4") caster			12.7cm (5") caster			15.24cm (6") caster			17.78cm (7") caster			20.32cm (8") caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
		STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"	STD	+ 2.54 cm 1"	+ 5.08cm 2"
4" 10.2 cm	P1	33.02cm 13"	35.56cm 14"	38.10cm 15"	34.29cm 13.5"	36.83cm 14.5"	39.37cm 15.5"												
	P2	34.29cm 13.5"	36.83cm 14.5"	39.37cm 15.5"	35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"									
5" 12.7 cm	P1																		
	P2				35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"									
	P3	35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"						
	P4	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"	40.64cm 16"	43.18cm 17"	45.72cm 18"	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"			
7" 17.8 cm	P1										41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"			
	P2							41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"			
	P3							43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"
	P4	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20"	48.26cm 19.5"	50.80cm 20.5"	53.3cm 21"
	P5	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	48.26cm 19"	50.80cm 20"	53.3cm 21"	49.53cm 19.5"	52.07cm 20.5"	54.6cm 21.5"

Height available +/- 0,65cm (1/4") (see stem bolt and fork combination at the end of the order form)

CASTERS				Available		
CA01-03 to CA01-08	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6", 7", 8"		
CA10 to CA11	☐	Pneumatic casters (1 ¼")		6", 8"		
CA14 to CA16	☐	Newton UltraCaster Composite (1 ½ ")		4", 5", 6"		
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 ½ ")		4", 5", 6"		
CASTER OPTIONS						
CA12	☐	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) ¹				
CA13	☐	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)				

1- Not available with WC-19

C2 - REAR Seat-to-floor Height

WHEEL SIZES		Height available		Specify desired height: _____
☐	20" (451) wheel	31.75cm / 12.5"	to 43.81cm / 17.25"	
☐	22" (501) wheel	34.29cm / 13.5"	to 46.35cm / 18.25"	
☐	24" (540) wheel	36.83cm / 14.5"	to 48.89cm / 19.25"	
☐	25" (559) wheel	38.10cm / 15"	to 50.16cm / 19.75"	
☐	26" (590) wheel	39.37cm / 15.5"	to 51.34cm / 20.25"	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1cm" of selected seat to floor height"

Footrest						
SWING-IN/SWING-OUT FOOTREST						
FR03-01	☐	60° Swing-in/swing-out footrest - 34.29 - 43.18cm (13.5" - 17")		Set at*:		
FR03-02	☐	60° Swing-in/swing-out footrest - 36.83 - 49.53cm (14.5" - 19.5")		Set at*:		
FR03-03	☐	60° Swing-in/swing-out footrest - 39.29 - 52.07cm (15.5" - 20.5")		Set at*:		
FR04-01	☐	70° Swing-in/swing-out footrest - 34.29 - 43.18cm (13.5" - 17")		Set at*:		
FR04-02	☐	70° Swing-in/swing-out footrest - 36.8 3- 49.53cm (14.5" - 19.5")		Set at*:		
FR04-03	☐	70° Swing-in/swing-out footrest - 39.37 - 52.07cm (15.5" - 20.5")		Set at*:		
FR15-01	☐	80° Swing-in/swing-out footrest - 34.29 - 40.64cm (13.5" - 16") ⁵		Set at*:		
FR15-02	☐	80° Swing-in/swing-out footrest - 36.83 - 49.53cm (14.5" - 19.5") ⁵		Set at*:		
FR15-03	☐	80° Swing-in/swing-out footrest - 39.37 - 52.07cm (15.5" - 20.5") ⁵		Set at*:		
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension) ¹	19.05 - 27.94cm (7.5" - 11")	Set at*:		
FR07	☐	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}	24.13 - 33.02cm (9.5" - 13")	Set at*:		
FR08	☐	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}	24.13 - 36.83cm (9.5" - 14.5")	Set at*:		
FOOTREST OPTIONS ³						
FR09	☐	High mount (footplate attached to the hanger with clamps) ³	12.70 - 34.29cm (5" - 13.5")	Set at*:		
FR10	☐	Elevating legrest with calf pad, length adjustable ⁴ Left ☐ Right ☐	35.56 - 55.88 cm (14"-22")	Set at:		
FR99	☐	Omit Footrest				

1- Available only with two piece adjust. footplates; 2- May interfere with front casters; 3- Not available with front mount and the elevating legrest;

4- The elevating legrest replace the footrest side selected; Not available with one piece footplate; 5- Not available with 8" caster

FOOTPLATES						
FP01	☐	Standard - Composite ¹				
FP02	☐	Newton adjustable angle - Composite ^{*VP01}				
FP03	☐	Newton adjustable angle - Aluminum				
FP04	☐	Newton adjustable angle - Carbon fiber with rubber inserts				
FP06	☐	One piece adjustable angle flip-up - Aluminum				
FP07	☐	One piece adjustable angle flip-up - Carbon fiber with rubber inserts				
FOOTPLATE ACCESSORIES						
FP08	☐	Length adjustable heel loop ²				
FP09	☐	Calf strap				
FP30	☐	Bodypoint Aeromesh Padded calf strap				

1- Not available with front mount; 2- not available with one piece footplate

SWING AWAY RESIDUAL LIMB SUPPORT						
FR27	☐	Residual limb support - flat 15.24x20.32cm (6"x8")	Left ☐ Right ☐			
FR28	☐	Residual limb support - flat 20.32x24.4cm (8"x10")	Left ☐ Right ☐			
FR29	☐	Residual limb support - flat 25.4x25.4cm (10"x10")	Left ☐ Right ☐			
FR30	☐	Residual limb support - flat 25.4x35.56cm (10"x14")	Left ☐ Right ☐			
FR32	☐	Residual limb support - contour 15.24x20.32cm (6"x8")	Left ☐ Right ☐			
FR33	☐	Residual limb support - contour 20.32x24.4cm (8"x10")	Left ☐ Right ☐			

Backrest									
FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)									
BC06-912		Straight Carbon fiber back cane	22.86 - 30.48cm (9" - 12")	increment of 1/2"		Set at ² :			
BC06-1215		Straight Carbon fiber back cane	30.48 - 38.10cm (12" - 15")	increment of 1/2"		Set at ² :			
BC06-1518		Straight Carbon fiber back cane	38.10 - 45.72cm (15" - 18")	increment of 1/2"		Set at ² :			
BC06-1821		Straight Carbon fiber back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"		Set at ² :			
BC08-1315		8° Bend Carbon fiber back cane	34.29 - 38.1cm (13.5" - 15")	increment of 1/2"		Set at ² :			
BC08-1518		8° Bend Carbon fiber back cane	38.1 - 45.72cm (15" - 18")	increment of 1/2"		Set at ² :			
BC08-1821		8° Bend Carbon fiber back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"		Set at ² :			
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)									
BC09-1215	☐	Straight Carbon fiber back cane	30.48 - 38.10cm (12" - 15")	increment of 1/2"		Set at ² :			
BC09-1518	☐	Straight Carbon fiber back cane	38.10 - 45.72cm (15" - 18")	increment of 1/2"		Set at ² :			
BC09-1821	☐	Straight Carbon fiber back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"		Set at ² :			
BC10-1315	☐	8° bend Carbon fiber back cane	34.29 - 38.1cm (13.5" - 15")	increment of 1/2"		Set at ² :			
BC10-1518	☐	8° bend Carbon fiber back cane	38.1 - 45.72cm (15" - 18")	increment of 1/2"		Set at ² :			
BC10-1821	☐	8° bend Carbon fiber back cane	45.72 - 53.34cm (18" - 21")	increment of 1/2"		Set at ² :			
		Select back cane angle:	85° ¹ ☐	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐	
PUSH HANDLES OPTIONS (ONE option at a time)									
PH01	☐	Push handles							
PH05	☐	Fold down push handles ^{3, 5}							
PH06-920	☐	clamp on height adjustable push handles 22.86 - 50.8cm (9" - 20") ^{3,4, 5}							
PH06-923	☐	clamp on height adjustable push handles 22.86 - 58.4cm (9" - 23") ^{3,4, 5}							
PH02	☐	Omit push handles							
STROLLER HANDLE AND STABILIZER BAR (only available with push handle)									
PH04	☐	Stroller handle (not available with the HD option)							
PH03	☐	Newton folding stabilizer bar (not available with stroller handle)							

1- Not available with 8 degree bend cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 120kg (265 lb) weight limit

4- Not compatible with swing away armrests, 8° bend back canes and amputee axle plates; 5- Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

BACK UPHOLSTERY									
BH07	☐	Slip-on Active-Mesh back upholstery							
BH08	☐	Adjustable tension Active-Mesh back upholstery							
BH99	☐	Omit Back upholstery							
SEAT CUSHION									
SC01	☐	5.08 cm (2") Cushion							
SC02	☐	7.62 cm (3") Cushion							
SC99	☐	Omit seat cushion							

Rear wheels												
ONE ARM DRIVE ¹			Only available with Newton One rear wheels, permanent axles, 0° camber									
HR22	☐	One arm drive - aluminum inner handrim ²	Left ☐	Right ☐								
HR60	☐	One arm drive - plastic coated inner handrim ²	Left ☐	Right ☐								
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)										
REAR WHEELS					Available							
RW01	☐	Mag ³								20", 22" 24"		
RW02	☐	Newton One - spoke wheel ⁴								20", 22" 24"		
RW03	☐	Newton Gravity - ultralight wheel								20", 22", 24", 25", 26"		
RW04	☐	Spinergy Spox	Blk STD ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐			22", 24", 25", 26"		
RW05	☐	Spinergy LX	Blk STD ☐	Pur ⁵ ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	Pink ⁵ ☐	Grn ⁵ ☐	Org ⁵ ☐	22", 24", 25", 26"
RW07	☐	Spinergy CLX ⁶	Blk STD ☐	Pur ⁵ ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	Pink ⁵ ☐	Grn ⁵ ☐	Org ⁵ ☐	24", 25"
TIRES					Available							
AIRLESS												
TI01	☐	Soft Urethane 3.49cm (1 3/8") - Low Tread, black ⁹								20", 22" 24"		
TI02	☐	Soft Urethane 3.49cm (1 3/8") - Medium Tread, black ^{7,9}								22", 24", 25", 26"		
TI04	☐	Soft Urethane 2.54cm (1") - No tread, black - <i>Shox</i>								22", 24", 25", 26"		
TI08	☐	Pneumatic Airless insert 3.49cm (1 3/8") - Medium Tread ⁸								22", 24", 26"		
PNEUMATIC												
TI05	☐	Pneumatic 3.49cm (1 3/8") Medium Tread ⁹								20", 22", 24", 25", 26"		
TI06	☐	Pneumatic HP 2.54cm (1") - Low Tread, Puncture resistant, black - <i>SpeedLite</i>								20", 22", 24", 25", 26"		
TI07	☐	Pneumatic HP 2.54cm (1") - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>								24", 25", 26"		
TI03	☐	Pneumatic 5.08cm (2") - Large Tread, black ¹⁰								22", 24", 25"		

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; *Only available with 24" wheel size; Not available with 21" and 22" wide*; 2- Only standard handrim position available;

3- Not compatible with HP tires; *Only available with aluminum or plastic coated handrim, standard position*; 4- Not compatible with Shox tires; 5- color spokes require an additionnal 10 days lead-time;

6- Only available with 1" tires; 7- 22" is only available with Mags and Newton One rear wheels; 8- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels

9- Not available with flip-back armrest when the 6° camber is selected; 10- Not available with armrests when the 3° or 6° camber are selected;

Handrims					
HANDRIM				Available	
HR01	☐	Aluminum anodized		20", 22", 24", 25", 26"	
HR61	☐	Aluminum black hard anodized ²		20", 22", 24", 25"	
HR02	☐	Plastic coated ¹		20", 22", 24", 25", 26"	
HR03	☐	High friction coated ^{1,2,3}		22", 24", 25", 26"	
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ² <i>See on page 6 for color choice and price</i>		20", 22", 24", 25", 26"	
HR28	☐	Surge ^{1,2,3} (1 1/2" Oval handrim w/ Gripton high friction strip)		22", 24", 25", 26"	
HR05	☐	Surge LT ^{1,2,3} (1 1/8" Oval handrim w/ Gripton high friction strip)		22", 24", 25", 26"	
HR04	☐	Natural Fit ^{1,2,3}		22", 24", 25", 26"	
HR62	☐	Simi H ^{1,2,3} (1 3/8" Oval handrim with high friction strip)		22", 24", 25", 26"	
HR63	☐	Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)		22", 24", 25", 26"	
HR64	☐	Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim)		22", 24", 25", 26"	
HR99	☐	Omit handrim			

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;

Handrims option					
ASSEMBLY POSITION					
HR29	☐	Standard Handrim position	Full tab		
HR030	☐	Narrow Handrim position no cut	Full tab		
HR30	☐	Narrow Handrim position with cut ³	Short tab 2 positions		
HR31	☐	Super Narrow Handrim ³	Short tab 2 positions STD ☐ Short tab 1 position ☐		
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS					
HR23	☐	Standard grip			
HR27	☐	Super grip			
HR099	☐	Omit thumb piece			


3- Only available with Aluminum and Air Grip handrims;

AXLE ¹					
AX01	☐	Quick release axle ¹			
AX02	☐	Permanent axle			
AX03	☐	Quad release axle ¹			
AXLE PLATE					
AX04	☐	Amputee adjustable axle plate ²			

1- Not required with one arm drive; 2 - May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down

REAR WHEEL CAMBER					
CM00	☐	0° camber			
CM03	☐	3° camber			
CM06	☐	6° camber (not available with amputee axle plate)			

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable													
AMPUTEE	2,5 1"	3,2 1/4"	3,8 1/2"	4,4 3/4"	5,1 2"	5,7 1/4"	6,3 1/2"	7 3/4"	7,6 3"	8,3 1/4"	8,9 1/2"	9,5 3/4"	10,1 4"
☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 3,8 cm/1.5"

Possible conflict between a seat depth under 40.6 cm (16") and a center of gravity over 6,3 cm (2,5").

Possible conflict between a 66 cm (26") rear wheel and a center of gravity over 6,3 cm (2,5").

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

1-Only available with the amputee axle plate

WHEEL LOCK					
WL08	☐	Newton Wheel lock - push to lock			
WL16	☐	Aluminum push to lock with extensions			
WL02	☐	Pull to lock			
WL04	☐	Newton grade aid push to lock (works best with pneumatic or high pressure tire)			
WL13 WL14	☐	6" removable extension handle ¹ *VP01			
WL05	☐	Scissors Lock			
WL06	☐	Attendant Lock (additional set of wheel lock mounted to the back)			
WL15	☐	Unilateral wheel lock - pull to lock ²	Left ☐ Right ☐		
WL99	☐	Omit wheel locks			

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04); *Only the extension for the Newton push to lock is available with the vision package (VP01)*

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

Armrests								
HEIGHT ADJUSTABLE FLIP BACK ARMREST <i>(Convertible to "T" armrest)</i> ^{1,4} *VP01								
AR01-812	☐	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")	Locked ☐	Unlocked STD ☐			
AR01-1014	☐	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")	Locked ☐	Unlocked STD ☐			
AR02-812	☐	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")	Locked ☐	Unlocked STD ☐			
AR02-1014	☐	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")	Locked ☐	Unlocked STD ☐			
HEIGHT ADJUSTABLE "T" ARMREST ¹ *VP01								
AR13-610	☐	25.4cm (10") Desk length armpad ⁵	15.25 - 25.40cm (6" - 10")					
AR03-812	☐	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")					
AR03-1014	☐	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")					
AR14-610	☐	35.56cm (14") Full length armpad ⁵	15.25 - 25.40cm (6" - 10")					
AR04-812	☐	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")					
AR04-1014	☐	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")					
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²								
AR05	☐	Tubular armrest - short pad	20.32 - 27.94cm (8 - 11")					
AR23	☐	Tubular armrest - long pad	20.32 - 27.94cm (8 - 11")					
AR06	☐	25.4cm (10") Desk Length flat armpad	20.32 - 27.94cm (8 - 11")					
AR17	☐	35.56cm (14") Full Length flat armpad	20.32 - 27.94cm (8 - 11")					
ARMREST OPTIONS								
AR12	☐	35.56cm (14") Full Length Armpad - Gel Ovation (pair) ³						
AR99	☐	Omit armrest						
¹⁻ May have 1/2" adjustment restriction on both limits of the range; ²⁻ Could interfere with some rigid backrests; Not available with attendant lock and fixed carbon fiber side guards;								
³⁻ Only available with full length armrest; Not available with swing away and desk length armrests; ⁴⁻ Not available with attendant lock; ⁵⁻ May have some restriction with seat-to-floor height and rear wheel								
SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>								
SG01	☐	Plastic side guards <i>(removable)</i>						
SG02	☐	Carbon fiber side guards <i>(fixed)</i>						
SG03	☐	Carbon fiber side guards <i>(removable)</i>						
SG06	☐	Plastic fender side guards - straight bracket <i>(removable)</i>						
SG08	☐	Plastic fender side guards - offset bracket <i>(removable)</i> ¹						
SG04	☐	Carbon fiber fender side guards <i>(removable)</i>						
SG99	☐	Omit side guard						
SIDE GUARDS OPTION								
SG07	☐	Side guards clip (pair) ²						
¹⁻ only compatible with 0° camber; Increase distance between fender side guards by 1 1/8"; ²⁻ Not available with fixed carbon fiber side guard								
POSITIONING BELT								
SB01-S	☐	Velcro adjustable - Small 88.9cm (35") ¹						
SB01-L	☐	Velcro adjustable - Long 152.40cm (60")						
SB02-S	☐	Auto buckle - Small 88.9cm (35") ¹						
SB02-M	☐	Auto buckle - Medium 121.92cm (48")			New			
SB02-L	☐	Auto buckle - Long 152.40cm (60")						
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 142.24cm (56")						
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 154.94cm (61")						
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 129.54cm (51")						
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 142.24cm (56")						
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 154.94cm (61")						
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 129.54cm (51")						
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 142.24cm (56")						
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 154.94cm (61")						
SB03	☐	Hardware attachment for 2 extra points (1 pair)						
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 73.66 (18")						
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 93.98cm (24")						
SB99	☐	Omit positioning belt						
¹⁻ Not recommended for seat width of 45.72cm (18") and more								
HEADREST ¹								
AC01-S	☐	Profiled Headrest - Small 8.8x15.2cm (3.5"x6")						
AC01-M	☐	Profiled Headrest - Medium 11.4x17.8cm (4.5"x7")						
AC01-L	☐	Profiled Headrest - Large 14x2.3cm(5.5"x8")						
AC28	☐	Headrest mounting hardware with removable non folding stabilizer bar ¹						

¹⁻ Only available with push handles - Not available with a stroller handle or a folding stabilizer bar

ACCESSORIES													
AC03	☐	Anti-tippers											
AC29	☐	Swing-away anti-tipper ¹	Left ☐	Right ☐									
AC27	☐	Tip assist ²											
AC06	☐	Cane and crutch holder	Left ☐	Right ☐									
AC07	☐	Oxygen tank holder ⁴ - D and E Type - installed on the left											
AC08	☐	Spoke guards ³	Available: 20", 22", 24", 26"										
AC09	☐	Footrest neoprene impact guard (pair) - Black											
AC10	☐	Seat pouch											
AC30	☐	Newton multi-purpose pouch											
AC26	☐	BodyPoint Mobility Bag											
AC23	☐	Newton Backpack											
AC12	☐	Velcro Closing Strap											
AC13	☐	Tool kit											
AC16	☐	Motion Composites Polo Shirt	W/S ☐	W/M ☐	W/L ☐								
			M/S ☐	M/M ☐	M/L ☐	M/XL ☐	M/XXL ☐	M/XXXL ☐					

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 34,3cm (13 1/2") and 50,2cm (19 3/4"); 2- Cannot be combined wi


3- Not available with Mag wheels and Spinergy; 4- Not compatible with swing away armrest or amputee axle plate and WC-19 transit tie-down

Colors																				
Options ¹	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)
																				
AirGrip	☐		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Frame ²		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
FC99	☐	Custom Color	Specify PPG Paintit color: additional 10 days lead-time																	

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

2- With carbon fiber crossbrace;

Fork and Stem bolt combinations



Vision Package Accents



Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
	16	1.0 to 2.75
22 in	17	1.0 to 3.25
	18+	1.0 to 3.75
24 in	16	1.0 to 2.0
	17	1.0 to 2.75
	18+	1.0 to 3.75

Notes :